




























Centrifuges-Large Size C

Models: CN-650R, CN-650RD, CN-650, CN-5101, CN-5301, Optional Rotors & Adapters For Centrifuges

Applicable Model	Rotor		Max. Speed (rpm)	Max. RCF (xg)	Capacity (ml x pcs)	Tube Size (φ x l)mm		Cat. No.		
						Max.	Min.			
CN-650R CN-650RD CN-650 CN-5301	RS-250		3,500	2,250	1.5x96	10.7x100	8x41	2-6000-21		
					6x112	12x110	8.5x60	2-6000-22		
					10x80	15.8x100	12x60	2-6000-23		
					15x40TC	16.5x122	14x60	2-6000-25		
					15x48	17.2x122	14.5x60	2-6000-24		
					50x16TC	29.5x125	27.5x50	2-6000-26		
					250x4	62x130	58x70	2-6000-27		
									2-6000-21	
	RS-1548		3,500	2,430	15x48	17.5x107	14x80	2-6000-12		
	RS-2504		3,500	2,480	250x4	65.5x130	61x90	2-6000-13		
CN-650R CN-650RD CN-650 CN-5101 CN-5301	RS-9602		2,500	1,060	96micro x4	129x86	126x83	2-6000-14		
	RS-100		3,500	2,090	6x56	13x105	10x50	2-5400-21		
			3,500	2,150	10x48	16x107	13x50	2-5400-22		
			3,500	2,150	15x16TC	16.7x135	13.5x86	2-3400-32		
			3,500	2,150	15x32	17.5x115	14x80	2-5400-23		
			4,000	2,770	50x4	35x135	32x90	2-3400-22		
			5,000	4,640	50x4TC	29.5x137	26.5x90	2-3400-23		
			3,500	2,120	50x8	35x120	32x90	2-5400-24		
			3,500	2,120	50x8TC	29.5x127	26.5x90	2-5400-25		
			4,000	2,930	100x4	45.5x130	42x95	2-3400-21		
										2-5400-21
	RS-50		4,500	3,350	15x16	16.7x115	13.5x86	2-3400-32		
			4,000	2,600	50x4	35x117	32x90	2-3400-22		
			5,500	5,270	50x4TC	29.5x120	26.5x90	2-3400-23		
			5,000	4,200	100x4	45x112	42x95	2-3400-21		
									2-3400-32	
	RA-5008		6,000	4,500	50x8	29.5x118	27.5x90	2-6000-16		
	RA-1512H		6,000	4,500	15x12	16.7x130	13.5x86	2-6000-15		