

## C Centrifuges-Medium Size

**Model: CN-3101, Economy Centrifuge, up to 5000rpm, Max. Capacity 400ml**



**Features:**

LCD speed display • Regressive 30 min. timer or continuous operation • Stable and stepless speed change • Stainless steel inside bowl • Breakdown indication and fuse safety device • Auto current-restraining function • Overheat and overload safety protection • Hidden auto-balance device • Safety door lock and cover open indication • Auto accelerating & decelerating functions. • Auto-cooling air circulation • See-through window • Optional rotor selections for various tube size (See in page 19 in the full product catalog) • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~5,000	3,790	400	0~30 or "ON"	20.5	36×48×34

**Model: CN-3301, Brushless General Purpose Centrifuges, up to 6500rpm, Max. Capacity 400ml**



**Features:**

LCD speed display • Regressive 30 min. timer or continuous operation • Quiet brushless motor • Stable and stepless speed change • Stainless steel inside bowl • Breakdown indication & fuse safety device • Auto current-restraining function • Overheat and overload safety protection • Hidden auto-balance device • Safety door lock and cover open indication • Auto accelerating & decelerating functions • Auto-cooling air circulation • See-through window • Optional rotor selections for various tube size • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~6,500	4,500	400	0~30 or "ON"	18	36×48×34

**Model: CN-350, Medium Size Centrifuge, up to 6500rpm, Max. Capacity 400ml**



**Features:**

LCD, speed and time display with RCF setting • Quiet brushless motor • Microprocessor control • 10 memory programs for preset parameters • Auto power off when open • Auto cover lock with motor running • Overheat and overload safety protection • Auto power off for rotor imbalance • Auto buzzer when finish • Hidden auto-balance device • Power pilot light and fuse safety device • Auto accelerating & decelerating functions • See-through window • Optional rotor selections for various tube size • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~6,500	4,500	400	99:59	20.5	36×48×34

**Model: CN-350RD, Medium Size Refrigerated Centrifuge, up to 5000rpm, Max. Capacity 480ml**














**Features:**

LCD, temperature, speed and RCF setting also display • Temp. range : 0°C~40°C • 10 memory programs for preset parameters • Overheat and overload safety protection • Stainless steel inside bowl • Microprocessor control • Auto power off for rotor imbalance • Auto buzzer when finish • High function & high effective control system • Breakdown indication & fuse safety device • Auto accelerating & decelerating functions • Auto-cooling air circulation • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~5,000	4,640	480	0~30 or "ON"	40	45×56×38

## Optional Rotors & Adapters For Centrifuges CN-350RD/350/3101/3301

Rotor		Max. Speed (rpm)	Max. RCF (x g)	Capacity (ml x pcs)	Tube Size (x 1)mm Max. Min.		Cat. No.
RS-50		3,500	1,990	100x4	45x112	42x95	2-3400-21
		4,000	2,600	50x4	35x117	32x90	2-3400-22
		4,000	2,600	50x4TC	29.5x120	26.5x90	2-3400-23
		3,500	1,990	15x16	16.7x115	13.5x86	2-3400-32
		3,500	1,990	15x12TC	17x120	15x85	2-3400-24
		3,500	1,990	10x24	15.5x113	14x85	2-3400-25
		3,500	1,700	6x40	12.7x88	11x63	2-3400-26
		3,500	1,700	6x36	12x106	11x85	2-3400-27
		3,500	1,330	1.5x48	11x56	10x36	2-3400-28
							
RS-1508		3,500	1,990	15x8	17x122	15x90	2-3000-14
RA-1508		6,000	3,940	15x8	16.7x122	15x90	2-5052-11
RA-1512H		6,000	4,500	15x12	16.7x122	15x90	2-6000-15
RA-1516		4,000	2,250	15x16	16.7x115	15x90	2-3000-12
RA-1024		3,500	1,840	10x24	16.7x107	15x90	2-3000-13
RA-5004		6,000	3,980	50x4	29.5x130	27.5x94	2-5052-12
RA-5006		6,000	3,980	50x6	29.5x130	27.5x90	2-5052-22
RA-5008		6,000	4,500	50x8	29.5x117	27.5x90	2-1000-12
REA-1520		6,000	3,140	1.5x20	11x50	10x40	2-5052-13