



Carbon Dioxide (CO2)

\$82.00

[View Product](#)



VGM-10

- 100PPM Acetylene
- 100PPM Carbon Dioxide
- 100PPM Carbon Monoxide
- 100PPM Ethane
- 100PPM Ethylene
- 100PPM Hydrogen
- 100PPM Methane
- 100PPM Nitrogen
- 100PPM Oxygen
- 100PPM Propane
- Balance: Argon

VGM-10

\$485.00

[View Product](#)

VGM-9

30.07% Hydrogen
1.24% Propane
1.24% Propylene
1.25% Isobutane
1.25% Butane
1.25% Trans-2-Butene
1.25% Cis-2-Butene
1.25% Ethylene
1.26% Ethane
1.25% Argon
1.24% Nitrogen
1.25% Methane
1.25% Carbon Monoxide
Balance: Carbon Dioxide

VGM-9

\$655.00

[View Product](#)

VGM-8

0.25% Nitrogen
0.03% Carbon Dioxide
5.58% Ethane
0.03% Propane
0.006% Isobutane
0.01% Butane
0.01% Isopentane
0.01% Neopentane
0.05% Hexane
Balance: Methane

VGM-8

\$540.00[View Product](#)

VGM-7

0.10% Oxygen	0.05% Pentane
1.00% Nitrogen	0.025% Hexane
1.00% Carbon Dioxide	0.025% Heptane
2.52% Ethane	Balance: Methane
1.00% Helium	
3.00% Hydrogen	
0.30% Propane	
0.10% Isobutane	
0.10% Butane	
0.05% Isopentane	

VGM-7

\$585.00[View Product](#)

VGM-5

141.5 PPM	1,1-Dichloroethylene
144.5 PPM	1,1-Dichloroethane
141.5 PPM	1,2-Dichloroethane
258.9 PPM	1,3-Butadiene
183.0 PPM	2-Chloropropene
144.5 PPM	Cis-1,2-Dichloroethylene
144.5 PPM	Trans-1,2-Dichloroethylene
131.0 PPM	Vinyl Bromide
538.0 PPM	Vinyl Acetylene
868.4 PPM	Ethyl Chloride
1664.0 PPM	Methyl Chloride
	Balance: Nitrogen

VGM-5

\$570.00

[View Product](#)



VGM-4

0.10% Methane
0.10% Ethane
0.10% Propane
0.10% Butane
0.10% Isobutane
0.10% Pentane
0.10% Isopentane
20.88% Oxygen
Balance: Nitrogen

VGM-4

\$395.00

[View Product](#)



VGM-3

- 0.15% Acetylene**
- 0.15% Carbon Dioxide**
- 0.15% Propylene**
- 0.15% Hydrogen**
- 0.15% Carbon Monoxide**
- 0.15% Ethane**
- 0.15% 1-Butene**
- 0.15% Propane**
- 0.15% Methane**
- 0.15% Etylene**
- Balance Nitrogen**

VGM-3

\$490.00

[View Product](#)



VGM-2

- | | |
|-------------------------------|----------------------------------|
| 100 PPM Methane | 100 PPM Cis-2-Butene |
| 100 PPM Ethane | 100 PPM 1,3-Butadiene |
| 100 PPM Ethylene | 100 PPM Isopentane |
| 100 PPM Acetylene | 100 PPM Pentane |
| 100 PPM Propane | 100 PPM 1-Pentene |
| 100 PPM Propylene | 100 PPM Trans-2-Pentene |
| 100 PPM Isobutane | 100 PPM Cis-2-Pentene |
| 100 PPM Propadiene | 100 PPM 2-Methyl-2-Butene |
| 100 PPM Butane | 100 PPM Hexane |
| 100 PPM 1-Butene | Balance: Helium |
| 100 PPM Isobutylene | |
| 100 PPM Trans-2-Butene | |

VGM-2

\$980.00

[View Product](#)



VGM-1

0.25% Methane
0.25% Ethane
0.25% Propane
0.25% Butane
0.25% Isobutane
0.10% Pentane
0.10% Isopentane
0.10% Hexane
Balance: Nitrogen

VGM-1

\$395.00

[View Product](#)



Cylinder Devalving Tool

\$152.60

[View Product](#)



Tedlar Span Bag 9x9 Double Valve

\$69.00

[View Product](#)



Gas Tubing

\$5.10

[View Product](#)



**Compressed
Gas Hazmat fee
paid to carrier**

A one-time \$60 hazardous material fee payable to the carrier is required on any shipment containing compressed gas.

\$60.00

[View Product](#)