


**Membrane Solutions**



## **MS® Serological Pipette-GRADUATED**

*ISO9001: 2000; FDA registered; CE certification*



**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd  
*Australian Distributors; Importers & Manufacturers*

MS® Serological Pipette fulfil all requirements set for disposable serological pipettes: **crystal clear** transparency, bidirectional graduations that provide a choice of ascending and descending volume reading, and accurate anti-drip delivery.

The cotton wadding sits firmly in the tube hinders an over-filling of the pipette sucking apparatus. The **colour**-coding on the pipettes assists in avoiding confusion during routine examinations and make it simple to determine the size/volume of each pipette.

Pipettes are pyrogen free, non cytotoxic, non haemolytic and come single wrapped in a convenient blister package. All the features guarantee them can be safely used in cell and tissue culture, bacteriological, and clinical research.

**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd  
*Australian Distributors; Importers & Manufacturers*



Sterile Serological Pipette



The same from their back

Items	Capacity/ Subdivisions (ml)	Color Code	Sterile
LBGSP000001	1 / 0.01	YELLOW	Yes
LBGSP000002	2 / 0.01	GREEN	Yes
LBGSP000005	5 / 0.1	BLUE	Yes
LBGSP000010	10 / 0.1	RED	Yes
LBGSP000025	25 / 0.2	PURPLE	Yes
LBGSP000050	50 / 1.0	BLACK	Yes

Certified as having been sterilized by gamma radiation

Items	Capacity/ Subdivisions (ml)	Color Code	Sterile
LBGSP001001	1 / 0.01	YELLOW	NO
LBGSP001002	2 / 0.01	GREEN	NO
LBGSP001005	5 / 0.1	BLUE	NO
LBGSP001010	10 / 0.1	RED	NO
LBGSP001025	25 / 0.2	PURPLE	NO
LBGSP001050	50 / 1.0	BLACK	NO

## Membrane Solutions, LLC