

Crimping Tools for Vials : **Crimp & Decappers**

Hand Crimper

Fixed Head Standard and Flip Off Cap – Type
Standard

Adjustable Head (for height NOT Cap Diameter)



All Stainless Steel Hand Crimper

Autoclavable; Removable Rubber Grips



New 2012 mini **Pneumatic Hand Air Crimper**

- requires Air compressor



Bench Press Crimper

- Lever operated



All Size heads available

Crimping and Decapping

Serum Vials :13mm, 20mm also 32mm

Chromatography : 8, 11, 13, 20mm Caps

Air & Bench Press Interchangeable Heads

Hand Crimp : Fixed Head Only

Hand Crimpers - Non Adjustable (Standard)

13mm, 20mm in stock Melbourne : for Serum Vials

also 8mm, 11mm, & 32mm for chromatography vials

To order . . .

Adjustable Hand Crimper 13mm, 20mm; Pneumatic Hand Air Crimper

**Adjustable
Hand
Crimper**



**Air
Crimper**
with
**detachable
Heads**
including
Decapper
Heads



The **Adjustable Hand Crimper** has an adjustment feature which will allow it to be used with a variety of septa, discs, stoppers or liners with a thickness range of 0.25 to 3.0mm. The unit looks like the standard hand crimper but it is supplied with a wrench and special instructions for making adjustments.

The user may adjust the stroke or movement for each application as needed.

Available for 13mm and 20mm Diameter Crimp Vials both Standard Caps and Flip-off Cap types

ALL Stainless Steel Hand Crimper

- Autoclavable
- Removable Rubber Grips



Ordering ;

CT-13/20Crimp-SA and CT-13/20Crimp-FOA adjustable

CT-13/20Crimp-S and CT-13/20Crimp-FO NON-adjustable : **NEW Pricing \$259.00 each**

NEW 2012 All Stainless Steel Hand : Crimper CT-13/20Crimp-S and FO : Pricing \$327.00 each
Autoclavable . . . with removable Rubber Grips

Hand Pneumatic Crimper / Detachable Heads CT-13/20AirCrimp \$625.00 each

Order Heads separately CT-13/20AirHead-S/FO \$250.00 each

: Requires separate 70psi Air Supply and Pressure regulator / hose



Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA