

Gas Syringes

Series C

- ▶ **Optional Chaney adapter** – much lower dead volume than luer-style connectors
- ▶ **Perfect syringe for everyday use** – all our standard features, but without a valve for pressurized sample storage



	100 µl	250 µl	500 µl	1 ml	2 ml	5 ml	10 ml
Standard:	030025	030031	030032	030033	030034	030035	030036
Std. w/Chaney adapter:	101430	101440	101448	101455	101463	101470	101478
Luer:	030045	030051	030052	030053	030054	030055	030056

Replacement needles (Luer needles on page 11)

	Standard	Chaney adapter
Size:	.029" x .005" x 2"	.018" x .008" x .75"
Type:	Bevel, open end	Side port, taper
Pkg. of 3:	943050	101935

Add sample storage capability to any syringe with a Mininert syringe valve