


Series C-120

Fixed needle: .019" x .005" x 2", bevel, open end

- ▶ **Fixed needle** – zero dead volume connection
- ▶ **Leak-tight to 250 psi** – liquids and gases
- ▶ **Extended syringe barrel** – attach a bulb or suction line for easier cleaning and filling



In addition to standard Pressure-Lok features such as PTFE plunger tip and needle seat, the C-120 offers as standard equipment a barrel extension beyond the rear flange which permits attachment of a vacuum bulb or line for easier filling or cleaning.

5 µl	10 µl	25 µl	50 µl	100 µl
120021	120022	120023	120024	120025