


Liquid Syringes

Purge and Trap Syringe

Accept standard Luer hub needles

- ▶ **Luer lock** – for use with purge and trap analysis
- ▶ **Frosted glass on syringe barrel** – allows easy volume reading and accepts writing
- ▶ **Plunger tip of stress-formed virgin PTFE** – self-lubricating and durable
- ▶ **Heavy duty rear flange** – flat edges keep the syringe from rolling



The heavy duty construction of our purge and trap syringe helps eliminate costs of replacement and injuries due to breakage. All of our environmental syringes are repairable.

For this syringe, Luer needles must be ordered separately. (All of our other Luer syringes include a needle.)

5 ml	10 ml
030039	030040

Luer needles

Size:	.028" x .016" x 2"
Type:	Bevel, open end
Pkg. of 3:	943061