

# SmartSite



80220  
Needleless Injection Site,  
Swabbable Female Luer...  
13 mm



90406  
Needleless Injection  
Site,  
Swabbable Female  
Luer...



80151  
Needleless Injection Site,  
Swabbable Female Luer...  
0.142 inch ID (3.6mm)




80211  
Needleless Injection Site,  
Swabbable Female Luer...  
fits 1/4 ID Tubing (0.25 inch  
6.4mm)



**80150**  
Needleless Injection Site,  
Swabbable Female Luer...



**80134**  
Needleless Injection Site,  
Swabbable Female Luer...  
both ports fit 0.268 inch OD  
Tubing (6.8 mm)



**17629**  
Needleless Access T Connector,  
Female Luer Lock Fits 0.143 inch  
OD Tubing (3.6 mm)



**80212**  
Needleless Injection Site,  
Neutral Pressure



**80167**  
Needleless Injection Site,  
Swabbable Female Luer...



**80209**  
Needleless Injection Site,  
Swabbable Female Luer...



**80221**  
Needleless Injection Site,  
Swabbable Female Luer...  
20 mm



**80210**  
Needleless Injection Site,  
Swabbable Female Luer...  
Fits 1/9 inch ID Tubing (0.125 inch,  
3.2mm)



91600  
Needleless Injection Site,  
Swabbable Female Luer...



80189  
Needleless Injection Site,  
Swabbable Female Luer...



90403  
Swabbable Needleless One-  
Way Sampling Valve



Q60001S  
Injection Site,  
Sample Assortment Kit Organiser 10  
1/2 inch x 6 3/4 inch x 1 7/8 inch

**94629**

Urinary Connector with Non-Lated  
Sampling Port and . . .  
Fits 0.394 inch -0.511 inch Tubing x  
0.394 inch ID (10mm - 13mm x10 mm,  
Side Arm Air Vent



**99778**

4-Way Stopcock, Femal Luer Lock  
Swabbabe Luer



**99777**

4-Way Stopcock, Female Luer Lock  
Swabbabe Luer



**80365**

T-Connector, Luer Activated, ,  
Female Luer Lock 0.156 inch ID x  
0.25 inch )D (3.9 mm x 6.4mm)





**80364**  
Y-Connector, Luer Activated, Female  
Luer Lock  
0.156 inch ID x 0.25 inch OD (3.9 mm  
x 6.4 mm)



**80363**  
Luer Activated Valve, Female Luer  
Lock, Male Luer

**80362**  
Luer Activated, Female Luer Lock  
0.295 inch OD (7.5 mm )

