



Metal Luer Connectors

Nickel Plated Brass

Male Luer to I.D. Tubing Connectors

Connectors Are Also Available in Stainless Steel & Customized Configurations

Metal Luer Connectors



Luer Slip

99154

Fits 1/16"-3/32" I.D. Tubing
(1.6mm-2.4mm)



Luer Lock

99156

Fits 1/16"-3/32" I.D. Tubing
(1.6mm-2.4mm)



Luer Slip

99152

Fits 1/8"-3/16" I.D. Tubing
(3.2mm-4.8mm)



Luer Slip

99153

Fits 1/8"-3/16" I.D. Tubing
(3.2mm-4.8mm)



Luer Lock

99155

Fits 1/8"-3/16" I.D. Tubing
(3.2mm-4.8mm)



Luer Slip

99150

Fits 1/4" I.D. Tubing
(6.4mm)



Luer Slip

99151

Fits 1/4" I.D. Tubing
(6.4mm)



Luer Lock

99160

Fits 1/4"-5/16" I.D. Tubing
(6.4mm-8mm)

Female Luer to I.D. Tubing Connectors



Luer Lock

99158

Fits 1/16"-3/32" I.D. Tubing
(1.6mm-2.4mm)



Luer Lock

99157

Fits 1/8"-3/16" I.D. Tubing
(3.2mm-4.8mm)



Luer Lock

99159

Fits 3/16"-1/4" I.D. Tubing
(4.8mm-6.4mm)



99161

Female Luer Lock
To Male Luer Slip



99162

Female Luer Lock
To Male Luer Slip



99163

Female Luer Lock
To Male Record



99164

Female Record
To Male Slip