

Direct Replacement FID Collector Assembly Kit for Agilent 5890 GCs

- Constructed of high-quality stainless steel.
- Meets or exceeds manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#
FID Collector Assembly Kit (includes insulators)	19231-60690	kit	23010
FID Collector Assembly Kit w/Siltek® Ignitor Castle	—	kit	21131

Replacement FID Parts for Agilent 5890 GCs

- Meets or exceeds manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#
A) FID Collector (includes insulators)	19231-20970 19231-20960 19231-20950 19231-20940	ea.	21138
B) FID Collector Nut and Washer	5181-3311	set	21136
C) FID Ignitor*	19231-60680	ea.	21001
D) FID Ignitor Castle	19231-20910	ea.	21137
Siltek® FID Ignitor Castle	—	ea.	21135

*Also fits OI Analytical 4410 detector (similar to OI part # 191833).



SILTEK®
versions available.

Direct Replacement FID Collector Assembly Kit for Agilent 6890/6850/7890 GCs

- Constructed of high-quality stainless steel.
- Meets or exceeds manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#
FID Collector Assembly Kit (includes insulator)	G1531-60690	kit	21699
FID Collector Assembly Kit w/Siltek® Ignitor Castle	—	kit	21132

Replacement FID Parts for Agilent 6890/6850/7890 GCs

- Meets or exceeds manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#
A) FID Collector (includes insulators)	G1531-20690 G1531-20700 19231-20940	ea.	21139
B) FID Collector Nut and Washer	5181-3311	set	21136
C) FID Ignitor*	19231-60680	ea.	21001
D) FID Ignitor Castle	19231-20910	ea.	21137
Siltek® FID Ignitor Castle	—	ea.	21135

*Also fits OI Analytical 4410 detector (similar to OI part # 191833).

**Direct Replacement FID Flow Measuring Adaptor for Agilent 5890/6890/6850/7890 GCs**

- Makes setting flows easy.
- Meets or exceeds manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#
FID Flow Measuring Adaptor	19301-60660	ea.	21000

Rugged design!

**FID Gauge Pack for Agilent 5890 GCs**

Pressure regulators and gauges for air & hydrogen. The 1/8-inch bulkhead allows easy hookup to instrument. Rated for inlet pressures to 250psi (1724kPa), outlet pressures of 0 to 60psi (0-414kPa).

Description	qty.	cat.#
FID Gauge Pack for Agilent 5890 GCs	ea.	22071