

Supplies for Varian Instruments: Liner Seals, Septa, Septum Nuts



Inlet Liner Seals for Varian 1177 Injectors

- Meets or exceeds manufacturer's performance.

Description	Max. temp.	Similar to Varian part #	10-pk. cat.#	50-pk. cat.#
6.35mm ID Graphite O-rings for split liners	450°C	—	20296	20297
6.5mm ID Graphite O-rings for splitless liners	450°C	39-26119-40	20298	20299



Liner Seals for Varian 1078/1079

Description	Max. temp.	Similar to Varian part #	qty.	cat.#
5mm Graphite Liner Seals for Varian 1078/1079 GCs	450°C	392611919 392534201	10-pk.	22683

Now precision molded!

Septa for Varian GCs

- Precision molding assures consistent, accurate fit.
- Preconditioned and ready to use.
- Do not adhere to hot metal surfaces.
- Packaged in precleaned glass jars.

Septum Diameter	25-pk.	50-pk.	100-pk.
Thermolite® Septa (usable to 340°C inlet temp.)			
9mm	27132	27133	27134
9.5mm (9/16")	27135	27136	27137
10mm	27138	27139	27140
11mm (7/16")	27141	27142	27143

IceBlue™ Septa (usable to 250°C inlet temp.)			
9mm	—	27156	27157
9.5mm (9/16")	—	27158	27159
10mm	—	27160	27161
11mm (7/16")	—	27162	27163

BTO® Septa (usable to 400°C inlet temp.)			
9mm CenterGuide™	—	27104	27105
9.5mm (9/16")	—	27106	27107
10mm	—	27108	27109
11mm (7/16") CenterGuide™	—	27110	27111

HANDY septum size chart

Injector Type	Septum Diameter (mm)
Packed Column	9.5/10
Split/Splitless 1078/1079	10/11
1177	9
1075/1077	11

Stainless Steel Septum Nuts for Varian GCs

- Maintain a sharper knife edge longer than the original equipment aluminum design.
- Greatly minimize septum fragmentation and increase septum lifetime.

Description	Similar to Varian part #*	qty.	cat.#
26 Gauge Needle Guide Septum Nut	03-949666-01	ea.	21310
Standard Stainless Steel Septum Nut	03-949666-01	ea.	21304

*Similar to Varian part # 03-949666-01, except Restek nuts are not equipped with remote-start capabilities.

