

Graphite O-Rings for PerkinElmer Auto SYS™ XL PSS

| Description | Max. temp. | Similar to PE part # | qty. | cat.# |
|---|------------|----------------------|--------|-------|
| Graphite O-Rings for PerkinElmer Auto SYS™ XL PSS | 450°C | N610-1751 | 10-pk. | 21475 |
| Graphite O-Rings for PerkinElmer Auto SYS™ XL PSS | 450°C | N610-1751 | 25-pk. | 21476 |



Viton® O-Rings for PerkinElmer PSS

| Description | Max. temp. | Similar to PE part # | qty. | cat.# |
|------------------------------------|------------|----------------------|--------|-------|
| Viton® O-Rings for PerkinElmer PSS | 250°C | N6101747 | 10-pk. | 20366 |



Silicone O-Rings for PerkinElmer Auto SYS™ GCs

| Description | Max. temp. | Similar to PE part # | qty. | cat.# |
|--|------------|----------------------|--------|-------|
| Silicone O-Rings for PerkinElmer Auto SYS™ GCs | 250°C | N6101374 | 10-pk. | 20262 |



11mm Septa for PerkinElmer GCs

- Preconditioned and ready to use.
- Packaged in precleaned glass jars.

Now precision molded!



Thermolite® Septa

(usable to 340°C inlet temp.)

| qty. | cat.# |
|---------|-------|
| 25-pk. | 27141 |
| 50-pk. | 27142 |
| 100-pk. | 27143 |



IceBlue™ Septa

(usable to 250°C inlet temp.)

| qty. | cat.# |
|---------|-------|
| 50-pk. | 27162 |
| 100-pk. | 27163 |



BTO® Septa

(usable to 400°C inlet temp.)

| qty. | cat.# |
|---------|-------|
| 50-pk. | 27110 |
| 100-pk. | 27111 |

HANDY septum size chart

| PE Instrument | Septum Diameter (mm) |
|---------------|----------------------|
| Sigma series | 11 |
| 900,990 | 11 |
| 8000 series | 11 |
| Auto SYS™ | 11 |
| Auto SYS™ XL | 11 |

Septum Cap for PerkinElmer Auto SYS™ XL

- Made of clear anodized aluminum and high-quality stainless steel.
- Meets or exceeds manufacturer's performance.



| Description | Similar to PE part # | qty. | cat.# |
|---|----------------------|------|-------|
| Septum Cap for PerkinElmer Auto SYS™ XL | N6100153 | ea. | 22322 |

Injector Adaptor for PerkinElmer Auto SYS™ XL

- Made of high-quality stainless steel.
- Meet or exceed manufacturer's performance.
- Siltek® treated version available for increased inertness.

| Description | Similar to PE part # | qty. | cat.# |
|--|----------------------|------|-------|
| For use with PE style capillary nuts | | | |
| A) Injector Adaptor for PerkinElmer Auto SYS™ XL | N6100157 | ea. | 22318 |
| B) Siltek® Treated Injector Adaptor for PerkinElmer Auto SYS™ XL | — | ea. | 22320 |
| For use with 1/8" compression style nuts | | | |
| C) Injector Adaptor for PerkinElmer Auto SYS™ XL | — | ea. | 22319 |
| D) Siltek® Treated Injector Adaptor for PerkinElmer Auto SYS™ XL | — | ea. | 22321 |

