

Supplies for Thermo Scientific Instruments: Replacement Parts



A



B



C



D



Adapters for Capillary Columns on Thermo Scientific TRACE™ and Focus SSL

- Use same installation distance as manufacturer's adapters.
- Made of high-quality stainless steel.

Description	Similar to TS part #	qty.	cat.#
For use with standard 1/16" ferrules.			
A) Adapter for Capillary Column on Detector Base	—	ea.	24916
B) Adapter for Capillary Column on Split/Splitless Injector	—	ea.	24917
For use with M4 ferrules. new!			
C) Adapter for Capillary Column on Detector Base	347 25 436	2-pk.	24969
D) Adapter for Capillary Column on Split/Splitless Injector	347 05 451	ea.	24970

Graphite Ferrules—For M4 Fittings for Thermo Scientific TRACE™ and Focus SSL

Ferrule ID	Fits Column ID	Similar to TS part #	Graphite	
			2-pk.	10-pk.
0.3mm	0.10/0.15mm	—	22221	22222
0.4mm	0.18/0.28mm	290-13488	20280	20281
0.5mm	0.32mm	290-13487	20282	20283
0.8mm	0.45/0.53mm	290-13486	20284	20285

new!

Fixing Nut for Capillary Column for Split/Splitless Injector on Thermo Scientific TRACE™ and Focus SSL

- Made of high-quality stainless steel.

Description	Similar to TS part #	qty.	cat.#
Fixing Nut for Capillary Column for Split/Splitless Injector on Thermo Scientific TRACE™ and Focus SSL	350 32 423	5-pk.	24973

new!

Nut for Terminal Fitting for Thermo Scientific TRACE™ GCs



Description	Similar to TS part #	qty.	cat.#
Nut for Terminal Fitting for Thermo Scientific TRACE™ GCs	350 221 25	2-pk.	24896