

Tube & Pipe Fittings

also available

For Siltek®/Sulfinert® and Silcosteel®-CR-treated fittings, see pages 392-393.

Tube & Pipe Fittings

The Instrumentation Group of Parker Balston® Corporation designs and manufactures a top-quality line of components and systems for use in process instrumentation, semiconductor manufacturing, and analytical instrumentation. Parker's product quality and delivery have made them a world-wide leader—and this is the level of quality and service Restek wants to deliver to you.

Fitting Type	Parker #	Similar to Swagelok® #	Size	qty.	Brass cat.#	316 Grade Stainless Steel	
						qty.	cat.#
nut 	1 Nu 1	102-1	1/16"	20-pk.	21800	5-pk.	21900
	2 Nu 2	202-1	1/8"	40-pk.	21801	10-pk.	21901
	4 Nu 4	402-1	1/4"	40-pk.	21802	10-pk.	21902
front ferrule 	1 FF 1	103-1	1/16"	20-pk.	21803	10-pk.	21903
	2 FF 2	203-1	1/8"	40-pk.	21804	20-pk.	21904
	4 FF 4	403-1	1/4"	40-pk.	21805	20-pk.	21905
back ferrule 	1 BF 1	104-1	1/16"	20-pk.	21806	10-pk.	21906
	2 BF 2	204-1	1/8"	40-pk.	21807	20-pk.	21907
	4 BF 4	404-1	1/4"	40-pk.	21808	20-pk.	21908
nut & ferrule set 	—	—	1/16"	10-pk.	21809	2-pk.	21909
	—	—	1/8"	20-pk.	21810	5-pk.	21910
	—	—	1/4"	20-pk.	21811	5-pk.	21911
plug 	1 BLP 1	100-P	1/16"	5-pk.	21815	2-pk.	21915
	2 BLP 2	200-P	1/8"	10-pk.	21816	4-pk.	21916
	4 BLP 4	400-P	1/4"	10-pk.	21817	4-pk.	21917
union 	1 SC 1	100-6	1/16"	3-pk.	21818	ea.	21918
	2 SC 2	200-6	1/8"	5-pk.	21819	2-pk.	21919
	4 SC 4	400-6	1/4"	5-pk.	21820	2-pk.	21920
reducing union 	2 RU 1	200-6-1	1/8" to 1/16"	5-pk.	21823	ea.	21923
	4 RU 1	400-6-1	1/4" to 1/16"	5-pk.	21824	2-pk.	21924
	4 RU 2	400-6-2	1/4" to 1/8"	5-pk.	21825	2-pk.	21925

Parker's (A-Lok) two-piece ferrules and NPT fittings are ideal for installing new equipment, modifying existing instrumentation, or replacing worn connections.

Size	Parker #	Similar to Swagelok®	qty.	Brass cat.#	316 Grade Stainless Steel		Fitting Type
					qty.	cat.#	
1/16"	1 ET 1	100-3	2-pk.	21826	ea.	21926	 tee
1/8"	2 ET 2	200-3	2-pk.	21827	ea.	21927	
1/4"	4 ET 4	400-3	2-pk.	21828	ea.	21928	
1/8"	2 ECR 2	200-4	2-pk.	21829	ea.	21929	 cross
1/4"	4 ECR 4	400-4	2-pk.	21830	ea.	21930	
^A 1/8" tube to ^B 1/16"	2 TUR 1	100-R-2	5-pk.	21831	2-pk.	21931	 A B tube end reducer
1/4" tube to 1/16"	4 TUR 1	100-R-4	5-pk.	21832	2-pk.	21932	
1/8" tube to 1/4"	2 TUR 4	400-R-2	5-pk.	21833	2-pk.	21933	
1/4" tube to 1/8"	4 TUR 2	200-R-4	5-pk.	21834	2-pk.	21934	
1/8"	2 PC 2	201-PC	5-pk.	21835	2-pk.	21935	 port connector
1/4"	4 PC 4	401-PC	10-pk.	21836	2-pk.	21936	
1/8" tube to 1/4"	2 PC 4	401-PC-2	5-pk.	21837	2-pk.	21937	
^A 1/8" to ^B 1/8" NPT	2 MSC 2N	200-1-2	10-pk.	21841	2-pk.	21941	 A B male connector
1/4" to 1/4" NPT	4 MSC 4N	400-1-4	10-pk.	21842	2-pk.	21942	
1/16" to 1/8" NPT	1 MSC 2N	100-1-2	5-pk.	21843	2-pk.	21943	
1/8" to 1/4" NPT	2 MSC 4N	200-1-4	10-pk.	21844	2-pk.	21944	
1/4" to 1/8" NPT	4 MSC 2N	400-1-2	10-pk.	21845	2-pk.	21945	
^A 1/8" to ^B 1/8" NPT	2 FSC 2N	200-7-2	5-pk.	21846	2-pk.	21946	 A B female connector
1/4" to 1/4" NPT	4 FSC 4N	400-7-4	5-pk.	21847	2-pk.	21947	
1/4" to 1/8" NPT	4 FSC 2N	400-7-2	5-pk.	21848	2-pk.	21948	
1/8" male*	2A-Q4VN	QC4D-200	—	—	ea.	21957	 male & female quick couplings
1/8" male	2A-Q4P	QC4S-200	—	—	ea.	21958	
1/8" female*	2A-Q4CN	QC4B-200	—	—	ea.	21959	
1/4" male*	4A-Q4VN	QC4D-400	—	—	ea.	21960	
1/4" male	4A-Q4P	QC4S-400	—	—	ea.	21961	
1/4" female*	4A-Q4CN	QC4B-400	—	—	ea.	21962	

*Includes self-sealing shut-off valve.