

Manufacturer	Instrument/ Model #	11mm Crimp-Top (p. 264-265)	9mm Short-Cap, Screw-Thread (p. 269)	8mm Crimp-Top, 8 x 30mm (p. 267)	Headspace (p. 268)	8mm Screw-Thread (p. 270-271)	10mm Screw-Thread (p. 271)	WTSP 48 Style (p. 272-273)	WTSP 96 Style (p. 273)
A.I.M.	CPS -100, 200	X				X			
A.I.	42 vial tray	X							
A.I.	60 vial tray			X					
A.I.	Headspace				X				
Agilent	1042, 1050, 1080, 1082, 1084, 1090, 1100, 5890, 6850, 6890, 7670A, 7671A, 7672, 7673A/B, 7683, 7890, 8042	X	X						
Agilent	Headspace				X				
Alcott	738	X				X	X		
Alltech	570	X				X	X		
Alltech	580 (standard tray)	X							
Altex		X					X		
AMS	42	X							
Antek	736 Unisampler, 738	X				X	X		
ASC		X				X	X		
Beckman	501, 502, 507					X	X		
Bruker	LC 51							X	
Carlo Erba/Fisons	42 vial tray	X							
Carlo Erba/Fisons	60 vial tray, AS 105 tray			X					
Carlo Erba/Fisons	AS 800	X	X			X			
Carlo Erba/Fisons	8000					X			
Carlo Erba	HS 250, HS 500, HS 800				X				
Carnegie	CMA-250/200	X							
Chrompack	CP 9000 GC Series	X							
CTC	CTC A105S	X							
CTC	Headspace				X				
Cueni	Headspace				X				
Dani	ALS 39.80, 86.80	X	X			X			
Dani	HSS 39.50, HSS 86.50, SPT 37.50				X				
Delsi		X							
Dynatech	42 vial tray	X	X			X			
Dynatech	60 vial tray			X					
Dynatech	LC2000		X						
Fisons	42 vial tray	X							
Fisons	60 vial tray			X					
GBC		X							
Gerstel		X							
Gilson	231 XL, 232 XL, 233 XL, Aspec XL	X	X			X			
Gilson	Asted XL	X				X			
Gynkotec	Gina, others	X				X			
Hitachi	AS-2000, AS-6000	X	X						
Hitachi	S6551					X			
Hitachi	L-7200					X	X	X	
Hitachi	L-7250					X	X		
IBM		X				X			
ICI	other	X				X			
Infochroma		X				X			
Jasco	851-AS, AS-950	X				X			
Jasco	LC800 & 900 series					X			
Kipp		X				X			
Kontron	360	X				X			
Kontron	460					X			
LDC	Marathon, Promis	X	X						
LDC	other	X							
LDC	713			X					
LDC	all except 713					X			
L.E.A.P. Technologies	CTC A200S	X	X			X	X		
L.E.A.P. Technologies	CTC A105S	X							

Manufacturer	Instrument/ Model #	11mm Crimp-Top (p. 264-265)	9mm Short-Cap, Screw-Thread (p. 269)	8mm Crimp-Top, 8 x 30mm (p. 267)	Headspace (p. 268)	8mm Screw-Thread (p. 270-271)	10mm Screw-Thread (p. 271)	WTSP 48 Style (p. 272-273)	WTSP 96 Style (p. 273)
L.E.A.P. Technologies	HS 500				X				
Magnus Scientific		X				X			
PerkinElmer	AI-1	X	X	X					
PerkinElmer	ISS 100, LC 600 42 vial tray	X	X						
PerkinElmer	ISS 200	X	X		X				
PerkinElmer	Integral 4000	X			X				
PerkinElmer	Autosystem	X			X	X			
PerkinElmer	420/B, 4900	X				X			
PerkinElmer	F40, F45, HS 6, HS 40				X				
PerkinElmer	LC 600 (60), AS100/B, 300, 8300, AS2000, LC 30, 2000, B 2000			X					
Pharmacia LKB	2157-010	X				X			
Precision Sampling	GC111, GC 311, LC 241-60	X		X		X			
Phillips	4247, 4710	X				X			
Phillips	LC-XP	X				X		X	
Polymer-Labs	GPC 110/210	X							
Pye	LCXP	X							
S.G.E.	M280D	X				X			
Sedere		X				X			
Shimadzu	SIL-10A, 10Ai, 10Axl	X	X			X	X		
Shimadzu	AOC-8B/9, SIL-6A, SIL-6B, -9A, -8A					X			
Shimadzu	AOC-14/1400, AOC-17					X		X	
Shimadzu	SIL-2AS							X	
Shimadzu	HSS-2B				X				
Siemens	AS 32, AS 200	X							
Sigma	2000			X					
Spark Holland	SPH 125	X	X						
Spark Holland	Marathon, Promis	X	X			X			
Spark Holland	Triathlon, Midas	X				X			
Spectra-Physics	8875, 8880	X	X			X			
Talbot		X				X			
Tekmar	7000, 7000/7050				X				
Thermo Scientific	TRACE GC 2000, AS2000	X	X	X					
TOA	ICA5450					X			
TosoHaas	TSK-6080, AS-8010, AS-8020	X				X		X	
Tosca 1		X							
Tracor	770, 771, 772	X							
TSP	8875, 8880, AS 100/1000, AS 300/3000	X				X			
Unicam	4710	X							
Unicam	4700LC, 4247	X				X			
Unicam	LC-XP	X				X		X	
United Technologies		X				X			
Varian/Rainin	Dynamax AI-1A, AI-200, Dynamax AI-3	X							
Varian/Rainin	DynamaxAI-3AI-1A, AI-2W					X			
Varian	8100/8200	X	X			X	X		
Varian	Marathon	X	X			X			
Varian	9100/9090/9095	X				X			
Varian	8035, 8000 series, 8055, 8085, 3800 GC					X			
Varian	Genesis				X				
Waters	Alliance 2690	X	X			X	X		
Waters	710, 717+							X	X
Waters	48-pos. M700							X	
Waters	96-pos. M700								X