

0.8mL, 8mm Crimp-Top Vials, VOA Sampling Vials

0.8mL, 8mm Crimp-Top Vials, 8 x 30mm

0.8mL Crimp-Top Vials (8 x 30mm, 8mm Crimp)

Description	100-pk.	1000-pk.
0.8mL Clear Flat Bottom Vials	24695	24696

0.8mL, 8mm Aluminum Seals w/Septa

Description	100-pk.	1000-pk.
Silver Seals w/Clear PTFE/Natural Rubber Septa	21757	21758



0.8mL, 8mm Crimper and Decapper

Description	qty.	cat.#
8mm Crimper	ea.	21735
8mm Decapper	ea.	21736



instrument reference

For 0.8mL, 8mm Crimp-Top Vials, 8 x 30mm

The following instruments are compatible with 0.8mL, 8mm crimp-top vials. For further information or questions, please contact your Restek sales representative or technical service.

Manufacturer	Instrument/Model #
A.I.	60 vial tray
Carlo Erba/Fisons	60 vial tray, AS 105 tray
Dynatech	60 vial tray
Fisons	60 vial tray
LDC	713
PerkinElmer	LC 600 (60), AS100/B, 300, 8300, AI-1, AS2000, LC 30, 2000, B 2000
Precision Sampling	GC 111, GC 311, LC241-60
Sigma	2000
Thermo Scientific	TRACE GC 2000, AS2000

also available

See page 274 for electronic crimpers and decappers.

Pre-Cleaned Volatile Organic Analyte (VOA) Sampling Vials

- Container, liner and closure cleaned and assembled.
- Clear or amber.
- Open top caps.
- Teflon® faced 0.125" silicone septa.
- Each case lot numbered.

Description	Screw-Thread	qty.	cat.#
	Size		
20mL Clear Pre-Cleaned VOA Vials	24mm-400	72-pk.	21798
20mL Amber Pre-Cleaned VOA Vials	24mm-400	72-pk.	21799
40mL Clear Pre-Cleaned VOA Vials	24mm-400	72-pk.	21796
40mL Amber Pre-Cleaned VOA Vials	24mm-400	72-pk.	21797
60mL Collection Vials, Clear Glass, for ASE® 200	24mm-400	72-pk.	26121
60mL Collection Vials, Amber Glass, for ASE® 200	24mm-400	72-pk.	26122
Replacement Septa for 20, 40 & 60mL Vials and 250mL Bottles (24mm Teflon®-lined 0.125" silicone septa)		100-pk.	24694

