



Headspace Autosampler Vials

| Description | 100-pk. | 1000-pk. |
|------------------------------------|---------|----------|
| a) 6mL Clear Vial | 21166 | 21167 |
| b) 10mL Clear Vial, Flat Bottom | 24683 | 24684 |
| c) 10mL Clear Vial, Rounded Bottom | 21164 | 21165 |
| d) 20mL Clear Vial, Flat Bottom | 24685 | 24686 |
| e) 20mL Clear Vial, Rounded Bottom | 21162 | 21163 |
| f) 27mL Clear Vial | 21160 | 21161 |

a) - f) vial to instrument compatibility are designated in instrument reference chart below.

instrument reference

For Headspace Vials

The following instruments are compatible with the headspace vials listed above (as letter designated above). For further information or questions, please contact your Restek sales representative or technical service.

| Manufacturer | Instrument/Model # | Manufacturer | Instrument/Model # |
|-----------------------|--|-----------------------|--|
| a) Carlo Erba | HS 500, HS 800 | d) A.I. | Headspace |
| Agilent | Headspace | Carlo Erba | HS 500, HS 800 |
| PerkinElmer | ISS 200 | CTC | Headspace |
| | | Cueni | Headspace |
| b) A.I. | Headspace | Dani | HSS 39.50, HSS 86.50, SPT 37.50 |
| Carlo Erba | HS500, HS 800 | Agilent | Headspace |
| Cueni | Headspace | | |
| Dani | HSS 39.50, HSS 86.50 | e) A.I. | Headspace |
| Agilent | Headspace | Cueni | Headspace |
| | | Dani | HSS 39.50, HSS 86.50, SPT 37.50 |
| c) A.I. | Headspace | L.E.A.P. Technologies | HS 500 |
| Cueni | Headspace | PerkinElmer | Autosystem GC, F40, F45, HS 6, Integral 4000, HS 40 |
| Dani | HSS 39.50, HSS 86.50 | Tekmar | 7000, 7000/7050 |
| L.E.A.P. Technologies | HS 500 | Varian | Genesis |
| PerkinElmer | Autosystem GC, F40, F45, HS 6, Integral 4000, HS 40 | | |
| Tekmar | 7000, 7000/7050 | f) Shimadzu | HSS-2B |
| Varian | Genesis | | |

20mm Aluminum Seals w/Septa, Assembled

| Description | 100-pk. | 1000-pk. |
|---|---------|----------|
| Silver Seal w/ PTFE/Gray Butyl Rubber | 21761 | 21762 |
| Silver Seal w/ PTFE/Silicone | 21763 | 21764 |
| Pressure Release Silver Seal w/ PTFE/Gray Butyl Rubber Septum | 21765 | 21766 |
| Pressure Release Silver Seal w/ PTFE/Silicone Septum | 21767 | 21768 |
| 20mm Magnetic Crimp Cap w/Teflon®/Silicone Septum (tan/white) with 8mm Hole | 22831 | 22832 |
| 20mm Magnetic Crimp Cap w/Teflon®/Silicone Septum (tan/white) with 5mm Hole | 22833 | 22834 |
| 20mm Magnetic Crimp Cap w/PTFE/Butyl Septum with 8mm Hole | 22835 | 22836 |
| 20mm Magnetic Crimp Cap w/PTFE/Butyl Septum with 5mm Hole | 22837 | 22838 |

5mm hole is compatible with the following systems: Carlo Erba HS500/HS800, CTC 500, Fisons HS500/HS800 and Varian/Chrompack 9020/25. 8mm hole is compatible with the Combi Pal.

also available

See page 274 for 20mm electronic crimper and decapper.

tech tip

Compatibility concerns?
Refer to the Septum Selection Guide on page 263.

