


Autosampler Syringes for Agilent GCs

Replacement Teflon®-Tipped Plunger Assembly for Agilent 7673, 7683, and 6850 Autosampler Syringes




SGE

Syringe Volume	Needle Term.	SGE cat.#	qty.	Restek cat.#
10 μ L	F	031808	2-pk.	21218
10 μ L	R	031809	2-pk.	21284

Needle Gauge for Agilent 7673, 7683, 6850 Syringes


Single Gauge Needle

- The most popular gauge for Agilent 7673.
- Best for Merlin Microseal septum and standard septum-equipped GCs.
- Packed column injection ports.
- Split/splitless injection ports.




Single Gauge Needle

- On-column injection ports.
- Split/splitless injection ports.



23s-26s—Dual Gauge (tapered)

- Durability of a 23s gauge needle.
- Ability of a 26s gauge needle to perform split/splitless and on-column injections.



CROSS-reference

For Agilent Syringes

Listed below are some Hamilton and SGE syringes, cross-referenced to Agilent part numbers for your convenience. For further information or questions, contact your Restek sales representative or technical service.

Similar to Agilent Part #	Restek Part #	Description	Page
5181-1267	21212, 24596	10 μ L Standard Syringe, Fixed Needle, 23-26 Gauge	278
5181-1273	21210, 24593	5 μ L Standard Syringe, Fixed Needle, 23-26 Gauge	278
5181-3319	20166	Needle for 10 μ L Standard Syringe, 23-26 Gauge	279
5181-3321	21213, 24598	10 μ L Standard Syringe, Removable Needle, 23-26 Gauge	278
5181-3356	21218	Replacement Teflon®-Tipped Plunger for 10 μ L Fixed Needle Syringe	280
5181-3358	24919	Teflon®-Tipped Plunger for 10 μ L Standard Syringe, Removable Needle	290
5181-3360	21215, 24600	10 μ L Standard Syringe, Fixed Needle, 23-26 Gauge (6-pk.)	278
5181-3365	21218	Teflon®-Tipped Plunger for 10 μ L Standard Syringe, Fixed Needle	280
5181-8809	24894, 24789	10 μ L Gas-tight Syringe, Fixed Needle, 23 Gauge	279
5181-8811	20164, 24825	Needle for 10 μ L Gas-tight Syringe, 23 Gauge	279
5181-8813	21224, 24797	10 μ L Gas-tight Syringe, Removable Needle, 23 Gauge	279
5182-0830	24821	Needle for 0.5 μ L Standard Syringe, 23 Gauge	279
5182-0832	20166	Needle for 5-10 μ L Autosampler Syringe, 23-26 Gauge	279
5182-0834	20172, 24793	5 μ L Standard Syringe, Removable Needle, 23 Gauge	278
5182-0835	21223, 21211	5 μ L Standard Syringe, Removable Needle, 23-26 Gauge	278
5182-0875	20170, 24783	5 μ L Standard Syringe, Fixed Needle, 23 Gauge (6-pk.)	278
9301-0713	20167, 24785	10 μ L Standard Syringe, Fixed Needle, 23 Gauge	278
9301-0714	24595, 24784	10 μ L Standard Syringe, Fixed Needle, 26 Gauge	278
9301-0725	20169, 24787	10 μ L Standard Syringe, Fixed Needle, 23 Gauge (6-pk.)	278
9301-0891	24592, 24780	5 μ L Standard Syringe, Fixed Needle, 26 Gauge	278
9301-0892	20168, 24781	5 μ L Standard Syringe, Fixed Needle, 23 Gauge	278