

## Autosampler Syringes for PerkinElmer Autosystem™ and Clarus 500 GCs

- Hamilton syringes are designed and tested to meet critical autosampler performance.
- Needle point styles are developed to withstand multiple and rapid injections through a septum.



Hamilton 5µL Autosampler Syringe

### Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton			Restek cat.#
					Model	cat.#	qty.	
5µL	ASN	23	2.756"	3	75ASN/PE	88035	ea.	24953
5µL	ASN	26	2.756"	3	75ASN/PE	88040	ea.	24952



SGE 5µL Autosampler Syringe

### SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE			Restek cat.#
					Model	cat.#	qty.	
0.5µL	R	23	70mm	Cone	0.5BR-PE-0.63	000478	ea.	24811
0.5µL	R	26	70mm	Cone	0.5BR-PE-0.47	000475	ea.	24810
5µL	F	23	70mm	Cone	5F-PE-0.63	001954	ea.	24813
5µL*	F	23	70mm	Cone	5F-PE-GT-0.63	001957	ea.	24815
5µL	F	26	70mm	Cone	5F-PE-0.47	001953	ea.	24812
5µL*	F	26	70mm	Cone	5F-PE-GT-0.47	001955	ea.	24814
50µL	F	23	70mm	Cone	50F-PE-0.63	004670	ea.	24816

\*Gas-tight.

### Guide to Needle Termination Codes

#### Hamilton:

(N) Cemented Needle  
 (RN) Removable Needle  
 (ASN) Autosampler Cemented Needle  
 (ASRN) Autosampler Removable Needle  
 (TLL) Teflon® Luer Lock  
 (KH) Knurled Hub  
 (LT) Luer Tip  
 (LTN) Luer Tip Cemented Needle


#### SGE:

(F) Fixed Needle  
 (R) Removable Needle  
 (LL) Luer Lock  
 (LT) Luer Tip

## Replacement Needle and Plunger Kit for PerkinElmer Autosampler Syringes

### SGE Needle and Plunger Kit

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	qty.	Restek cat.#
0.5µL	23	70mm	Cone	033765	kit	24818



Needle and Plunger Kit

also available

### Looking for septa?

Restek carries a full line of septa for every type of application and instrument. See pages 134-135.

