

## Autosampler Syringes for CTC/LEAP GCs

- Hamilton syringes are designed and tested to meet critical autosampler performance.
- Needle point styles are developed to withstand multiple and rapid injections through a septum.



Hamilton 10µL Autosampler Syringe

### Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton			Restek cat.#
					Model	cat.#	qty.	
1.2µL	N	26	2"/51mm	AS	7701.2N	203185	ea.	22755
5µL	N	26s	2"/51mm	AS	75N	203189	ea.	22756
10µL	N	23s	2"/51mm	AS	701N	203361	ea.	22757
10µL	N	23s-26s	2"/51mm	AS	701N	203362	ea.	22758
10µL	N	23s	2"/51mm	2	701N	203363	ea.	22759
10µL	N	26s	2"/51mm	2	701N Slim Line	203072	ea.	22760
10µL	N	26s	2"/51mm	AS	701N	203205	ea.	22761
10µL	N	26s	2"/51mm	2	701N	203269	6-pk.	22762
25µL*	N	26s	2"/51mm	AS	1702N	203043	ea.	22763
25µL*	N	26s	2"/51mm	AS	1702N Slim Line**	203074	ea.	22764
100µL*	N	26s	2"/51mm	AS	1710N Slim Line**	203076	ea.	22765

\*Gas-tight syringe.

\*\*Barrel OD = 6.7mm; all other 25µL and 100µL syringes have a 7.9mm barrel OD.



SGE 10µL Autosampler Syringe

### SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE			Restek cat.#
					Model	cat.#	qty.	
0.5µL	R	23	50mm	Cone	0.5BNR-C/F-0.63	000492	ea.	24928
0.5µL	R	26	50mm	Cone	0.5BNR-C/F-0.47	000490	ea.	24926
5µL	F	23	50mm	Cone	5F-C/F-5/0.63	001981	ea.	24921
5µL	F	26	50mm	Cone	5F-C/F-5/0.47	001982	ea.	24920
10µL	F	23	50mm	Cone	10F-C/F-0.63	002981	ea.	24923
10µL	R	23	50mm	Cone	10R-C/F-0.63	002984	ea.	24932
10µL	F	26	50mm	Cone	10F-C/F-5/0.47C	002980	ea.	24922
10µL	F	26	50mm	Cone	SK10F-C/F-5/0.47C	002986	6-pk.	24925
10µL	R	26	50mm	Cone	10R-C/F-5/0-0.47C	002982	ea.	24930
10µL*	F	23	50mm	Cone	10F-C/T-GT-5/0.63C	002987	ea.	22749
25µL	F	26	50mm	Cone	25F-C/T-5/0.47C	003980	ea.	22750
25µL*	F	23	50mm	Cone	25F-C/T-GT-5/0.63C	003987	ea.	22751
25µL*	R	26	50mm	S-hole	25R-C/T-GT-5/0.47H	003988	ea.	22752
100µL*	F	23	50mm	Cone	100F-C/T-GT-5/0.63C	005335	ea.	22753
100µL*	R	26	50mm	Cone	100R-C/T-GT-5/0.47C	005333	ea.	22754

\*Gas-tight syringe.

also available

Replacement plunger assemblies for these syringes are available on request.

## Replacement Needles for CTC/LEAP Autosampler Syringes



### Needles for SGE Syringes

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	qty.	Restek cat.#
0.5µL	23	50mm	Cone	033772	ea.	24929
10µL	23	50mm	Cone	037787	2-pk.	24933
10µL	26	50mm	Cone	037785	2-pk.	24931