

Syringes for Waters U6K Valves

- Reinforced plungers.
- Removable needles.
- Replaceable plunger/barrel assemblies.
- Barrel stop prevents plunger blowout.



Hamilton Syringes

Hamilton 50µL Syringe for Waters U6K Valves

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		Restek
					Model	cat.#	
10µL	RN	25s	1.97"	3	801	84815	ea. 21255
25µL	RN	25s	1.97"	3	802	84816	ea. 21256
50µL	RN	25s	1.97"	3	805	84817	ea. 21257
100µL	RN	25s	1.97"	3	810	84818	ea. 21258
250µL	RN	25s	1.97"	3	825	84819	ea. 21259



SGE Syringes

SGE 50µL Syringe for Waters U6K Valves

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		Restek
					Model	cat.#	
25µL	R	25	50mm	3	25R-WLC-RA8	003375	ea. 24873
50µL	R	25	50mm	3	50R-WLC-RA8	004375	ea. 24874
100µL	R	25	50mm	3	100R-WLC-RA8	005375	ea. 24875
250µL	R	25	50mm	3	250R-WLC-RA8	006375	ea. 24876
500µL	R	25	50mm	3	500R-WLC-RA8	007375	ea. 24877

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) Teflon® Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip

Replacement Needles for Syringes for Waters U6K Valves

Needles for Hamilton Syringes



Syringe Volume	Needle Gauge	Needle Length	Point Style	Hamilton cat.#	Restek cat.#
10-100µL	25s	1.97"	3	8647-01	6-pk. 21270
250µL-10mL	25s	1.97"	3	8648-01	6-pk. 21271