

Syringes for CTC/LEAP LC Autosamplers



Hamilton 100µL LC Autosampler Syringe

Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		qty.	Restek cat.#
					Model	cat.#		
10µL	N	22s	2"/51mm	3	701N	203073	ea.	22743
25µL	N	22s	2"/51mm	3	1702N	203075	ea.	22744
100µL	N	22s	2"/51mm	3	1710N Slim Line*	203077	ea.	22745
100µL	N	22	2"/51mm	3	1710N	203235	ea.	22746
250µL	N	22	2"/51mm	3	1725N	203079	ea.	22747
500µL	N	22	2"/51mm	3	1750N	203349	ea.	22748

*Barrel OD = 6.7mm; all other 25µL and 100µL syringes have a 7.9mm barrel OD.

also **available**

Replacement plunger assemblies for these syringes are available on request.



SGE 25µL LC Autosampler Syringe

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		qty.	Restek cat.#
					Model	cat.#		
10µL	F	22s	2"/51mm	LC	10F-CTC-LC	002710	ea.	22737
25µL	F	22s	2"/51mm	LC	25F-C/TLC	003984	ea.	22738
25µL*	R	22s	2"/51mm	LC	25R-C/T-GT-LC	003985	ea.	22739
50µL*	F	22s	2"/51mm	LC	50F-C/T-GT-LC	004810	ea.	22740
100µL*	R	22s	2"/51mm	LC	100R-C/T-GT-LC	005330	ea.	22741
500µL*	F	22s	2"/51mm	LC	500F-CTC-GT-LC(0.41)	007720	ea.	22742

*Gas-tight syringe.

did you **know**?

Restek sells nitrogen generators from Parker for LC/MS systems. See [page 245](#).



Guide to Needle Termination Codes

Hamilton:

(N) Cemented Needle
(RN) Removable Needle
(ASN) Autosampler Cemented Needle
(ASRN) Autosampler Removable Needle
(TLL) Teflon® Luer Lock
(KH) Knurled Hub
(LT) Luer Tip
(LTN) Luer Tip Cemented Needle

SGE:

(F) Fixed Needle
(R) Removable Needle
(LL) Luer Lock
(LT) Luer Tip