


Bottle Tops



Hub-Cap (assembly of the bottle cap and plug)

Hub-Cap 4 Liter Bottle Tops

Most bottles use a GL45 cap. New Hub-Cap bottle tops are a great way to neatly keep your mobile phase lines where they belong. Use instead of parafilm, aluminum foil, or tape on your mobile phase reservoirs.



Description	qty.	cat.#
Hub-Cap (assembly of the bottle cap and plug)	kit	26541
Hub-Cap Multi-pack	3-pk.	26542



Hub-Cap Adapter

Hub-Cap Adapter and Opti-Cap™

Hub-Cap Adapters

Allow the use of the Opti-Cap™ with 4-liter solvent bottles.

Description	qty.	cat.#
Hub-Cap Adapter	ea.	26538
Hub-Cap Adapter Multi-pack	3-pk.	26539
Hub-Cap Adapter and Opti-Cap™	kit	26540

new!

Hub-Cap Filter Kit

Transfer and filter mobile phase in a single step!

Kit includes: bottle adapter, bottle adapter nut, filter inlet cap, grid support, vacuum hose barb, tube compression fitting, 47mm grid, 47mm .22µm filter membrane, 47mm .45µm filter membrane, 1/4" OD x 1/8" ID ultra chemical resistant, Teflon® FEP lined Tygon® tubing (3'), 6" x 6" box with shrink wrap insert



Hub-Cap Filter Kit



Unscrew and lift off top.



Place membrane filter on top of grid.



Reattach top.



Connect vacuum line to side port.

Description	qty.	cat.#
Hub-Cap Filter Kit	kit	26395
Replacement Membrane Filters	qty.	cat.#
Polypropylene Membrane Filters, 47mm, 0.45µm	100-pk.	26396
Polypropylene Membrane Filters, 47mm, 0.22µm	100-pk.	26397
Nylon Membrane Filters, 47mm, 0.45µm	100-pk.	26398
Nylon Membrane Filters, 47mm, 0.22µm	100-pk.	26399



Polypropylene Membrane Filters