



**Restek Pack in a Box Kit**

Restek's Pack in a Box Kit is a complete column-packing system. Everything you need to pack your own analytical HPLC columns is included: pump, reservoir, column hardware, packing material, and system control software (computer not included). Detailed packing instructions and a DVD video also are included. The packing pump is electric, so no pressurized air is required. AC input requirement: 85-265 VAC.

Description	qty.	cat.#
Pack in a Box Kit	kit	26408

**HPLC Organizers\***



**Bill Grove**  
Instrument Innovations  
2+ years of service!



LCLocker™ HPLC Organizer



Deluxe BenchBooster™ Organizer



Mini pHPerch™ Storage Unit



TopLoader™ BalanceBank™ Storage Unit



HPLC 30-Column Storage Cabinet



Book Holders



Open Supply Bins, 13-bin unit



Glove Box Dispensers

Description	dimensions	qty.	cat.#
LCLocker™	24 x 12 x 6"	ea.	25149
BenchBooster™	24 x 7 x 12"	ea.	25150
Mini pHPerch™	13 x 12 x 6"	ea.	25147
TopLoader™	12 x 12 x 7"	ea.	25148
30 Column Cabinet	17 <sup>7</sup> / <sub>8</sub> x 15 x 2 <sup>7</sup> / <sub>8</sub> "	ea.	25159
Book Holder, Small	0.75" ID	ea.	25151
Book Holder, Large	1.5" ID	ea.	25152
Open Supply Bin, 4-bin Unit	24 x 12 x 10"	ea.	25153
Open Supply Bin, 5-bin Unit	12 x 16 x 10.5"	ea.	25154
Open Supply Bin, 13-bin Unit	12 x 7.5 x 12"	ea.	25155
Glove Box Dispenser, Single	5 <sup>7</sup> / <sub>8</sub> x 11 <sup>1</sup> / <sub>2</sub> x 4 <sup>1</sup> / <sub>2</sub> "	ea.	25156
Glove Box Dispenser, Double	10 <sup>3</sup> / <sub>8</sub> x 11 <sup>1</sup> / <sub>2</sub> x 4 <sup>1</sup> / <sub>2</sub> "	ea.	25157
Glove Box Dispenser, Triple	15 <sup>3</sup> / <sub>8</sub> x 11 <sup>1</sup> / <sub>2</sub> x 4 <sup>1</sup> / <sub>2</sub> "	ea.	25158

\*Please note: Accessories and supplies in photographs are not included.