



- Choose stainless steel or Siltek® deactivation for improved inertness.
- Caps include frit, PEEK™ seal, and PTFE O-ring.



please note

Resprep™ designed parts are not interchangeable with Dionex designed parts.

Replacement Parts for ASE® 200 Systems, Restek Enhanced Design

Resprep™ Extraction Cell Parts for ASE® 200 Systems, Restek Enhanced Design

- Simpler design with fewer parts.
- Faster installation.
- Easier cleaning.
- More cost effective than original equipment design.

To simplify replacing the frit and PEEK™ washer in the ASE® 200 extraction cell cap, we developed an enhanced cell design that eliminates the snap ring. The frit presses snugly into a groove in a PEEK™ washer design that fits into a recessed cavity in the enhanced extraction cell.

These Restek replacement parts will keep your extraction system running smoothly. All parts are economically priced to save you money and are designed to meet or exceed the performance of the original manufacturer's parts.

Description	Stainless Steel		Siltek® Treated	
	qty.	cat.#	qty.	cat.#
Extraction Cell Kit, Resprep™ for ASE® 200, 1mL	kit	25980	kit	25981
Extraction Cell Kit, Resprep™ for ASE® 200, 5mL	kit	25982	kit	25983
Extraction Cell Kit, Resprep™ for ASE® 200, 11mL	kit	25984	kit	25985
Extraction Cell Kit, Resprep™ for ASE® 200, 22mL new!	kit	25976	kit	25977
Extraction Cell Kit, Resprep™ for ASE® 200, 33mL	kit	25986	kit	25987
Extraction Cell Body, Resprep™ for ASE® 200, 1mL	ea.	25960	ea.	25961
Extraction Cell Body, Resprep™ for ASE® 200, 5mL	ea.	25962	ea.	25963
Extraction Cell Body, Resprep™ for ASE® 200, 11mL	ea.	25964	ea.	25965
Extraction Cell Body, Resprep™ for ASE® 200, 22mL new!	ea.	25978	ea.	25979
Extraction Cell Body, Resprep™ for ASE® 200, 33mL	ea.	25966	ea.	25967
Extraction Cell Caps, Resprep™ for ASE® 200	2-pk.	25968	2-pk.	25969
PEEK™ Seal/Frit Assembly, Resprep™ for ASE® 200	2-pk.	25970	2-pk.	25971
Frit, Resprep™ for ASE® 200	12-pk.	25972	12-pk.	25973

Description	qty.	cat.#
PEEK™ Seal, Resprep™ for ASE® 200	12-pk.	25974
PEEK™ Seal, Resprep™ for ASE® 200	48-pk.	25975
Funnel, Resprep™ for ASE® 200 new!	ea.	24099
PTFE O-Rings for ASE® 200 & ASE® 300 Caps	100-pk.	26187
Viton® O-Rings for ASE® 200 & ASE® 300 Caps	50-pk.	26188

20mm Filters for ASE® 200 Extraction Cells

- Cellulose or glass fiber construction.
- Cellulose filters available in economical 1000-packs.



Description	Similar to Dionex part #	qty.	cat.#
Cellulose Filters for ASE® 200	049458	100-pk.	26118
Cellulose Filters for ASE® 200	049458	1000-pk.	26190
Glass Fiber Filters for ASE® 200	047017	100-pk.	26119

60mL Sample Collection Vials for ASE® 200 Systems

- Cleaned, assembled, and ready to use.
- Clear or amber glass.
- Caps and Teflon®-lined septa included.



Description	Similar to Dionex part #	qty.	cat.#
60mL Collection Vials, Clear Glass, for ASE® 200	048784	72-pk.	26121
60mL Collection Vials, Amber Glass, for ASE® 200	048781	72-pk.	26122
Replacement Septa for 20, 40 & 60mL Vials and 250mL Bottles (24mm Teflon®-lined 0.125" silicone septa)	—	100-pk.	24694