

Sulfinert® Treated Swagelok® Sample Cylinders

- Stable storage of samples containing ppb levels of sulfur compounds.
- Manufactured by Swagelok®; D.O.T. rated to 1,800psi (12,411kPa) at room temperature.
- 304 grade stainless steel with 1/4" female NPT threads on both ends.

Description	Size	qty.	cat.#
Sulfinert® Sample Cylinder	75cc	ea.	24130
Sulfinert® Sample Cylinder	150cc	ea.	24131
Sulfinert® Sample Cylinder	300cc	ea.	24132
Sulfinert® Sample Cylinder	500cc	ea.	24133
Sulfinert® Sample Cylinder	1000cc	ea.	24134
Sulfinert® Sample Cylinder	2250cc	ea.	21394



Sulfinert® Treated Alta-Robbins Sample Cylinder Valves

- All wetted parts are Sulfinert® treated for inertness.
- Compatible with Sulfinert® treated Swagelok® sample cylinders.
- Large, durable, Kel-F® seat ensures leak-free operation; temperature range: -40°C to 120°C.

Description	qty.	cat.#
1/4" NPT Exit	ea.	21400
1/4" Compression Exit	ea.	21401
1/4" NPT with Dip Tube*	ea.	21402
1/4" NPT with 2850psi Rupture Disc	ea.	21403
1/4" NPT Male Inlet x 1/4" Female Outlet with 2850psi Rupture Disc	ea.	21404

*To order catalog #21402 (Sulfinert Alta-Robbins Sample Cylinder Valve, 1/4" NPT with Dip Tube), please call Customer Service at 800-356-1688, ext. 3, or contact your Restek representative. Specify dip tube length or % outage when ordering (maximum length = 5.25" / 13.3cm). Note: End of part will not be treated after cutting tube to length.



Sulfinert® Treated Rupture Disc Tee

2850psig (19,650kPa) rating; 1/4" NPT connections.

Description	qty.	cat.#
Sulfinert® Treated Rupture Disc Tee (1/4" NPT connections)	ea.	21396
Replacement Rupture Disc (not Sulfinert® treated)	ea.	24298



Sulfinert® Treated Gas Sampling Valves and Sample Loops

- Ideal for samples containing low concentrations of sulfur or other active compounds.
- Sample loop sizes from 5µL to 5cc.

Sulfinert® Treated Gas Sampling Valves (1/16" fittings, 0.40mm port diameter; "W Type" valve)

Description	qty.	cat.#
Sulfinert® Gas Sampling Valve; 4-Port	ea.	20584
Sulfinert® Gas Sampling Valve; 6-Port	ea.	20585
Sulfinert® Gas Sampling Valve; 10-Port	ea.	20586



Replacement Rotors (Not Coated)

Description	qty.	cat.#
Replacement Rotor for 4-Port Sulfinert® Gas Sampling Valve	ea.	20587
Replacement Rotor for 6-Port Sulfinert® Gas Sampling Valve	ea.	20588
Replacement Rotor for 10-Port Sulfinert® Gas Sampling Valve	ea.	20589

Sulfinert® Treated Gas Sample Loops (1/16" fittings, for "W Type" valves)

Description	Size	qty.	cat.#
Sulfinert® Sample Loops	5µL	ea.	22840
Sulfinert® Sample Loops	10µL	ea.	22841
Sulfinert® Sample Loops	20µL	ea.	22842
Sulfinert® Sample Loops	25µL	ea.	22843
Sulfinert® Sample Loops	50µL	ea.	22844
Sulfinert® Sample Loops	100µL	ea.	22845
Sulfinert® Sample Loops	250µL	ea.	22846
Sulfinert® Sample Loops	500µL	ea.	22847
Sulfinert® Sample Loops	1mL	ea.	22848
Sulfinert® Sample Loops	2mL	ea.	22849
Sulfinert® Sample Loops	5mL	ea.	22850

