

Treated Straight, 6-Foot Length Tubing

Silcosteel®-CR Treated Straight Seamless 316L Grade Stainless Steel Tubing

- Individual 6-foot pieces.

In response to customer requests, we offer 6-foot straight lengths of 1/8-, 1/4-, and 3/8-inch Silcosteel®-CR treated tubing. This tubing can be cut to your exact requirements using a standard tubing cutter.

6-foot Length ($\pm 1/2"$)

OD	ID	qty.	cat.#
1/8" (3.18mm)**	0.055" (1.40mm)	ea.	22898
1/4" (6.35mm)**	0.180" (4.57mm)	ea.	22899
3/8" (9.52mm)***	0.277" (7.04mm)	ea.	22900

ordering note

Required length in meters
x 3.2808 = length in feet.

Siltek®/Sulfinert® Treated Straight Seamless 316L Grade Stainless Steel Tubing

- Individual 6-foot pieces.

In response to customer requests, we offer 6-foot straight lengths of 1/8-, 1/4-, and 3/8-inch Siltek®/Sulfinert® treated tubing. This tubing can be cut to your exact requirements using a standard tubing cutter.

6-foot Length ($\pm 1/2"$)

OD	ID	qty.	cat.#
1/8" (3.18mm)**	0.055" (1.40mm)	ea.	22901
1/4" (6.35mm)**	0.180" (4.57mm)	ea.	22902
3/8" (9.52mm)***	0.277" (7.04mm)	ea.	22903

**0.035" wall thickness

***0.049" wall thickness



frequently asked question

Which treated tubing should I use?

Restek Performance Coatings offers several grades of tubing with varied treatments. This chart will help you determine the tubing best suited to your application with respect to pressure, sensitivity of your analysis, and acidic environment exposure. For more frequently asked questions on treated tubing, see [page 391](#) or contact the Restek Coatings experts.

Restek Tubing Reference Guide

