

Neat & Single Analyte Solutions



Top: Matt Hepfer, Analytical Reference Materials
Manufacturing Technician/Analyst

Bottom: Diane Thompson, Shipping Technician



Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	cat.#
abamectin	71751-41-2	Neat	100mg	PS-2068
Abate®	3383-96-8	Neat	100mg	PS-665
Acarol	18181-80-1	Neat	250mg	PS-1062
Accent®	111991-09-4	Neat	100mg	PS-2005
acenaphthene	83-32-9	M	1,000	31267
acenaphthylene	208-96-8	M	1,000	31268
acephate	30560-19-1	Neat	1g	PS-738
acetaldehyde-2,4-DNPH	1019-57-4	ACN	100	33074
acetamidrid	135410-20-7	Neat	100mg	PS-2219
acetochlor	34256-82-1	Neat	100mg	PS-2040
acetochlor oa	—	Neat	25mg	MET-2040B
acetone	67-64-1	PTM	5,000	30245
acetone-2,4-DNPH	1567-89-1	ACN	100	33075
acetonitrile	75-05-8	DMSO	2.05mg/mL	36281
acetonitrile	75-05-8	PTM	1,000	30495
acetophenone	98-86-2	PTM	5,000	30621
acibenzolar-S-methyl	135158-54-2	Neat	250mg	PS-2212
acifluorfen	50594-66-6	M	1,000	32255
acifluorfen	50594-66-6	Neat	1g	PS-1016
acifluorfen, amino	—	Neat	10mg	MET-1016A
acifluorfen methyl ester	50594-67-7	M	1,000	32256
acifluorfen methyl ester	—	Neat	100mg	PS-1111
acrolein	107-02-8	PTM	10,000	30499
acrolein	107-02-8	W	10,000	30478
acrolein-2,4-DNPH	888-54-0	ACN	100	33076
acrylamide	79-06-1	M	1,000	30494
acrylonitrile	107-13-1	PTM	2,000	30246
agribrom	16079-88-2	Neat	500mg	PS-2089
Akton	1757-18-2	Neat	1g	PS-675
alachlor	15972-60-8	M	1,000	32204
alachlor	15972-60-8	Neat	1g	PS-357
alachlor oa	—	Neat	25mg	MET-357E
aldicarb	116-06-3	Neat	1g	PS-734
aldicarb sulfoxide	1646-87-3	Neat	10mg	PS-1054
aldoxycarb	1646-88-4	Neat	10mg	PS-1055
aldrin	309-00-2	M	1,000	32205
aldrin	309-00-2	Neat	100mg	PS-69
di-allate	2303-16-4	Neat	100mg	PS-507
tri-allate	2303-17-5	Neat	1g	PS-506
allethrin (5mL)	584-79-2	D	0.1mg/mL	PS-96
alloydim - sodium (5mL)	66003-55-2	ACN	100	PS-1200
allyl chloride	107-05-1	PTM	2,000	30248
4-allyl-2-methoxyphenol	97-53-0	Neat	1g	PS-940
allyl trimethylammonium bromide	—	Neat	1g	PS-360
alprazolam	28981-97-7	PTM	1,000	34042
amdro	67485-29-4	Neat	100mg	PS-1001
ametryne	834-12-8	Neat	500mg	PS-383
2-aminobenzimidazole	934-32-7	Neat	250mg	MET-222A
aminocarb	2032-59-9	Neat	250mg	PS-751
3-amino-2,5-dichlorobenzoic acid	133-90-4	Neat	100mg	PS-314
2-amino-4,6-dinitrotoluene	35572-78-2	ACN	1,000	31670
4-amino-2,6-dinitrotoluene	19406-51-0	ACN	1,000	31671
2-amino-N-isopropyl benzamide	—	Neat	50mg	MET-1011C
aminomethyl phosphonic acid	1066-51-9	Neat	100mg	MET-1051A
aminomethyl phosphonic acid (AMPA)	1066-51-9	W	100	32428
aminopyralid	150114-71-9	Neat	100mg	PS-2315
aminotriazole	61-82-5	Neat	1g	PS-382
amitraz	33089-61-1	Neat	500mg	PS-1005
ammonium fluosilicate	16919-19-0	Neat	1g	PS-291
ammonium picrate	131-74-8	ACN	2,000	31890
ammonium sulfamate	7773-06-0	Neat	1g	PS-344
Amobam	3566-10-7	Neat	1g	PS-228
amobarbital	64-43-7	PTM	1,000	34028
d-amphetamine	51-63-8	PTM	1,000	34020
tert-amyl alcohol	75-85-4	PTM	10,000	30631
tert-amyl ethyl ether (TAE)	919-94-8	PTM	2,000	30617

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	cat.#
tert-amyl methyl ether (TAME)	994-05-8	PTM	2,000	30629
ancymidol	12771-68-5	Neat	1g	PS-1008
5-α-androstane	438-22-2	D	2,000	31065
anilazine	101-05-3	Neat	250mg	PS-388
aniline	62-53-3	M	1,000	31470
anthracene	120-12-7	A	1,000	31269
anthracene-d10	1719-06-08	D	2,000	31037
9,10-anthraquinone	84-65-1	Neat	1g	PS-926
antifoam agent for purge & trap	N/A	Neat	1mL	31822
4-ap	504-24-5	Neat	1g	PS-1010
aprobarbital	77-02-1	PTM	1,000	34029
Aramite	140-57-8	H	2,000	31624
Aramite (5mL)	140-57-8	H	0.1mg/mL	PS-850
Aroclor 1016	12674-11-2	H	1,000	32006
Aroclor 1016	12674-11-2	I	200	32064
Aroclor 1016	12674-11-2	TO	500mg/kg	32076
Aroclor 1016	12674-11-2	TO	50mg/kg	32075
Aroclor 1221	11104-28-2	H	1,000	32007
Aroclor 1221	11104-28-2	I	200	32065
Aroclor 1221	11104-28-2	TO	500mg/kg	32078
Aroclor 1221	11104-28-2	TO	50mg/kg	32077
Aroclor 1232	11141-16-5	H	1,000	32008
Aroclor 1232	11141-16-5	I	200	32066
Aroclor 1232	11141-16-5	TO	500mg/kg	32080
Aroclor 1232	11141-16-5	TO	50mg/kg	32079
Aroclor 1242	53469-21-9	H	1,000	32009
Aroclor 1242	53469-21-9	I	200	32067
Aroclor 1242	53469-21-9	TO	500mg/kg	32082
Aroclor 1242	53469-21-9	TO	50mg/kg	32081
Aroclor 1248	12672-29-6	H	1,000	32010
Aroclor 1248	12672-29-6	I	200	32068
Aroclor 1248	12672-29-6	TO	500mg/kg	32084
Aroclor 1248	12672-29-6	TO	50mg/kg	32083
Aroclor 1254	11097-69-1	H	1,000	32011
Aroclor 1254	11097-69-1	I	200	32069
Aroclor 1254	11097-69-1	TO	500mg/kg	32086
Aroclor 1254	11097-69-1	TO	50mg/kg	32085
Aroclor 1260	11096-82-5	H	1,000	32012
Aroclor 1260	11096-82-5	I	200	32070
Aroclor 1260	11096-82-5	TO	500mg/kg	32088
Aroclor 1260	11096-82-5	TO	50mg/kg	32087
Aroclor 1262	37324-23-5	H	1,000	32409
Aroclor 1268	11100-14-4	H	1,000	32410
Asana®	66230-04-4	Neat	100mg	PS-2004
Aspon	3244-90-4	Neat	50mg	PS-663
asulam	3337-71-1	Neat	1g	PS-1009
atraton	1610-17-9	Neat	500mg	PS-395
atrazine	1912-24-9	A	1,000	32208
atrazine	1912-24-9	Neat	500mg	PS-380
atrazine desethyl	6190-65-4	Neat	50mg	MET-380B
atrazine desethyl des-isopropyl-2-hydroxy	—	Neat	50mg	MET-380E
atrazine desethyl-2-hydroxy	—	Neat	50mg	MET-380D
atrazine desisopropyl	1007-28-9	Neat	50mg	MET-58A
atrazine-2-hydroxy	2163-68-0	Neat	50mg	MET-380F
aviation gas (5mL)	8006-69-1	PTM	50,000	30208
aviation gas	8006-69-1	PTM	2,500	30094
aviation gas	8006-69-1	PTM	50,000	30207
azadirachtin	11141-17-6	Neat	10mg	PS-2075
azamethiphos (5mL)	35575-96-3	D	100	PS-1202
azinphos-methyl oxon (5mL)	—	T	100	MET-666A
azobenzene	103-33-3	D	1,000	31496
azobenzene	103-33-3	Neat	1g	PS-851
azocyclofen	41083-11-8	Neat	250mg	PS-2063
azoxystrobin	131860-33-8	Neat	250mg	PS-2141

*Volume is 1mL/ampul unless otherwise noted. Concentration is µg/mL unless otherwise noted.

Solvent code:

A = acetone
ACN = acetonitrile
C = carbon disulfide
Cy = cyclohexane

D = methylene chloride
DMSO = dimethyl sulfoxide
EA = ethyl acetate
H = hexane

I = isoctane
Ip = isopropanol
M = methanol
MTBE = methyl *tert*-butyl ether

PTM = purge & trap grade methanol
T = toluene
TO = transformer oil
W = water (DI)

Cat.#s starting with PS- & MET- are ChemService pesticides that are subject to size and price changes without notification. They are not available in Europe.

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL* cat.#	Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL* cat.#
banol	671-04-5	Neat	100mg PS-1085	bromoacetyl bromide	598-21-0	Neat	500mg PS-273
barban	101-27-9	Neat	100mg PS-540	<i>p</i> -bromoaniline	106-40-1	Neat	1g MET-374B
barbital	57-44-3	PTM	1,000 34030	bromobenzene	108-86-1	PTM	2,000 30250
Baycarb	3766-81-2	Neat	1g PS-1046	2-bromobutanoic acid	80-58-0	MTBE	2,000 31881
Baycor	55179-31-2	Neat	100mg PS-1093	2-bromobutyrate	3196-15-4	MTBE	2,000 31882
Baygon®	114-26-1	Neat	1g PS-551	2-bromochlorobenzene	694-80-4	PTM	2,000 30228
Baytan®	55219-65-3	Neat	250mg PS-1064	4-bromochlorobenzene	106-39-8	PTM	2,000 30230
Baythroid®	68359-37-5	Neat	250mg PS-1090	1-bromo-2-chloroethane	107-04-0	PTM	2,000 30469
BDMC	—	Neat	100mg PS-450	bromochloromethane	74-97-5	PTM	2,000 30225
Beam	41814-78-2	Neat	100mg PS-1083	2-bromo-1-chloropropane	3017-95-6	PTM	2,000 30226
benalaxyl	71626-11-4	Neat	250mg PS-2142	bromodichloromethane	75-27-4	PTM	2,000 30251
benazolin-ethyl	25059-80-7	Neat	250mg PS-2193	4-bromo-3,5-dimethylphenyl-N- methylcarbamate (BDMC)	3766-81-2	M	100 32274
bendiocarb	22781-23-3	Neat	250mg PS-739	(2-bromoethyl)trimethyl ammonium bromide	2758-06-7	Neat	1g PS-361
benfluralin	1861-40-1	Neat	1g PS-389	1-bromo-4-fluorobenzene	460-00-4	A	1,000 31854
benomyl	17804-35-2	Neat	100mg PS-222	4-bromofluorobenzene	460-00-4	PTM	2,000 30026
bensulfuron-methyl	83055-99-6	Neat	100mg PS-1082	4-bromofluorobenzene	460-00-4	PTM	2,500 30067
bensulide	741-58-2	Neat	250mg PS-425	4-bromofluorobenzene	460-00-4	PTM	5,000 30003
bensultap	17606-31-4	Neat	1g PS-2106	4-bromofluorobenzene	460-00-4	PTM	10,000 30082
bentazon	25057-89-0	M	1,000 32257	bromoforn	75-25-2	PTM	2,000 30252
bentazon	25057-89-0	Neat	250mg PS-1011	bromomethane	74-83-9	PTM	2,000 30253
bentazon methyl derivative	—	Neat	50mg PS-1112	bromomethane solution (5mL)	74-83-9	M	10,000 PS-155
benthiocarb	28249-77-6	Neat	1g PS-1017	1-bromo-2-nitrobenzene	577-19-5	A	1,000 32279
benzaldehyde	100-52-7	D	2,000 33017	<i>p</i> -bromophenoxy acetic acid	—	Neat	1g PS-312
benzaldehyde-2,4-DNPH	1157-84-2	ACN	100 33077	bromophos ethyl	4824-78-6	Neat	100mg PS-2008
benzene	71-43-2	DMSO	10mg/mL 36282	bromophos methyl	2104-96-3	Neat	50mg PS-2009
benzene	71-43-2	D	250 35262	2-bromopropionic acid	598-72-1	MTBE	1,000 31653
benzene	71-43-2	PTM	2,000 30249	bromoxynil	1689-84-5	Neat	1g PS-426
benzene-d6	1076-43-3	PTM	2,000 30025	bromoxynil-heptanoate	—	Neat	50mg MET-426A
benzidine	92-87-5	M	1,000 31441	bromoxynil-methyl ether	—	Neat	100mg MET-426B
benzo(a)anthracene	56-55-3	M	1,000 31270	bromoxynil octanoic acid ester	1689-99-2	Neat	1g PS-434
benzo(a)pyrene	50-32-8	A	1,000 31271	bromuconazole	116255-48-2	Neat	100mg PS-2301
benzo(b)fluoranthene	205-99-2	A	1,000 31272	bronopol	52-51-7	Neat	1g PS-2128
benzo(ghi)perylene	191-24-2	D	1,000 31273	bupropfen	69327-76-0	Neat	100mg PS-2067
benzo(k)fluoranthene	207-08-9	A	1,000 31274	Busan® 40	51026-28-9	Neat	100mg PS-2019
benzoic acid	65-85-0	D	2,000 31879	Busan® 85	128-03-0	Neat	100mg PS-2020
benzoic acid	65-85-0	M	1,000 31415	butabarbital	125-40-6	PTM	1,000 34031
benzoylcegonine	519-09-5	PTM	1,000 34016	butachlor	23184-66-9	Neat	250mg PS-348
benzoylprop ethyl	22212-55-1	Neat	100mg PS-2079	1,3-butadiene	106-99-0	PTM	1,000 30622
benzphetamine	5411-22-3	PTM	1,000 34022	butalbitol	77-26-9	PTM	1,000 34032
benzthiazuron (5mL)	1929-88-0	D	100 PS-1201	1,4-butanediol	110-63-4	M	1,000 34078
benzyl benzoate	120-51-4	H	5,000 31847	(s)-(-)-1,2,4-butanetriol	42890-76-6	pyridine	1,000 33024
benzyl benzoate	120-51-4	Neat	1g PS-904	(s)-(-)-1,2,4-butanetriol (5mL)	42890-76-6	pyridine	1,000 33032
<i>o</i> -benzyl- <i>p</i> -chlorophenol	120-32-1	Neat	1g PS-170	1-butanol	71-36-3	D	250 35260
benzyl thiocyanate	3012-37-1	Neat	1g PS-98	1-butanol	71-36-3	PTM	50,000 30474
α -BHC	319-84-6	M	1,000 32206	<i>tert</i> -butanol	75-65-0	PTM	50,000 30470
α -BHC	319-84-6	Neat	100mg PS-692	<i>tert</i> -butanol-d9	25725-11-5	PTM	20,000 30618
β -BHC	319-85-7	A	1,000 32209	2-butanone (MEK)	78-93-3	PTM	5,000 30254
β -BHC	319-85-7	Neat	100mg PS-691	2-butanone-2,4-DNPH	958-60-1	ACN	100 33078
δ -BHC	319-86-8	M	1,000 32217	2-[2-butoxyethoxy]ethyl thiocyanate	112-56-1	Neat	250mg PS-730
δ -BHC	319-86-8	Neat	100mg PS-693	butoxypolypropylene glycol	9003-13-8	Neat	1g PS-908
ϵ -BHC	6108-10-7	T	100 PS-2222	butralin	33629-47-9	Neat	500mg PS-1004
γ -BHC (lindane)	58-89-9	M	1,000 32226	butylate	2008-41-5	Neat	1g PS-503
BHC mixed isomers	608-73-1	Neat	1g PS-70	4- <i>tert</i> -butyl-2-chloro-phenol	—	Neat	10mg MET-607A
bifenazate	149877-41-8	Neat	250mg PS-2194	2- <i>sec</i> -butyl-4,6-dinitro-phenol triethanolamine salt	6420-47-9	Neat	1g PS-436
bifenox	42576-02-3	Neat	1g PS-311	2-butyl-2-ethyl-1,3-propanediol	115-84-4	Neat	1g PS-101
bifenthrin	82657-04-3	Neat	100mg PS-2003	butylmesityl oxide oxalate	532-34-3	Neat	1g PS-907
binapacryl	485-31-4	Neat	100mg PS-2087	2-(4- <i>tert</i> -butylphenoxy) cyclohexanol	—	Neat	100mg MET-858A
bioallethrin	28057-48-9	Neat	250mg PS-791	<i>n</i> -butyraldehyde-2,4-DNPH	1527-98-6	ACN	100 33079
<i>s</i> -bioallethrin	28434-00-6	Neat	250mg PS-790	γ -butyrolactone (GBL)	96-48-0	ACN	1,000 34077
bioresmethrin	28434-01-7	Neat	500mg PS-2080	cadmium succinate	141-00-4	Neat	1g PS-280
biphenyl	92-52-4	Neat	1g PS-2032	cadusafos	95465-99-9	Neat	100mg PS-2199
bis(<i>p</i> -chlorophenyl)urea	—	Neat	500mg MET-59A	caffeine	74051-80-2	W	5 31800
N,N'-bis(3,4-dichloro-phenyl) urea	—	Neat	100mg MET-60A	caffeine	74051-80-2	W	25 31801
bis(2-ethylhexyl)adipate	103-23-1	M	1,000 31449	caffeine	74051-80-2	W	125 31802
bis(2-ethylhexyl)phthalate	117-81-7	D	1,000 31420	caffeine	74051-80-2	W	250 31803
bis[tri- <i>n</i> -butyltin]oxide	56-35-9	Neat	1g PS-282	caffeine	74051-80-2	W	500 31804
Bitrex®	3734-33-6	Neat	1g PS-2084	caffeine	58-08-2	M	1,000 34084
bladex	21725-46-2	Neat	1g PS-387				
brodifacoum	56073-10-0	Neat	50mg PS-1056				
bromacil	314-40-9	Neat	1g PS-62				
bromazepam	1812-30-2	PTM	1,000 34043				

free data

Available on Our Website: Lot Certificates, Data Packs, and MSDSs

For complete information detailing manufacturing and testing for Restek inventoried reference standards, visit our website at www.restek.com.

To view lot certificates and/or an MSDS, enter the catalog number of the product in the Search feature. For a free data pack (Adobe® PDF file), enter the catalog number and lot number of the product.

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	Individual cat.#
cannabidiol	13956-24-1	PTM	1,000	34011
cannabinol	521-35-7	PTM	1,000	34010
ε-caprolactam	105-60-2	D	2,000	31833
captafol	2425-06-1	Neat	1g	PS-260
captan	133-06-2	Neat	1g	PS-25
carbaryl	63-25-2	Neat	1g	PS-84
carbazole	86-74-8	D	1,000	31836
carbazole	86-74-8	M	1,000	31430
carbendazim	10605-21-7	Neat	100mg	PS-1077
carbofuran	1563-66-2	Neat	1g	PS-754
carbofuran-3-keto	—	Neat	100mg	MET-754E
carbofuranphenol-3-keto	—	Neat	50mg	MET-754F
carbon disulfide	75-15-0	PTM	2,000	30258
carbon tetrachloride	56-23-5	DMSO	20mg/mL	36283
carbon tetrachloride	56-23-5	PTM	2,000	30259
carbophenothion	786-19-6	Neat	1g	PS-88
carbophenothion methyl (5mL)	953-17-3	M	100	MET-88A
carbophenothion methyl-O-analog	—	Neat	50mg	MET-88B
carbophenothion methyl sulfone (5mL)	—	EA	10ng/µL	MET-88C
carbophenothion methyl sulfoxide (5mL)	—	EA	10ng/µL	MET-88D
carbosulfan	55285-14-8	Neat	250mg	PS-1063
carboxin	5234-68-4	Neat	1g	PS-294
carfentrazone-ethyl	128639-02-1	Neat	100mg	PS-2170
cartap hydrochloride	22042-59-7	Neat	500mg	PS-2069
chloral hydrate	302-17-0	ACN	1,000	30609
α-chloralose	15879-93-3	Neat	250mg	PS-2112
β-chloralose	16376-36-6	Neat	1g	PS-2113
chloramben methyl ester	7286-84-2	Neat	100mg	PS-1113
chloranil	118-75-2	Neat	1g	PS-11
chloranacryl	2164-09-2	Neat	1g	PS-63
chlorbenside	103-17-3	Neat	100mg	PS-853
chlorbromuron	13360-45-7	Neat	1g	PS-368
chlordan	57-74-9	Neat	1g	PS-75
chlordan (technical)	57-74-9	H	1,000	32021
chlordan (technical)	57-74-9	I	5,000	32072
chlordan (technical)	57-74-9	M	2,000	32016
α-chlordan	5103-71-9	M	1,000	32207
cis-chlordan	5103-71-9	Neat	10mg	PS-75-1
trans-chlordan	5103-74-2	Neat	10mg	PS-75-2
γ-chlordan	5566-34-7	M	1,000	32227
chloridiazepoxide	438-41-5	PTM	1,000	34044
chloridimeform	6164-98-3	Neat	1g	PS-840
chlorfenapyr	122453-73-0	Neat	100mg	PS-2171
chlorfenvinphos	470-90-6	Neat	100mg	PS-1087
chlorfluorecol-methyl ester	2536-31-4	Neat	100mg	PS-1022
chlorimuron ethyl	90982-32-4	Neat	100mg	PS-1081
chlormephos	24934-91-6	Neat	250mg	PS-2209
chlormequat chloride	999-81-5	Neat	1g	PS-363
chloroacetic acid	79-11-8	Neat	1g	PS-31
2-chloroallyl diethyldithiocarbamate	95-06-7	Neat	250mg	PS-19
4-chloroaniline	106-47-8	D	2,000	31211
chlorobenzene	108-90-7	DMSO	1.8mg/mL	36284
chlorobenzene	108-90-7	PTM	2,000	30261
chlorobenzene-d5	3114-55-4	PTM	2,000	30223
chlorobenzilate	510-15-6	M	1,000	32211
chlorobenzilate	510-15-6	Neat	100mg	PS-854
2-chlorobenzonitrile	873-32-5	Neat	250mg	MET-2015A
2-chloro-4,6-diamino- 1,3,5-triazine	3397-62-4	Neat	100mg	MET-58C
2-chloro-2',6'-diethyl- acetanilide	—	Neat	100mg	MET-357A
chloroethane	75-00-3	PTM	2,000	30263

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	Individual cat.#
2-chloroethanol	107-07-3	PTM	2,000	30264
(2-chloroethyl)phosphonic acid	16672-87-0	Neat	1g	PS-325
2-chloroethyl vinyl ether	110-75-8	PTM	2,000	30265
1-chloro-2-fluorobenzene	348-51-6	PTM	2,000	30040
1-chloro-4-fluorobenzene	352-33-0	PTM	2,500	30066
chloroform	67-66-3	DMSO	0.3mg/mL	36285
chloroform	67-66-3	PTM	2,000	30266
p-chloromandelic acid	492-86-4	Neat	1g	PS-310
chloromethane	74-87-3	PTM	2,000	30267
4-chloro-2-methylaniline	95-69-2	Neat	1g	MET-840A
4-chloro-2-methylphenol	1570-64-5	Neat	1g	MET-40A
2-chloronaphthalene	91-58-7	M	1,000	31284
chloroneb	2675-77-6	Neat	1g	PS-255
1-chloro-3-nitrobenzene	121-73-3	H	1,000	31875
4-chloro-3-nitrobenzotrifluoride	121-17-5	A	1,000	32282
2-chloro-4-nitrophenol	619-08-9	Neat	1g	MET-669A
1-chlorooctadecane	3386-33-2	D	10,000	31098
1-chlorooctane	111-85-3	PTM	10,000	30084
chlorophacinone	3691-35-8	Neat	500mg	PS-2085
o-chlorophenoxy acetic acid	614-61-9	Neat	1g	PS-38
p-chlorophenoxy acetic acid	122-88-3	Neat	1g	PS-39
2-[3-chlorophenoxy] propionamide	—	Neat	1g	PS-359
2-(3-chlorophenoxy) propionic acid	101-10-0	Neat	1g	PS-317
3-(4-chlorophenyl)methyl urea	—	Neat	100mg	MET-60B
p-chlorophenylurea	140-38-5	Neat	100mg	MET-1028A
6-chloro-2-picolinic acid methyl ester (5mL)	—	H	100	MET-419B
chloroprene	126-99-8	PTM	5,000	30238
chlorothalonil	1897-45-6	Neat	1g	PS-1020
4-chlorotoluene	106-43-4	PTM	2,000	30269
chlorotoluron	15545-48-9	Neat	250mg	PS-2078
4-chloro-σ-tolylxyacetic acid	94-74-6	Neat	1g	PS-40
4-chloro-σ-tolylxyacetic acid methyl ester	2436-73-9	Neat	100mg	PS-1105
chloroxuron	1982-47-4	Neat	1g	PS-376
chloroxynil	1891-95-8	Neat	500mg	PS-2090
chlorpropham	101-21-3	Neat	1g	PS-52
chlorpyrifos	2921-88-2	M	1,000	32212
chlorpyrifos	2921-88-2	Neat	1g	PS-674
chlorpyrifos methyl	5598-13-0	Neat	250mg	PS-418
chlorpyrifos methyl-O-analog	—	Neat	50mg	MET-418B
chlorpyrifos oxon	5598-15-2	Neat	100mg	MET-674B
chlorsulfuron	64902-72-3	Neat	100mg	PS-1065
chlorsulfuron-5-hydroxy (5mL)	—	EA	10ng/µL	MET-1065A
chlorthal	1861-32-1	Neat	1g	PS-33
chrysene	218-01-9	A	1,000	31275
clethodim	99129-21-2	Neat	100mg	PS-2045
clobazam	22316-47-8	PTM	1,000	34045
clofentezine	74115-24-5	Neat	100mg	PS-2015
clonazepam	1622-61-3	PTM	1,000	34046
cloransulam-methyl	147150-35-4	Neat	100mg	PS-2165
cocaethylene	529-38-4	ACN	1,000	34066
cocaine	53-21-4	PTM	1,000	34015
codeine	76-57-3	PTM	1,000	34000
command	81777-89-1	Neat	100mg	PS-2001
continine	486-56-6	M	1,000	34086
copper naphthenates	—	Neat	1g	PS-2028
copper oxychloride	1332-65-6	Neat	1g	PS-292
copper-8-quinolinolate	10380-28-6	Neat	500mg	PS-2091
coumaphos	56-72-4	Neat	100mg	PS-656
coumaphos-O-analog	321-54-0	Neat	50mg	MET-656D
coumatetralyl	5836-29-3	Neat	250mg	PS-2093
creatine	57-00-1	Neat	1g	MET-250A

*Volume is 1mL/ampul unless otherwise noted. Concentration is µg/mL unless otherwise noted.

Solvent code:

A=acetone
ACN=acetonitrile
C=carbon disulfide
Cy=cyclohexane

D=methylene chloride
DMSO=dimethyl sulfoxide
EA=ethyl acetate
H=hexane

I=isooctane
Ip=isopropanol
M=methanol
MTBE=methyl *tert*-butyl ether

PTM=purge & trap grade methanol
T=toluene
TO=transformer oil
W=water (DI)

Cat.#s starting with PS- & MET- are ChemService pesticides that are subject to size and price changes without notification. They are not available in Europe.

Individual Compounds: Environmental, Petrochemical, Other

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	cat.#	Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	cat.#
creosote oil	8001-58-9	D	50,000	31838	<i>p,p'</i> -DDMU	1022-22-6	Neat	10mg	MET-699B
crotonaldehyde-2,4-DNPH	1527-96-4	ACN	100	33080	2,4'-DDT	789-02-6	M	1,000	32200
crotoxyphos	7700-17-6	Neat	50mg	PS-603	4,4'-DDT	50-29-3	M	1,000	32203
cyanazine	21725-46-2	A	1,000	32215	<i>o,p'</i> -DDT	789-02-6	Neat	50mg	PS-698
cyanophos	2636-26-2	Neat	100mg	PS-2034	<i>p,p'</i> -DDT	50-29-3	Neat	1g	PS-699
cyclanilide	113136-77-9	Neat	100mg	PS-2243	DDT (mixture <i>p,p'</i> & <i>o,p'</i>)	—	Neat	1g	PS-74
cycloate	1134-23-2	Neat	1g	PS-502	decachlorobiphenyl (BZ #209)	—	—	—	—
cyclohexane	110-82-7	DMSO	19.4mg/mL	36286	(5mL)	2051-24-3	A	200	32030
1,2-cyclohexanediol	—	—	—	—	decachlorobiphenyl (BZ #209)	2051-24-3	A	200	32029
(<i>cis</i> & <i>trans</i> mixture)	931-17-9	Neat	1g	MET-858B	decachlorobiphenyl (BZ #209)	2051-24-3	I	10	32289
cycloheximide	66-81-9	Neat	100mg	PS-1002	decafluorobiphenyl	434-90-2	A	1,000	31855
cycloprate	54460-46-7	Neat	500mg	PS-2092	decafluorobiphenyl	434-90-2	D	2,000	31041
cyclosulfamuron	136849-15-5	Neat	250mg	PS-2214	decafluorobiphenyl	434-90-2	ACN	1,000	31842
lambda-cyhalothrin	91465-08-6	Neat	100mg	PS-2018	decafluorotriphenylphosphine	—	—	—	—
cymoxanil	57966-95-7	Neat	100mg	PS-1067	(DFTPP)	5074-71-5	D	2,500	31001
cypermethrin	52315-07-8	Neat	100mg	PS-1068	<i>n</i> -decane	124-18-5	Neat	1mL	31858
alpha-cypermethrin	67375-30-8	Neat	250mg	PS-2083	deltamethrin	52918-63-5	Neat	100mg	PS-2071
beta-cypermethrin	—	Neat	250mg	PS-2203	demeton O	298-03-3	Neat	100mg	PS-661
cyprazine	22936-86-3	Neat	100mg	PS-2104	demeton S	126-75-0	Neat	100mg	PS-662
cyproconazole	113096-99-4	Neat	50mg	PS-2130	demeton O&S	8065-48-3	Neat	100mg	PS-661/662
cyprodinil	121552-61-2	Neat	500mg	PS-2144	desethyl-atrazine	6190-65-4	A	1,000	32445
cyromazine	66215-27-8	Neat	500mg	PS-2073	desisopropylatrazine	1007-28-9	A	1,000	32446
Cythioate	115-93-5	Neat	1g	PS-1025	desmedipham	13684-56-5	Neat	250mg	PS-1015
2,4-D (2,4-dichlorophenylacetic acid)	94-75-7	M	1,000	32239	Devrinol®	15299-99-7	Neat	250mg	PS-553
2,4-D	94-75-7	Neat	1g	PS-41	Dexon	140-56-7	Neat	1g	PS-262
2,4-D butoxyethyl ester	1929-73-3	Neat	1g	PS-315	dextromethorphan HBr monohydrate	125-69-9	M	1,000	34081
2,4-D butoxypropylene ester	1928-45-6	Neat	1g	PS-298	dextro-propoxyphene	1639-60-7	PTM	1,000	34008
2,4-D butyl ester	94-80-4	Neat	1g	PS-42	diafenthion	80060-09-9	Neat	250mg	PS-2221
2,4-D <i>sec</i> -butyl ester	—	Neat	1g	PS-301	N,N-diallyl-2-chloroacetamide	93-71-0	Neat	1g	PS-351
2,4-D dimethylamine salt	2008-39-1	Neat	1g	PS-337	diazepam	439-14-5	PTM	1,000	34047
2,4-D ethanolamine salt	3599-58-4	Neat	1g	PS-338	diazinon	333-41-5	Neat	1g	PS-90
2,4-D ethyl ester	533-23-3	Neat	1g	PS-304	diazinon-O-analog	—	Neat	100mg	MET-90A
2,4-D isooctyl ester	25168-26-7	Neat	1g	PS-319	dibenz(a,h)anthracene	53-70-3	D	1,000	31276
2,4-D isopropanol amine salt	—	Neat	1g	PS-332	4,4'-dibromobiphenyl	92-86-4	D	2,000	31039
2,4-D isopropylamine salt	—	Neat	1g	PS-299	4,4'-dibromobiphenyl	92-86-4	EA	500	32092
2,4-D isopropyl ester	94-11-1	Neat	1g	PS-43	dibromochloromethane (chlorodibromomethane)	124-48-1	PTM	2,000	30271
2,4-D methyl ester	1928-38-7	M	1,000	32240	1,2-dibromo-3-chloropropane (DBCP)	96-12-8	Neat	1g	PS-1
2,4-D methyl ester	1928-38-7	Neat	1g	PS-41-1	1,2-dibromo-3-chloropropane (DBCP)	96-12-8	PTM	2,000	30270
2,4-D propylene glycol ester	—	Neat	1g	PS-302	dibromofluoromethane	1868-53-7	Neat	100mg	30634
Dacthal® (DCPA dimethyl ester)	1861-32-1	M	1,000	32216	1,2-dibromoethane	106-93-4	Neat	1g	PS-154
dalapon	75-99-0	ACN	1,000	32432	1,2-dibromoethane	106-93-4	PTM	2,000	30272
dalapon	75-99-0	M	1,000	32253	dibromomethane	74-95-3	PTM	2,000	30430
dalapon	75-99-0	M	2,000	32056	4,4'-dibromooctafluorobiphenyl	10386-84-2	D	2,000	31040
dalapon	75-99-0	Neat	1g	PS-32	4,4'-dibromooctafluorobiphenyl	10386-84-2	H	250	32053
dalapon methyl ester	17640-02-7	H	2,000	32057	4,4'-dibromooctafluorobiphenyl	10386-84-2	MTBE	2,000	31856
dalapon methyl ester	17640-02-7	M	1,000	32254	2,3-dibromopropionic acid	600-05-5	MTBE	1,000	31655
dalapon methyl ester	—	Neat	100mg	PS-1100	2,5-dibromotoluene	615-59-8	PTM	1,000	30435
daminozide	1596-84-5	Neat	1g	PS-391	2,5-dibromotoluene	615-59-8	PTM	10,000	30453
dazomet	533-74-4	Neat	1g	PS-17	3-(2,2-dibromovinyl)-2,2-dimethyl-(1-cyclopropane) carboxylic acid	63597-73-9	Neat	100 $\mu\text{g/mL}$	MET-2071B
2,4-DB	94-82-6	M	1,000	32241	dibutyl succinate	141-03-7	Neat	1g	PS-903
2,4-DB	94-82-6	Neat	1g	PS-306	dibutylchlorendate	1770-80-5	A	200	32025
2,4-DB dimethylamine salt	—	Neat	1g	PS-316	dicamba	1918-00-9	M	1,000	32247
2,4-DB methyl ester	18625-12-2	Neat	100mg	PS-1101	dicamba	1918-00-9	Neat	1g	PS-346
DCPA diacid	2136-79-0	M	200	32261	dicamba methyl ester	6597-78-0	M	1,000	32248
<i>p,p'</i> -DDA	83-05-6	Neat	100mg	MET-698A	dicamba methyl ester	6597-78-0	Neat	100mg	PS-1102
2,4'-DDD	53-19-0	M	1,000	32098	dicapthion	2463-84-5	Neat	1g	PS-669
4,4'-DDD	72-54-8	M	1,000	32201	dichlofenthion	97-17-6	Neat	1g	PS-103
<i>o,p'</i> -DDD	53-19-0	Neat	100mg	PS-694	dichlofluanid	1085-98-9	Neat	1g	PS-2047
<i>p,p'</i> -DDD	72-54-8	Neat	1g	PS-72	dichlone	117-80-6	Neat	1g	PS-12
<i>o,p'</i> -DDD olefin	—	Neat	10mg	MET-694A	dichlormid	37764-25-3	Neat	100mg	PS-2054
2,4'-DDE	3424-82-6	M	1,000	32099	2,4-dichloroaniline	554-00-7	Neat	1g	MET-372A
4,4'-DDE	72-55-9	M	1,000	32202	2,6-dichlorobenzamide	—	Neat	1g	MET-392B
<i>o,p'</i> -DDE	3424-82-6	Neat	50mg	PS-695					
<i>p,p'</i> -DDE	72-55-9	Neat	100mg	PS-696					

*Volume is 1mL/ampul unless otherwise noted. Concentration is $\mu\text{g/mL}$ unless otherwise noted.

Solvent code:

A=acetone
ACN=acetonitrile
C=carbon disulfide
Cy=cyclohexane

D=methylene chloride
DMSO=dimethyl sulfoxide
EA=ethyl acetate
H=hexane

I=isooctane
Ip=isopropanol
M=methanol
MTBE=methyl *tert*-butyl ether

PTM=purge & trap grade methanol
T=toluene
TO=transformer oil
W=water (DI)

Cat.#s starting with PS- & MET- are ChemService pesticides that are subject to size and price changes without notification. They are not available in Europe.

Compound	Solvent	Individual	Compound	Solvent	Individual
Packaged 1mL/ampul*	CAS#	Code	$\mu\text{g/mL}^*$	cat.#	
1,2-dichlorobenzene	95-50-1	M	1,000	31442	
1,3-dichlorobenzene	541-73-1	M	1,000	31443	
1,4-dichlorobenzene	106-46-7	ACN	1,000	30498	
1,4-dichlorobenzene	106-46-7	M	1,000	31444	
<i>o</i> -dichlorobenzene	95-50-1	Neat	1g	PS-2	
<i>p</i> -dichlorobenzene	106-46-7	Neat	1g	PS-3	
1,2-dichlorobenzene-d4	2199-69-1	PTM	2,000	30049	
4,4'-dichlorobenzhydrol	90-97-1	Neat	1g	MET-854A	
3,3'-dichlorobenzidine	91-94-1	D	2,000	31835	
3,3'-dichlorobenzidine	91-94-1	M	2,000	31026	
3,3'-dichlorobenzidine-free base	91-94-1	Neat	100mg	31884	
2,6-dichlorobenzoic acid	50-30-6	Neat	1g	MET-392A	
3,5-dichlorobenzoic acid	51-36-5	MTBE	1,000	31652	
3,5-dichlorobenzoic acid	51-36-5	Neat	1g	PS-1107	
2,6-dichlorobenzoic acid methyl ester	—	Neat	250mg	MET-392C	
3,5-dichlorobenzoic acid methyl ester	2905-67-1	M	1,000	32264	
3,5-dichlorobenzoic acid methyl ester	2905-67-1	MTBE	1,000	31649	
3,5-dichlorobenzoic acid methyl ester	—	Neat	100mg	PS-1114	
2,6-dichlorobenzonitrile	1194-65-6	Neat	1g	PS-392	
4,4'-dichlorobenzophenone	90-98-2	Neat	100mg	MET-854B	
1,4-dichlorobutane	110-56-5	PTM	2,000	30227	
<i>trans</i> -1,4-dichloro-2-butene	110-57-6	PTM	2,000	30274	
dichlorodifluoromethane (CFC-12)	75-71-8	PTM	2,000	30275	
1,1-dichloroethane	75-34-3	PTM	2,000	30276	
1,2-dichloroethane	107-06-2	DMSO	25mg/mL	36288	
1,2-dichloroethane	107-06-2	PTM	2,000	30277	
1,2-dichloroethane-d4	17060-07-0	PTM	2,000	30027	
1,1-dichloroethene	75-35-4	DMSO	40mg/mL	36287	
1,1-dichloroethene	75-35-4	PTM	2,000	30278	
<i>cis</i> -1,2-dichloroethene	156-59-2	PTM	2,000	30279	
<i>trans</i> -1,2-dichloroethene	156-60-5	PTM	2,000	30280	
<i>cis</i> -1,2-dichloroethylene	156-59-2	DMSO	4.67mg/mL	36289	
<i>trans</i> -1,2-dichloroethylene	156-60-5	DMSO	4.67mg/mL	36290	
3,6-dichloro-2-hydroxy benzoic acid	—	Neat	50mg	MET-346B	
1,1-dichloro-1-nitro-ethane	594-72-9	Neat	1g	PS-5	
dichlorophen	97-23-4	Neat	1g	PS-6	
2,6-dichlorophenol	87-65-0	M	1,000	31409	
N-[(2,4-dichlorophenoxy) acetyl]-DL-methionine	—	Neat	10mg	PS-390	
2,4-dichlorophenylacetic acid	19719-28-9	M	200	32049	
2,4-dichlorophenyl acetic acid methyl ester	N/A	H	200	32050	
2,4-dichlorophenyl benzenesulfonate	97-16-5	Neat	1g	PS-856	
1,2-dichloropropane	78-87-5	Neat	1g	PS-153	
1,2-dichloropropane	78-87-5	PTM	2,000	30281	
2,2-dichloropropane	594-20-7	PTM	2,000	30283	
1,3-dichloropropane	542-75-6	Neat	1g	PS-152	
<i>cis</i> -1,3-dichloropropane	10061-01-5	PTM	2,000	30284	
<i>trans</i> -1,3-dichloropropane	10061-02-6	PTM	2,000	30285	
3',4'-dichloropropionanilide	709-98-8	Neat	500mg	PS-356	
2,3-dichloropropionic acid	565-64-0	MTBE	1,000	31650	
2,3-dichloropropionic acid methyl ester	3674-09-7	MTBE	1,000	31651	
1,2-dichlorotetrafluoroethane (CFC-114)	76-14-2	PTM	2,000	30476	
3-(2,2-dichlorovinyl)-2,2-dimethyl-(1-cyclopropane) carboxylic acid	55701-05-8	Neat	1mL	MET-758D	
dichlorprop	120-36-5	M	1,000	32249	
dichlorprop	120-36-5	Neat	250mg	PS-44	
dichlorprop methyl ester	57153-17-0	M	1,000	32250	
dichlorprop methyl ester	—	Neat	100mg	PS-1103	
dichlorvos	62-73-7	Neat	1g	PS-89	
diclofop acid	40843-25-2	Neat	100mg	MET-1036B	
diclofop methyl	51338-27-3	Neat	1g	PS-1036	
dicloran	99-30-9	Neat	1g	PS-293	
diclosulam	145701-21-9	Neat	250mg	PS-2208	
<i>o,p</i> -dicofol (5mL)	—	ACN	100	MET-82A	
dicrotophos	141-66-2	Neat	250mg	PS-602	
dicyclopentadiene	77-73-6	Neat	500mg	PS-2033	
dieldrin	60-57-1	M	1,000	32218	
dieldrin	60-57-1	Neat	250mg	PS-76	
diesel fuel #2 composite (5mL)	68334-30-5	D	50,000	31259	

Compound	Solvent	Individual	Compound	Solvent	Individual
Packaged 1mL/ampul*	CAS#	Code	$\mu\text{g/mL}^*$	cat.#	
diesel fuel #2 composite	68334-30-5	D	5,000	31093	
diesel fuel #2 composite	68334-30-5	D	50,000	31258	
diesel fuel #2: 25% weathered	68334-30-5	D	5,000	31234	
diesel fuel #2: 50% weathered	68334-30-5	D	5,000	31235	
diesel fuel #2: 75% weathered	68334-30-5	D	5,000	31236	
diesel fuel #2: unweathered	68334-30-5	D	5,000	31233	
diesel/biodiesel 80:20	67784-80-9	D	5,000	31880	
diethyl ethyl	38727-55-8	Neat	1g	PS-1006	
diethofencarb	87130-20-9	Neat	250mg	PS-2153	
2-diethylamino-6-methyl pyrimidin-4-ol	—	Neat	100mg	MET-644A	
2,6-diethylaniline	579-66-8	Neat	500mg	MET-357C	
diethyl dithiobis (thionoformate)	502-55-6	Neat	1g	PS-66	
diethyl ether (ethyl ether)	60-29-7	PTM	2,000	30286	
diethyl phosphate	598-02-7	Neat	100mg	MET-90C	
O,O-diethyl phosphorochloridothioate	2524-04-1	Neat	1g	PS-670	
N,N-diethyl- <i>m</i> -toluamide	134-62-3	Neat	1g	PS-902	
difenacoum	56073-07-5	Neat	100mg	PS-2143	
difenoconazole	119446-68-3	Neat	250mg	PS-2159	
difenoquat-desmethyl	—	Neat	100mg	MET-404C	
difenoquat methylsulfate	43222-48-6	Neat	1g	PS-404	
diflubenuron	35367-38-5	Neat	1g	PS-1028	
diflufenican	—	Neat	100mg	PS-2316	
diflufenzopyr	109293-97-2	Neat	250mg	PS-2148	
1,4-difluorobenzene	540-36-3	PTM	2,000	30032	
2,3-dihydro-2,2-dimethyl benzofuran-3,7-diol	—	T	10ng/ μL	MET-754D	
2,3-dihydro-2,2-dimethyl-benzofuran-7-ol	1563-38-8	Neat	1g	MET-754B	
3,5-diiodo-4-hydroxybenzonitrile	1689-83-4	Neat	250mg	PS-397	
diisopropyl ether (DIPE)	108-20-3	PTM	2,000	30627	
dikegulac acid	18467-77-1	Neat	250mg	PS-2190	
dimethypo	52207-48-4	Neat	250mg	PS-2184	
dimethachlon	24096-53-5	Neat	500mg	PS-2154	
dimethenamid	87674-68-8	Neat	250mg	PS-2161	
dimethenamid-P	163515-14-8	Neat	100mg	PS-2291	
dimethipin	55290-64-7	Neat	1g	PS-2052	
dimethoate	60-51-5	Neat	250mg	PS-659	
dimethomorph	110488-70-5	Neat	500mg	PS-2138	
1,2-dimethoxyethane	173201-80-4	DMSO	0.5mg/mL	36291	
dimethylarsinic acid	75-60-5	Neat	500mg	PS-51	
N,N-dimethylacetamide	127-19-5	DMSO	5.45mg/mL	36292	
2,5-dimethylbenzaldehyde-2,4-DNPH	152477-96-8	ACN	100	33081	
2,2-dimethyl-1,3-benzodioxole-4-ol	—	Neat	500mg	MET-739A	
dimethylchlorosilane (DMDCS) (20mL)	75-78-5	Neat	20mL	31840	
N,N-dimethylformamide	68-12-2	DMSO	4.4mg/mL	36293	
dimethyl- <i>p</i> -nitrophenylphosphate	950-35-6	Neat	100mg	PS-613	
2,6-dimethylphenol	576-26-1	M	1,000	31410	
1,1-dimethyl-3-phenylurea	101-42-8	Neat	1g	PS-61	
O,O-dimethyl phosphorochloridothioate	2524-03-0	Neat	1g	PS-91	
dimethylphthalate	131-11-3	Neat	1g	PS-901	
N,N-dimethylsuccinamic acid	2564-95-6	Neat	1g	PS-340	
di- <i>n</i> -butylphthalate	84-74-2	Neat	1g	PS-900	
diniconazole	83657-24-3	Neat	100mg	PS-2152	
3,5-dinitroaniline	618-87-1	ACN	1,000	31661	
1,2-dinitrobenzene	528-29-0	M	1,000	31453	
1,2-dinitrobenzene	99-65-0	ACN	1,000	31662	
1,4-dinitrobenzene	100-25-4	ACN	2,000	33205	
4,6-dinitro-2-sec-butyl-phenol ethanalamine salt	—	Neat	1g	PS-433	
4,6-dinitro- <i>o</i> -cresol	534-52-1	Neat	1g	PS-57	
2,4-dinitrophenol	51-28-5	M	1,000	31291	
2,4-dinitrotoluene	121-14-2	ACN	1,000	31663	
2,6-dinitrotoluene	606-20-2	ACN	1,000	31664	
3,4-dinitrotoluene	610-39-9	EA	2,000	33901	
3,4-dinitrotoluene	610-39-9	M	1,000	31452	
dinocap	39300-45-3	Neat	100mg	PS-27	

Individual Compounds: Environmental, Petrochemical, Other

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	Individual cat.#	Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	Individual cat.#
di-n-octyl phthalate	117-84-0	M	1,000	31426	ethidimuron	30043-49-3	Neat	100mg	PS-2061
dinoseb	88-85-7	M	1,000	32251	ethiofencarb	29973-13-5	Neat	100mg	PS-2058
dinoseb	88-85-7	Neat	100mg	PS-56	ethion	563-12-2	Neat	1g	PS-92
dinoseb methyl ether	6099-79-2	M	1,000	32252	ethofumesate	26225-79-6	Neat	1g	PS-1045
dinoseb methyl ether	6099-79-2	Neat	100mg	PS-1104	2-ethoxyethanol	110-80-5	DMSO	0.8mg/mL	36295
di-n-propyl isocin-chomeronate	136-45-8	Neat	1g	PS-906	ethoxyquin	91-53-2	Neat	250mg	PS-432
diolein (1,3-di[<i>cis</i> -octadecenyl] glycerol)	2465-32-9	pyridine	5,000	33022	ethoxysulfuron	126801-58-9	Neat	100mg	PS-2303
dioxacarb	6988-21-2	Neat	250mg	PS-2022	ethyl-3-amino- 2,5-dichlorobenzoate	—	Neat	1g	PS-326
1,4-dioxane	123-91-1	DMSO	1.9mg/mL	36294	ethylbenzene	100-41-4	DMSO	1.84mg/mL	36296
1,4-dioxane	123-91-1	D	2,000	31853	ethylbenzene	100-41-4	PTM	2,000	30290
1,4-dioxane	123-91-1	PTM	2,000	30287	ethylbenzene-d5	20302-26-5	PTM	2,000	30028
1,4-dioxane-d8	17647-74-4	PTM	2,000	30614	ethylbenzene-d10	25837-05-2	PTM	2,000	30029
dioxathion	78-34-2	Neat	100mg	PS-658	ethyl- <i>tert</i> -butyl ether (ETBE)	637-92-3	PTM	2,000	30628
dioxymeton-S-methyl	17040-19-6	Neat	100mg	PS-642	N-ethylcyclohexylamine	5459-93-8	Neat	1g	MET-502A
dipacinone	82-66-6	Neat	100mg	PS-288	S-ethyl dipropylthiocarbamate	759-94-4	Neat	500mg	PS-54
diphenamid	957-51-7	Neat	1g	PS-350	ethylenediamine	107-15-3	M	540	35222
diphenylacetone	86-29-3	Neat	1g	PS-64	ethylene glycol	3775-85-7	DMSO	3.1mg/mL	36297
diphenylamine	122-39-4	Neat	1g	PS-427	ethylene oxide	75-21-8	DMSO	500	36005
1,2-diphenylhydrazine	122-66-7	M	1,000	31497	2-ethyl-1,3-hexanediol	94-96-2	Neat	1g	PS-102
dipropetryn	4147-51-7	Neat	100mg	PS-411	2-ethylhexyl-2,4-dichlorophenoxy acetate	1928-43-4	Neat	1g	PS-319-1
diquat dibromide monohydrate	6385-62-2	Neat	1g	PS-365	ethylmercurichlorethiimide	2597-93-5	Neat	1g	PS-13
disodium methyl arsenate	144-21-8	Neat	500mg	PS-281	ethyl methacrylate	97-63-2	PTM	2,000	30289
disulfoton	298-04-4	Neat	1g	PS-652	2-[(2-ethyl-6-methyl-phenyl) amino]-1-propanol	—	Neat	50mg	MET-403A
disulfoton sulfone	2497-06-5	Neat	100mg	MET-652A	ethyl-1-naphthalene acetate	2122-70-5	Neat	1g	PS-335
disulfoton sulfoxide	2497-07-6	Neat	100mg	MET-652B	famphur	52-85-7	Neat	1g	PS-671
ditalimfos	5131-24-8	Neat	100mg	PS-2127	famphur-O-analog	—	Neat	50mg	MET-671A
ditalimfos (5mL)	5131-24-8	H	100	PS-1208	fenac	85-34-7	Neat	1g	PS-308
dithianon	3347-22-6	Neat	250mg	PS-2035	fenamiphos sulfone	—	Neat	100mg	MET-612B
dithiopyr	97886-45-8	Neat	100mg	PS-2173	fenamiphos sulfoxide	—	Neat	50mg	MET-612A
diuron	330-54-1	ACN	200	32450	fenarimol	60168-88-9	Neat	100mg	PS-1073
diuron	330-54-1	Neat	500mg	PS-60	fenbuconazole	114369-43-6	Neat	100mg	PS-2072
dodemorphan	1593-77-7	Neat	100mg	PS-1086	fenbutatin-oxide	13356-08-6	Neat	1g	PS-1060
dodine	2439-10-3	Neat	1g	PS-250	fenchlorphos	299-84-3	Neat	100mg	PS-657
Dyfonate®	944-22-9	Neat	100mg	PS-664	fenchlorphos-oxon	—	Neat	100mg	MET-657A
Dylox®	52-68-6	Neat	1g	PS-608	fenfluramine	16105-77-4	PTM	1,000	34023
ecgonine	5796-31-6	PTM	1,000	34017	fenhexamid	126833-17-8	Neat	250mg	PS-2202
ecgonine methyl ester	38969-40-3	PTM	1,000	34018	fenitrothion	122-14-5	Neat	500mg	PS-678
EDDP perchlorate	66729-78-0	M	1,000	34069	fenitrothion-O-analog (5mL)	2255-17-6	H	100	MET-678A
edifenphos	17109-49-8	Neat	500mg	PS-2066	fenoxaprop ethyl	66441-23-4	Neat	100mg	PS-1088
EGDN	628-96-6	M	1,000	31601	fenoxaprop-P-ethyl	71283-80-2	Neat	100mg	PS-2223
EGT	2514-53-6	Neat	1g	PS-328	fenoxycarb	79127-80-3	Neat	100mg	PS-2174
endosulfan	115-29-7	Neat	500mg	PS-81	fenpropathrin	64257-84-7	Neat	250mg	PS-2002
endosulfan I	959-98-8	M	1,000	32221	fenpropimorph (5mL)	67564-91-4	T	100	PS-1214
endosulfan II	33213-65-9	M	1,000	32222	fenpyroximate	134098-61-6	Neat	100mg	PS-2216
endosulfan alcohol	2157-19-9	Neat	50mg	MET-81A	fensulfthion	115-90-2	Neat	100mg	PS-667
endosulfan ether	3369-52-6	Neat	50mg	MET-81B	fentanyll	437-38-7	M	1,000	34082
endosulfan alpha isomer	959-98-8	Neat	100mg	PS-81-1	nor-fentanyll oxalate	1609-66-1	M	1,000	34083
endosulfan beta isomer	33213-65-9	Neat	100mg	PS-81-2	fenthion	55-38-9	Neat	500mg	PS-655
endosulfan lactone	3868-61-9	Neat	50mg	MET-81C	fenthion-O-analog (5mL)	6552-12-1	I	100	MET-655C
endosulfan sulfate	1031-07-8	M	1,000	32223	fenthion-sulfone	3761-42-0	Neat	100mg	MET-655A
endosulfan sulfate	1031-07-8	Neat	50mg	PS-81-3	fenthion sulfoxide	3761-41-9	Neat	100mg	MET-655D
endothol	145-73-3	Neat	250mg	PS-345	fenthion sulfoxide (5mL)	3761-41-9	D	100	MET-655B
endrin	72-20-8	M	1,000	32219	fentin hydroxide	76-87-9	Neat	500mg	PS-1031
endrin	72-20-8	Neat	250mg	PS-77	fenuron TCA	4482-55-7	Neat	1g	PS-364
endrin aldehyde	7421-93-4	M	1,000	32224	ferbam	14484-64-1	Neat	1g	PS-20
endrin aldehyde	7421-93-4	Neat	10mg	PS-77-1	fipronil	120068-37-3	Neat	100mg	PS-2136
endrin ketone	53494-70-5	M	1,000	32220	flamprop isopropyl	52756-22-6	Neat	100mg	PS-2037
endrin ketone	53494-70-5	Neat	10mg	PS-77-2	flamprop-methyl	52756-25-9	Neat	100mg	PS-2038
EPN	2104-64-5	Neat	1g	PS-93	flocoumafen	90035-08-8	Neat	10mg	PS-2149
ethalfuralin	55283-68-6	Neat	1g	PS-408	fluzifop (5mL)	69335-91-7	ACN	100	PS-1205
ethametsulfuron-methyl	97780-06-8	Neat	100mg	PS-2183	fluzifop-butyl	69806-50-4	Neat	250mg	PS-2240
ethanedial dioxime	557-30-2	Neat	1g	PS-412	fluzifop- <i>p</i> -butyl	79241-46-6	Neat	100mg	PS-1097
ethanol	64-17-5	PTM	2,000	30288					
ethanol	64-17-5	W	10,000	30466					

*Volume is 1mL/ampul unless otherwise noted. Concentration is $\mu\text{g/mL}$ unless otherwise noted.

Solvent code:

A=acetone	D=methylene chloride	I=isooctane	PTM=purge & trap grade methanol
ACN=acetonitrile	DMSO=dimethyl sulfoxide	Ip=isopropanol	T=toluene
C=carbon disulfide	EA=ethyl acetate	M=methanol	TO=transformer oil
Cy=cyclohexane	H=hexane	MTBE=methyl <i>tert</i> -butyl ether	W=water (DI)

Cat.#s starting with PS- & MET- are ChemService pesticides that are subject to size and price changes without notification. They are not available in Europe.

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	Individual cat.#
fluchloralin	33245-39-5	Neat	50mg	PS-1012
flucythrinate	70124-77-5	Neat	30mg	PS-1072
fludioxonil	131341-86-1	Neat	100mg	PS-2172
flufenacet	142459-58-3	Neat	250mg	PS-2205
flufenoxuron	101463-69-8	Neat	250mg	PS-2213
flumetralin	62924-70-3	Neat	250mg	PS-2168
flumetralin (5mL)	62924-70-3	A	100	PS-1211
flumetsulam	98967-40-9	Neat	500mg	PS-2111
flumiclorac-pentyl	87546-18-7	Neat	500mg	PS-2119
flunitrazepam	1622-62-4	PTM	1,000	34049
fluometuron	2164-17-2	Neat	1g	PS-369
fluoranthene	206-44-0	M	1,000	31277
fluorene	86-73-7	M	1,000	31278
2-fluoroacetamide	640-19-7	Neat	1g	PS-2029
fluorobenzene	462-06-6	PTM	2,000	30030
2-fluorobiphenyl	321-60-8	D	2,000	31091
2-fluorobiphenyl	321-60-8	D	10,000	31096
fluorodifen	15457-05-3	Neat	1g	PS-431
1-fluoronaphthalene	321-38-0	D	2,000	31092
2-fluorophenol	367-12-4	D	2,000	31047
flurazepam	1172-18-5	PTM	1,000	34050
fluroxypyr-meptyl	81406-37-3	Neat	100mg	PS-2191
flurprimidol	56425-91-3	Neat	100mg	PS-2131
flusilazole	85509-19-9	Neat	100mg	PS-2041
flutolanil	66332-96-5	Neat	1g	PS-2057
flutriaol	76674-21-0	Neat	100mg	PS-2177
Folicur®	107534-96-3	Neat	100mg	PS-1091
folpet	133-07-3	Neat	250mg	PS-26
fomesafen	72178-02-0	Neat	500mg	PS-2102
forchlorfenuron	68157-60-8	Neat	100mg	PS-2150
formaldehyde (in water)	50-00-0	Neat	1g	PS-2031
formaldehyde-DNPH	1081-15-8	ACN	500	31837
formaldehyde-2,4-DNPH	1081-15-8	ACN	100	33082
formaldehyde oxazoladine	—	T	2,000	33004
formamide	75-12-7	DMSO	1.1mg/mL	36298
formetanate hydrochloride	23422-53-9	Neat	1g	PS-400
formothion	2540-82-1	Neat	100mg	PS-2036
fosamine ammonium	25954-13-6	Neat	1g	PS-414
fosetyl-aluminum	39148-24-8	Neat	100mg	PS-2026
fluberidazole	3878-19-1	Neat	100mg	PS-2059
fuel oil # 4	N/A	D	5,000	31216
fuel oil # 4	N/A	D	50,000	31244
fuel oil # 5	N/A	D	5,000	31217
fuel oil # 5	N/A	D	50,000	31246
fuel oil # 6 (5mL)	N/A	D	50,000	31249
fuel oil # 6	N/A	D	5,000	31218
fuel oil # 6	N/A	D	50,000	31248
2-furanacrylic acid	539-47-9	Neat	1g	PS-341
gibberellic acid	77-06-5	Neat	100mg	PS-49
glufosinate-ammonium	77182-82-2	Neat	250mg	PS-2124
DL-glutethimide	18389-24-7	PTM	1,000	34058
glycerin	56-81-5	pyridine	500	33020
glycolaldehyde-2,4-DNPH	—	ACN	100	33091
glyodin	556-22-9	Neat	1g	PS-261
glyphosate	1071-83-6	Neat	1g	PS-1051
glyphosate (5mL)	1071-83-6	W	1,000	32427
glyphosate	1071-83-6	W	1,000	32426
glyphosate-isopropylammonium	38641-94-0	Neat	250mg	PS-2185
glyphosine	2439-99-8	Neat	1g	PS-2164
Guthion®	86-50-0	Neat	1g	PS-666
Guthion® ethyl	2642-71-9	Neat	100mg	PS-673
halofenozide	112226-61-6	Neat	500mg	PS-2146
halosulfuron-methyl	100784-20-1	Neat	100mg	PS-2108
1,6-HDIP	72375-27-0	DMSO	1,000	33002
heptachlor	76-44-8	M	1,000	32228
heptachlor	76-44-8	Neat	100mg	PS-78
heptachlor epoxide (isomer B)	1024-57-3	M	1,000	32230
heptachlor epoxide (isomer B)	1024-57-3	Neat	50mg	PS-700
trans-heptachlor epoxide	28044-83-9	Neat	50mg	PS-700-1

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	Individual cat.#
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ #180)	35065-29-3	I	10	32288
heptenophos	23560-59-0	Neat	250mg	PS-2225
hexachloroacetone	116-16-5	Neat	1g	PS-272
hexachlorobenzene	118-74-1	A	1,000	32231
hexachlorobenzene	118-74-1	Neat	1g	PS-690
2,2',3,4,4',5'-hexachlorobiphenyl (BZ #138)	35065-28-2	I	10	32286
2,2',4,4',5,5'-hexachlorobiphenyl (BZ #153)	35065-27-1	I	10	32287
hexachloro-1,3-butadiene	87-68-3	M	1,000	31435
hexachlorocyclopentadiene	77-47-4	M	1,000	32232
hexachlorodimethyl sulfone	—	Neat	1g	PS-171
hexachloroethane	67-72-1	M	1,000	31436
hexachlorophene	70-30-4	D	2,000	31811
hexachlorophene	70-30-4	Neat	1g	PS-1035
hexaconazole	79983-71-4	Neat	500mg	PS-2157
hexaflumuron	86479-06-3	Neat	250mg	PS-2074
hexaldehyde-2,4-DNPH	1527-97-5	ACN	100	33083
hexane	8031-34-3	DMSO	1.45mg/mL	36299
hexobarbital	56-29-1	PTM	1,000	34033
hexythiazox	78587-05-0	Neat	100mg	PS-1079
HMX	2691-41-0	ACN	1,000	31665
hydraulic oil	—	D	50,000	31839
hydrocodone	34195-34-1	PTM	1,000	34002
hydromorphone	71-68-1	PTM	1,000	34063
hydroprene	41096-46-2	Neat	50mg	PS-2134
γ-hydroxybutyrate (GHB)	502-85-2	M	1,000	34076
4-hydroxycoumarin	1076-38-6	Neat	1g	MET-104B
7-hydroxycoumarin	93-35-6	Neat	1g	MET-104A
5-hydroxy-diclofop methyl (5mL)	—	Cy	10ng/µL	MET-1036A
p-hydroxy-N,N-dimethyl-benzene sulfonamide	—	Neat	10mg	MET-671B
2-hydroxyethyl-n-octyl sulphide	3547-33-9	Neat	1g	PS-952
2-hydroxypyrazine	—	M	100	MET-668B
8-hydroxyquinoline	148-24-3	Neat	1g	PS-18
hymexazole	10004-44-1	Neat	25mg	PS-2186
imazalil	35554-44-0	Neat	100mg	PS-2043
imazalil sulfate	60534-80-7	Neat	100mg	PS-2043-1
imazamethabenz-methyl	81405-85-8	Neat	100mg	PS-2195
imazamox	114311-32-9	Neat	500mg	PS-2162
imazapic	104098-48-8	Neat	250mg	PS-2197
imazapyr	81334-34-1	Neat	100mg	PS-2016
imazaquin	81335-37-7	Neat	100mg	PS-2053
imazethapyr	81335-77-5	Neat	100mg	PS-2039
imidacloprid	138261-41-3	Neat	500mg	PS-2086
Imidan	732-11-6	Neat	1g	PS-653
2-imidazolidinethione	96-45-7	Neat	1g	PS-2180
iminobis (methyl-phosphonic acid)	6711-48-4	Neat	500mg	MET-2164A
indeno(1,2,3-cd)pyrene	193-39-5	D	1,000	31279
3-indoleacetic acid	87-51-4	Neat	1g	PS-320
3-indolebutyric acid	133-32-4	Neat	100mg	PS-48
3-indolepropionic acid	830-96-6	Neat	1g	PS-321
iodofenphos (5mL)	18181-70-9	A	100	PS-1206
iodomethane	74-88-4	PTM	2,000	30292
ioxynil-methyl	—	Neat	100mg	MET-397A
iprobefos	26087-47-8	Neat	500mg	PS-2065
iprodone	36734-19-7	Neat	1g	PS-1052
isazophos	42509-80-8	Neat	1g	PS-681
isobenzan	297-78-9	Neat	10mg	PS-704
isobornyl thiocyanacetate	115-31-1	Neat	1g	PS-731
isobutryldehyde-2,4-DNPH	2057-82-1	ACN	100	33084
isobutyl alcohol	78-83-1	PTM	2,000	30293
isodrin	465-73-6	Neat	1g	PS-2103
isofenphos	25311-71-1	Neat	100mg	PS-1003
isofenphos-des-N-isopropyl	—	Neat	100mg	MET-1003C
isofenphos-des-N-isopropyl	—	H	100	MET-1003A
isofenphos-des-N-isopropyl- O-analogue (5mL)	—	Cy	10ng/µL	MET-1003B

free data

Available on Our Website: Lot Certificates, Data Packs, and MSDSs

For complete information detailing manufacturing and testing for Restek inventoried reference standards, visit our website at www.restek.com.

To view lot certificates and/or an MSDS, enter the catalog number of the product in the Search feature. For a free data pack (Adobe® PDF file), enter the catalog number and lot number of the product.

Individual Compounds: Environmental, Petrochemical, Other

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	cat.#	Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	cat.#
isomalathion	3344-12-5	Neat	5ml	MET-86A	mecoprop-2,4,4-trimethyl-pentyl ester (5mL)	—	A	100	PS-1212
isoprocab	2631-40-5	Neat	500mg	PS-2070	mefenac	73250-68-7	Neat	500mg	PS-2140
isopropalin	33820-53-0	Neat	1g	PS-406	mefenoxam	70630-17-0	Neat	100mg	PS-2160
isopropylbenzene	98-82-8	Neat	1g	PS-2179	mefluide	53780-34-0	Neat	100mg	PS-2050
isopropylbenzene	98-82-8	PTM	2,000	30294	mepiridine	50-13-5	PTM	1,000	34004
isopropyl-4,4'-dichloro-benzilate	5836-10-2	Neat	1g	PS-857	mephobarbital	115-38-8	PTM	1,000	34034
2-isopropyl-6-methyl-4-pyrimidinol	2814-20-2	Neat	1g	MET-90B	mephosfolan	950-10-7	Neat	100mg	PS-1026
isoprothiolane	50512-35-1	Neat	250mg	PS-2101	mepiquat chloride	24307-26-4	Neat	10mg	PS-1047
isoproturon	34123-59-6	Neat	1g	PS-2000	meprobamate	57-53-4	PTM	1,000	34059
isovaleraldehyde-2,4-DNPH	2256-01-1	ACN	100	33085	mepronil	55814-41-0	Neat	100mg	PS-2227
2-isovaleryl-1,3-indanedione	83-28-3	Neat	250mg	PS-911	metalaxyl	57837-19-1	Neat	100mg	PS-1099
isoxaben	82558-50-7	Neat	1g	PS-2129	metalddehyde	9002-91-9	Neat	1g	PS-2030
isoxaflutole	141112-29-0	Neat	500mg	PS-2166	metamitron	41394-05-2	Neat	500mg	PS-2117
ivermectin	70288-86-7	Neat	100mg	PS-2220	Metasystox I®	919-86-8	Neat	50mg	PS-1096
jet fuel A (5mL)	64742-47-8	D	50,000	31243	metazachlor	—	Neat	100mg	PS-2312
jet fuel A	64742-47-8	D	5,000	31215	methabenzthiazuron	18691-97-9	Neat	100mg	PS-2048
jet fuel A	64742-47-8	D	50,000	31242	methacrifos	62610-77-9	Neat	100mg	PS-2218
JP-4 military fuel	8008-20-6	D	5,000	31219	methacrolein-2,4-DNPH	5077-73-6	ACN	100	33095
JP-4 military fuel	8008-20-6	D	50,000	31250	methacrylonitrile	126-98-7	PTM	2,000	30297
JP-4 military fuel	8008-20-6	PTM	50,000	30472	methadone	1095-90-5	PTM	1,000	34005
JP-5 military fuel (5mL)	8008-20-6	D	50,000	31253	metham sodium	137-42-8	Neat	500mg	PS-221
JP-5 military fuel	8008-20-6	D	5,000	31220	methamidophos	10265-92-6	Neat	100mg	PS-676
JP-5 military fuel	8008-20-6	D	50,000	31252	(+)-methamphetamine	51-57-0	PTM	1,000	34021
JP-8 military fuel	8008-20-6	D	5,000	31262	methanol	67-56-1	DMSO	15mg/mL	36401
JP-8 military fuel	8008-20-6	D	50,000	31254	methanol	67-56-1	W	10,000	30467
Kelthane	115-32-2	Neat	250mg	PS-82	methaqualone	340-56-7	PTM	1,000	34064
Kepon	143-50-0	Neat	50mg	PS-701	methidathion	950-37-8	Neat	1g	PS-679
kerosene composite (5mL)	84742-81-0	D	50,000	31257	methiocarb	2032-65-7	Neat	1g	PS-543
kerosene composite	84742-81-0	D	5,000	31094	methiocarb sulfone	—	Neat	50mg	MET-543A
kerosene composite	84742-81-0	D	50,000	31256	methiocarb sulfoxide (5mL)	—	T	100	MET-543B
kerosene: 25% weathered	84742-81-0	D	5,000	31230	methohexital	151-83-7	PTM	1,000	34035
kerosene: 50% weathered	84742-81-0	D	5,000	31231	methomyl	16752-77-5	Neat	1g	PS-775
kerosene: 75% weathered	84742-81-0	D	5,000	31232	Methoprene	40596-69-8	Neat	250mg	PS-1040
kerosene: unweathered	84742-81-0	D	5,000	31229	methoprotrotryne	841-06-5	Neat	100mg	PS-2122
kinoprene	42588-37-4	Neat	100mg	PS-2051	methoxychlor	72-43-5	M	1,000	32233
kresoxim-methyl	143390-89-0	Neat	250mg	PS-2182	methoxychlor	72-43-5	Neat	1g	PS-83
lactofen	77501-63-4	Neat	100mg	PS-2077	<i>o,p'</i> -methoxychlor	—	Neat	100mg	MET-83A
LAMPA	40158-98-3	ACN	1,000	34075	<i>p,p'</i> -methoxychlor-olefin	—	Neat	50mg	MET-83B
lead arsenate	7784-40-9	Neat	1g	PS-290	2-methoxyethanol	109-86-4	DMSO	0.25mg/mL	36402
lenacil	2164-08-1	Neat	1g	PS-415	1-(methylamino)anthraquinone	82-38-3	D	100	31823
leptophos	21609-90-5	Neat	1g	PS-677	methyl anthranilate	134-20-3	Neat	100mg	PS-2055
leptophos oxon	—	Neat	50mg	MET-677A	methyl arachidate	1120-28-1	Neat	35056	35056
levorphanol	5985-38-6	PTM	1,000	34003	methyl arachidonate	2566-89-4	Neat	35060	35060
lindane	58-89-9	Neat	1g	PS-71	methyl behenate	929-77-1	Neat	35062	35062
linuron	330-55-2	Neat	250mg	PS-372	methyl 2-bromopropionate	5445-17-0	MTBE	1,000	31654
Lontrel®	1702-17-6	Neat	250mg	PS-1069	methyl <i>tert</i> -butyl ether (MTBE)	1634-04-4	PTM	2,000	30402
lorazepam	846-49-1	PTM	1,000	34051	methylbutylketone	591-78-6	DMSO	0.25mg/mL	36400
LSD	50-37-3	ACN	25	34089	methyl caprate	110-42-9	Neat	35041	35041
LSD	50-37-3	ACN	100	34088	methyl caproate	106-70-7	Neat	35037	35037
malathion	121-75-5	Neat	500mg	PS-86	methyl caprylate	111-11-5	Neat	35039	35039
malathion monocarboxylic acid	35884-76-5	Neat	100mg	MET-86D	methylcyclohexane	108-87-2	DMSO	5.9mg/mL	36403
malathion monocarboxylic acid	35884-76-5	MTBE	100	MET-86B	methyl-2,3-dibromopropionate	1729-67-5	MTBE	1,000	31656
malathion-O-analog	1634-78-2	Neat	100mg	MET-86C	2-methyl-4,6-dinitrophenol	534-52-1	M	1,000	31292
maleic hydrazide	123-33-1	Neat	1g	PS-65	methyl eicosadienoate	2463-02-7	Neat	35058	35058
mancozeb	8018-01-7	Neat	500mg	PS-2095	methyl eicosatrienoate	55682-88-7	Neat	35059	35059
maneb	12427-38-2	Neat	1g	PS-22	methyl eicosenoate	2390-09-2	Neat	35057	35057
MCPA	97-74-6	M	1,000	32269	α -methylene- γ -butyrolactone (AMGBL)	547-65-9	ACN	1,000	34079
MCPB	94-81-5	Neat	100mg	PS-307	methylene chloride (dichloromethane)	75-09-2	DMSO	3mg/mL	36404
MCPB methyl ester	57153-18-1	Neat	500mg	PS-1084	methylene chloride (dichloromethane)	75-09-2	PTM	2,000	30401
MCPP	93-65-2	M	1,000	32271	methyl erucate	1120-34-9	Neat	35063	35063
3,4-MDA HCl	4764-17-4	M	1,000	34070	methyl heineicosanoate	6064-90-0	Neat	35061	35061
3,4-MDEA HCl	82801-81-8	M	1,000	34072	methyl heptadecanoate	1731-92-6	Neat	35050	35050
3,4-MDMA HCl	42542-10-9	M	1,000	34071	methyl heptanoate	106-73-0	Neat	35038	35038
4,4'-MDIP	72375-24-7	DMSO	1,000	33003	methyl isothiocyanate	556-61-6	Neat	1g	MET-221A
mecoprop	7085-19-0	Neat	250mg	PS-324					
mecoprop methyl ester	23844-56-6	Neat	100mg	PS-1106					
mecoprop-P	16484-77-8	Neat	500mg	PS-2116					

*Volume is 1mL/ampul unless otherwise noted. Concentration is $\mu\text{g/mL}$ unless otherwise noted.

Solvent code:

A=acetone
ACN=acetonitrile
C=carbon disulfide
Cy=cyclohexane

D=methylene chloride
DMSO=dimethyl sulfoxide
EA=ethyl acetate
H=hexane

I=isooctane
Ip=isopropanol
M=methanol
MTBE=methyl *tert*-butyl ether

PTM=purge & trap grade methanol
T=toluene
TO=transformer oil
W=water (DI)

Cat.#s starting with PS- & MET- are ChemService pesticides that are subject to size and price changes without notification. They are not available in Europe.

Compound	Solvent	Individual	Compound	Solvent	Individual				
Packaged 1mL/ampul*	CAS#	Code	μg/mL*	cat.#	Packaged 1mL/ampul*	CAS#	Code	μg/mL*	cat.#
methyl laurate	111-82-0	Neat		35043	morphine	6211-15-0	PTM	1,000	34006
methyl lignocerate	2442-49-1	Neat		35064	motor oil composite	64742-47-8	D	50,000	31464
methyl linoleate	112-63-0	Neat		35053	nabam	142-59-6	Neat	1g	PS-220
methyl linolenate	301-00-8	Neat		35054	naled	300-76-5	Neat	250mg	PS-605
methyl methacrylate	80-62-6	PTM	2,000	30299	naphthalene	91-20-3	M	1,000	31280
N-methyl-methamidophos	28167-49-9	Neat	100mg	PS-676-1	naphthalene-d8	1146-65-2	D	2,000	31043
methyl myristate	124-10-7	Neat		35045	1-naphthalene acetamide	86-86-2	Neat	1g	PS-28
methyl myristoleate	56219-06-8	Neat		35046	1-naphthaleneacetic acid	86-87-3	Neat	1g	PS-29
1-methylnaphthalene	90-12-0	M	1,000	31283	1,8-naphthalic anhydride	81-84-5	Neat	1g	PS-1048
2-methylnaphthalene	91-57-6	D	1,000	31285	2-naphthoxyacetic acid	120-23-0	Neat	1g	PS-35
methyl-1-naphthalene acetate	2876-78-0	Neat	1g	PS-30	N-1-naphthylacetamide	—	Neat	500mg	PS-354
N-methyl-N-1-naphthyl acetamide	—	Neat	1g	PS-352	N-1-naphthylphthalamic acid	132-66-1	Neat	250mg	PS-343
methyl nervonate	2733-88-2	Neat		35065	N-(1-naphthyl)phthalimide	—	Neat	1g	MET-343A
2-methyl-4-nitroaniline	99-55-8	M	1,000	31612	α-naphthylthiourea	86-88-4	Neat	1g	PS-106
3-methyl-4-nitrophenol	25081-34-2	Neat	1g	MET-678B	neburon	555-37-3	Neat	1g	PS-377
methyl nonadecanoate	1731-94-8	Neat		35055	Niclofen	1836-75-5	Neat	100mg	PS-394
methyl nonanoate	1731-84-6	Neat		35040	niclosamide	50-65-7	Neat	250mg	PS-1207
methyl oleate	112-62-9	Neat		35052	nicotine	54-11-5	M	1,000	34085
methyl palmitate	112-39-0	Neat		35048	nicotine	54-11-5	Neat	1g	PS-85
methyl palmitoleate	1120-25-8	Neat		35049	nitralin	4726-14-1	Neat	1g	PS-358
methyl parathion	298-00-0	Neat	250mg	PS-94	nitrapyrin	1929-82-4	Neat	100mg	PS-419
methyl pentadecanoate	7162-64-1	Neat		35047	nitrazepam	146-22-5	PTM	1,000	34053
4-methyl-2-pentanone (MIBK)	108-10-1	PTM	5,000	30400	nitrobenzene	99-95-3	ACN	1,000	31657
3-methylphenol	108-39-4	M	1,000	31403	nitrobenzene-d5	4165-60-0	D	2,000	31044
N-methylpiperidine	626-67-5	Neat	1g	MET-1047A	nitrobenzene-d5	4165-60-0	D	2,000	33904
N-methylpyrrolidone	872-50-4	DMSO	2.65mg/mL	36405	nitroglycerin	55-63-0	M	1,000	31498
methyl stearate	112-61-8	Neat		35051	nitroguanidine	556-88-7	M	1,000	31602
2-methylsulfonyl-4-tri-fluoromethyl benzoic acid	—	Neat	100mg	MET-2166B	2-nitromesitylene	603-71-4	M	2,000	33902
methyl tridecanoate	1731-88-0	Neat		35044	nitromethane	75-52-5	DMSO	0.25mg/mL	36406
methyl α-trifluoromethyl-benzoate	344-96-7	Neat	100mg	MET-2057A	4-nitrophenol	100-02-7	M	1,000	31296
methyl undecanoate	1731-86-8	Neat		35042	1-nitropropane	108-03-2	D	250	35263
methylpyrrolon	125-64-4	PTM	1,000	34060	2-nitropropane	79-46-9	PTM	2,000	30403
metobromuron	3060-89-7	Neat	1g	PS-374	N-nitrosodimethylamine	62-75-9	M	1,000	31427
metolachlor	51218-45-2	Neat	500mg	PS-403	N-nitrosodimethylamine-d6	17829-05-9	D	1,000	33910
S-metolachlor	87392-12-9	Neat	250mg	PS-403-1	N-nitrosodiphenylamine	86-30-6	M	1,000	31429
metolcarb	1129-41-5	Neat	500mg	PS-2096	N-nitroso-di-n-propylamine	621-64-7	M	1,000	31428
metoxuron	19937-59-8	Neat	250mg	PS-2178	N-nitrosodi-n-propylamine-d14	93951-96-3	D	1,000	33911
metribuzin	21087-64-9	A	1,000	32436	N-nitrosopendimethalin	—	Neat	100mg	MET-401A
metribuzin	21087-64-9	Neat	1g	PS-398	2-nitrotoluene	88-72-2	ACN	1,000	31659
metsulfuron methyl	74223-64-6	Neat	100mg	PS-1078	3-nitrotoluene	99-08-1	ACN	1,000	31660
cis-mevinphos	26718-65-0	Neat	100mg	PS-87-1	4-nitrotoluene	99-99-0	ACN	1,000	31658
MGK® 264	113-48-4	Neat	250mg	PS-962	cis-nonachlor	5103-73-1	Neat	25mg	MET-75B
MGK® R-11	126-15-8	Neat	1g	PS-910	trans-nonachlor	39765-80-5	Neat	25mg	MET-75C
mineral spirits:					nonatriacontane (C39) (10mL)	7194-86-7	C	3,000	31877
25% weathered	8030-30-6	D	5,000	31226	n-nonatriacontane (C39)	7194-86-7	C	3,000	31456
mineral spirits:					Norflurazon	27314-13-2	Neat	100mg	PS-1044
50% weathered	8030-30-6	D	5,000	31227	n-octacosane (C28)	630-02-4	D	1,000	31672
mineral spirits:					octamethylpyrophosphoramidate	—	Neat	50mg	PS-2338
75% weathered	8030-30-6	D	5,000	31228	omethoate	1113-02-6	Neat	100mg	PS-2017
mineral spirits: unweathered (5mL)	8030-30-6	D	50,000	31261	Omite®	2312-35-8	Neat	1g	PS-858
mineral spirits: unweathered	8030-30-6	D	5,000	31225	ormetoprin	6981-18-6	Neat	250mg	PS-2277
mineral spirits: unweathered	8030-30-6	D	50,000	31260	oryzalin	19044-88-3	Neat	1g	PS-410
Mirex	2385-85-5	Neat	100mg	PS-702	ovex	80-33-1	Neat	1g	PS-855
molinate	2212-67-1	Neat	250mg	PS-501	oxadiazon	19666-30-9	Neat	1g	PS-1050
Monceren®	66063-05-6	Neat	100mg	PS-1095	oxadiazon-hydroxy (5mL)	—	T	10ng/μL	MET-1050A
monocrotophos	6923-22-4	Neat	1g	PS-609	oxadixyl	77732-09-3	Neat	100mg	PS-2137
monolein (1-mono[cis-9-octadecenoyl]-rac-glycerol)	111-03-5	pyridine	5,000	33021	oxamyl	23135-22-0	Neat	250mg	PS-737
monolinuron	1746-81-2	Neat	250mg	PS-2210	oxazepam	604-75-1	PTM	1,000	34054
monomethyl tetrachloro-terephthalate	—	Neat	100mg	MET-33A	oxycarboxin	5259-88-1	Neat	1g	PS-295
monopalmitin	524-44-9	pyridine	5,000	33026	oxychlorane (5mL)	27304-13-8	M	100	MET-75A
monosodium acid methane arsonate sesquihydrate	2163-80-6	Neat	1g	PS-429	oxycodone	124-90-3	PTM	1,000	34007
monothiono TEPP	645-78-3	Neat	50mg	PS-2181	oxymeteton-methyl	301-12-2	Neat	50mg	PS-641
monuron	150-68-5	Neat	1g	PS-59	oxyfluorfen	42874-03-3	Neat	1g	PS-1033
monuron TCA	140-41-0	Neat	1g	PS-371	oxymorphone	76-41-5	PTM	1,000	34065
					oxytetracycline hydrochloride	2058-46-0	Neat	250mg	PS-2200
					paclobutrazol	76738-62-0	Neat	500mg	PS-2097
					paraoxon	311-45-5	Neat	100mg	PS-610
					paraquat Cl tetrahydrate	1910-42-5	Neat	1g	PS-366
					Parathion	56-38-2	Neat	500mg	PS-95
					penconazole	66246-88-6	Neat	100mg	PS-2207
					Pentac®	2227-17-0	Neat	1g	PS-862

free data

Available on Our Website: Lot Certificates, Data Packs, and MSDSs

For complete information detailing manufacturing and testing for Restek inventoried reference standards, visit our website at www.restek.com.

To view lot certificates and/or an MSDS, enter the catalog number of the product in the Search feature. For a free data pack (Adobe® PDF file), enter the catalog number and lot number of the product.

Individual Compounds: Environmental, Petrochemical, Other

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	cat.#	Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$	cat.#
pentachloroaniline	527-20-8	Neat	100mg	MET-150A	phenyl-2-thiourea	103-85-5	Neat	1g	PS-801
pentachloroanisole	1825-21-4	M	1,000	32268	phorate	298-02-2	Neat	100mg	PS-654
pentachloroanisole	1825-21-4	Neat	1g	PS-1109	phorate-oxon (5mL)	—	T	100	MET-654A
2,2',4,5,5'-pentachlorobiphenyl (BZ #101)	37680-73-2	I	10	32285	phorate sulfone (5mL)	—	H	100	MET-654B
2,3',4,4',5'-pentachlorobiphenyl (BZ #118)	31508-00-6	I	10	32293	phorate sulfoxide	2588-03-6	Neat	100mg	MET-654C
pentachloroethane	76-01-7	PTM	2,000	30404	phosalone	2310-17-0	Neat	1g	PS-682
pentachloronitrobenzene	82-68-8	EA	100	32091	Phosdrin	7786-34-7	Neat	250mg	PS-87
pentachloronitrobenzene	82-68-8	Neat	1g	PS-150	phosfolan	947-02-4	Neat	1g	PS-842
pentachlorophenol	87-86-5	M	1,000	31297	phosmet-oxon (5mL)	—	T	100	MET-653A
pentachlorophenol	87-86-5	Neat	1g	PS-7	phosphamidon	13171-21-6	Neat	250mg	PS-606
pentachlorothiophenol	133-49-3	Neat	100mg	MET-690B	phosphan	115-78-6	Neat	50mg	PS-67
<i>n</i> -pentacontane (C50)	6596-40-3	T	10	31685	phoxim	14816-18-3	Neat	1g	PS-2100
pentacosane (C25)	629-99-2	D	10,000	31487	phthalamic acid	—	Neat	100mg	MET-26A
pentafluorobenzene	363-72-4	PTM	2,000	30031	phthalanilic acid	4727-29-1	Neat	1g	PS-2114
pentafluorophenol	771-61-9	D	2,000	31048	picloram	1918-02-1	M	1,000	32265
2-pentanone	107-87-9	D	250	35261	picloram	1918-02-1	Neat	250mg	PS-274
pentazocine	64024-15-3	PTM	1,000	34062	picloram methyl ester	N/A	M	1,000	32266
pentobarbital	76-74-4	PTM	1,000	34036	picloram methyl ester	14143-55-6	Neat	100mg	PS-1108
perchloromethyl mercaptan	594-42-3	Neat	1g	PS-224	picric acid	88-89-1	M	1,000	31499
perfluorotributylamine (PFTBA)	311-89-7	Neat	1mL	30482	Pinnacle	79277-27-3	Neat	100mg	PS-2011
perfluorotributylamine (PFTBA)	311-89-7	Neat	1g	33027	pinoxaden	243973-20-8	Neat	100mg	PS-2319
permethrin	52645-53-1	Neat	250mg	PS-758	piperalin	3478-94-2	Neat	1g	PS-409
<i>cis</i> -permethrin	—	Neat	50mg	PS-758-1	piperonyl butoxide	51-03-6	Neat	100mg	PS-100
<i>trans</i> -permethrin (isomeric mix)	61949-77-7	Neat	50mg	PS-758-2	pirimicarb	23103-98-2	Neat	1g	PS-757
Perthane	72-56-0	Neat	1g	PS-73	pirimiphos-ethyl	23505-41-1	Neat	250mg	PS-643
PETN (pentaerythritol tetranitrate)	78-11-5	M	1,000	31600	pirimiphos-methyl	29232-93-7	Neat	250mg	PS-644
phenamiphos	22224-92-6	Neat	1g	PS-612	2-pivaloyl-1,3-indandione	83-26-1	Neat	250mg	PS-905
phenanthrene	85-01-8	M	1,000	31281	Polywax® 500	9006-42-2	Neat	1g	PS-225
phenanthrene-d10	1517-22-2	D	2,000	31045	Polywax® 655	9002-88-4	Neat	1g	36224
phencyclidine	956-90-1	PTM	1,000	34027	Polywax® 850	9002-88-4	Neat	1g	36225
phendimetrazine	50-58-8	PTM	1,000	34025	Polywax® 1000	9002-88-4	Neat	1g	36226
phenmedipham	13684-63-4	Neat	250mg	PS-1014	potassium cyanate	590-28-3	Neat	1g	PS-289
Phenmeryl nitrate	55-68-5	Neat	100mg	PS-202	prallethrin	23031-36-9	Neat	10mg	PS-2109
phenmetrazine	1707-14-8	PTM	1,000	34026	prazepam	2955-38-6	PTM	1,000	34055
phenobarbital	50-06-6	PTM	1,000	34037	primisulfuron-methyl	86209-51-0	Neat	500mg	PS-2169
phenol	108-95-2	M	1,000	31298	probenazole	27605-76-1	Neat	500mg	PS-2082
phenol-d6	13127-88-3	D	2,000	31049	prochloraz	67747-09-5	Neat	250mg	PS-2175
phenothiazine	92-84-2	Neat	1g	PS-733	procyclazine	32889-48-8	Neat	1g	PS-402
<i>d</i> -(<i>cis-trans</i>)-phenothrin	26002-80-2	Neat	50mg	PS-1030	procymidone	32809-16-8	Neat	250mg	PS-2126
phenoxyacetic acid	122-59-8	Neat	1g	PS-36	prodiamine	29091-21-2	Neat	500mg	PS-2133
3-phenoxybenzoic acid	3739-38-6	Neat	1g	MET-758C	profenofos	41198-08-7	Neat	250mg	PS-1024
<i>m</i> -phenoxybenzyl alcohol	13826-35-2	Neat	1g	MET-758B	Profluralin	26399-36-0	Neat	1g	PS-399
2-phenoxypropionic acid	940-31-8	Neat	1g	PS-37	prohexadione	88805-35-0	Neat	100mg	PS-2226-1
phentermine	1197-21-3	PTM	1,000	34024	prohexadione-calcium	127277-53-6	Neat	100mg	PS-2226
phenthoate	2597-03-7	Neat	100mg	PS-2025	promecarb	2631-37-0	Neat	100mg	PS-2027
phentioacetate	900-95-8	Neat	1g	PS-1021	prometon	1610-18-0	Neat	250mg	PS-386
phenyl mercuric acetate	62-38-4	Neat	1g	PS-14	prometryne	7287-19-6	A	1,000	32449
phenyl mercuric benzoate	—	Neat	1g	PS-15	prometryne	7287-19-6	Neat	1g	PS-384
phenyl mercuric chloride	100-56-1	Neat	1g	PS-214	propachlor	1918-16-7	M	1,000	32235
phenyl mercuric hydroxide	100-57-2	Neat	1g	PS-215	propachlor	1918-16-7	Neat	500mg	PS-355
phenyl mercuric lactate	—	Neat	1g	PS-216	propamocarb HCl	25606-41-1	Neat	250mg	PS-2167
phenyl mercuric oleate	104-60-9	Neat	100mg	PS-199	2-propanol	67-63-0	W	50,000	30473
phenyl mercuric phthalate	—	Neat	1g	PS-217	propazine	139-40-2	A	1,000	32448
phenylmercury borate	102-98-7	Neat	100mg	PS-203	propazine	139-40-2	Neat	500mg	PS-385
phenylmercury salicylate	28086-13-7	Neat	1g	PS-205	propazine-2-hydroxy (5mL)	—	ACN:W (75:25)	100	MET-385A
5-(phenylmethyl)- 3-furanmethanol	—	Neat	100mg	MET-1000A	propham	122-42-9	Neat	1g	PS-53
<i>o</i> -phenylphenol	90-43-7	Neat	1g	PS-420	prophos	13194-48-4	Neat	1g	PS-672
phenylpropanolamine HCl	154-41-6	M	1,000	34073	propineb	9016-72-2	Neat	250mg	PS-2062
					propionaldehyde-2,4-DNPH	725-00-8	ACN	100	33086
					propionitrile	107-12-0	PTM	2,000	30407
					<i>n</i> -propylbenzene	103-65-1	PTM	2,000	30408

*Volume is 1mL/ampul unless otherwise noted. Concentration is $\mu\text{g/mL}$ unless otherwise noted.

Solvent code:

A = acetone	D = methylene chloride	I = isooctane	PTM = purge & trap grade methanol
ACN = acetonitrile	DMSO = dimethyl sulfoxide	Ip = isopropanol	T = toluene
C = carbon disulfide	EA = ethyl acetate	M = methanol	TO = transformer oil
Cy = cyclohexane	H = hexane	MTBE = methyl <i>tert</i> -butyl ether	W = water (DI)

free data

Available on Our Website: Lot Certificates, Data Packs, and MSDSs

For complete information detailing manufacturing and testing for Restek inventoried reference standards, visit our website at www.restek.com.

To view lot certificates and/or an MSDS, enter the catalog number of the product in the Search feature. For a free data pack (Adobe® PDF file), enter the catalog number and lot number of the product.

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	Individual cat.#
S-propyl butylethylthio- carbamate	1114-71-2	Neat	1g	PS-55
propylene glycol dinitrate (PGDN)	6423-43-4	M	1,000	31821
propylzamide	23950-58-5	Neat	1g	PS-349
prosulfluron	94125-34-5	Neat	500mg	PS-2123
Prowl	40487-42-1	Neat	1g	PS-401
pymetrozine	123312-89-0	Neat	250mg	PS-2189
pyraclofos	77458-01-6	Neat	500mg	PS-2046
pyraflufen-ethyl	129630-17-7	Neat	100mg	PS-2306
pyrazon	1698-60-8	Neat	1g	PS-430
pyrazophos	13457-18-6	Neat	100mg	PS-2049
pyrazosulfuron-ethyl	93697-74-6	Neat	100mg	PS-2105
pyrazoxyfen	71561-11-0	Neat	1g	PS-2064
pyrene	129-00-0	M	1,000	31282
pyrethrum	8003-34-7	Neat	100mg	PS-97
pyridaben	96489-71-3	Neat	1g	PS-2125
pyridaphenthion	119-12-0	Neat	250mg	PS-2155
pyridate	55512-33-9	Neat	500mg	PS-2081
pyridine	85404-20-2	DMSO	1mg/mL	36407
pyridine	110-86-1	PTM	2,000	30409
pyridine-d5	7291-22-7	D	2,000	31046
pyriproxyfen	95737-68-1	Neat	100mg	PS-2110
pyrithiobac-sodium	123343-16-8	Neat	500mg	PS-2145
quinalphos	13593-03-8	Neat	250mg	PS-2010
quinclorac	84087-01-4	Neat	500mg	PS-2098
quinomethionate	2439-01-2	Neat	250mg	PS-860
quinoxifen	124495-18-7	Neat	100mg	PS-2284
quizalofop	76578-12-6	Neat	50mg	MET-1080A
quizalofop ethyl	76578-14-8	Neat	100mg	PS-1080
quizalofop-P-ethyl	100646-51-3	Neat	500mg	PS-2151
quizalofop methyl	—	Neat	500mg	MET-1080B
RDX	121-82-4	ACN	1,000	31666
resmethrin	10453-86-8	Neat	250mg	PS-1000
rimsulfuron	122931-48-0	Neat	100mg	PS-2118
rotenone	83-79-4	Neat	250mg	PS-99
Ruelene	299-86-5	Neat	50mg	PS-607
Safroin®	31218-83-4	Neat	100mg	PS-1076
secobarbital	29071-21-4	PTM	1,000	34038
Sesone	149-26-8	Neat	1g	PS-313
sethoxymim	74051-80-2	Neat	100mg	PS-2013
siduron	1982-49-6	Neat	1g	PS-373
Silvex	93-72-1	Neat	250mg	PS-47
Silvex isooctyl ester	—	Neat	1g	PS-329
Silvex methyl ester	—	Neat	250mg	PS-342
simazine	122-34-9	A	1,000	32236
simazine	122-34-9	Neat	100mg	PS-58
simazine hydroxy	—	Neat	25mg	MET-58B
simetryn	1014-70-6	Neat	50mg	PS-381
sodium arsenate	7778-43-0	Neat	1g	PS-428
sodium fluoroacetate	62-74-8	Neat	1g	PS-105
Sonar®	59756-60-4	Neat	100mg	PS-1070
spinosad	—	Neat	100mg	PS-2211
spiromesifen	283594-90-1	Neat	100mg	PS-2322
STB	132827-25-9	Neat	500mg	MET-222B
stearyl stearate (10mL)	2778-96-3	Cy	2,000	31636
stearyl stearate (10mL)	2778-96-3	H	2,000	31681
stearyl stearate	2778-96-3	Neat	100mg	31860
Stoddard solvent	8052-41-3	PTM	10,000	30487
Strobane	8001-50-1	Neat	100mg	PS-2023
styrene	100-42-5	PTM	2,000	30410
sulfadimethoxine	122-11-2	Neat	250mg	PS-2278
sulfaquinoxaline	59-40-5	Neat	10mg	PS-2056
sulfentrazone	122836-35-5	Neat	250mg	PS-2132
sulfolane	126-33-0	DMSO	0.8mg/mL	36413
sulfometuron methyl	74222-97-2	Neat	100mg	PS-1074
sulfosulfuron	141776-32-1	Neat	100mg	PS-2224
sulfotep	3689-24-5	Neat	50mg	PS-2024
sulfoxide	120-62-7	Neat	100mg	PS-976
sulprofos	35400-43-2	Neat	100mg	PS-1018
sulprofos sulfone (5mL)	—	T	100	MET-1018A
sulprofos-sulfoxide	34643-47-5	Neat	100mg	MET-1018C
super-rozol	28772-56-7	Neat	100mg	PS-1038
swep	1918-18-9	Neat	50mg	PS-1098
Sythane®	88671-89-0	Neat	100mg	PS-2006

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual µg/mL*	Individual cat.#
2,4,5-T	93-76-5	M	200	32243
2,4,5-T	93-76-5	Neat	250mg	PS-45
2,4,5-Tamine	—	Neat	1g	PS-296
2,4,5-T butoxyethyl ester	2545-59-7	Neat	1g	PS-309
2,4,5-T <i>n</i> -butyl ester	93-79-8	Neat	1g	PS-303
2,4,5-T 2-ethylhexyl ester	—	Neat	250mg	PS-297-1
2,4,5-T isooctyl ester	25168-15-4	Neat	100mg	PS-297
2,4,5-T isopropyl ester	—	Neat	100mg	PS-46
2,4,5-T methyl ester	1928-37-6	M	1,000	32244
talbutal	115-44-6	PTM	1,000	34039
tau-fluvalinate	102851-06-9	Neat	100mg	PS-1071
TCMTB	21564-17-0	Neat	100mg	PS-2021
2,4-TDIP	72375-21-4	DMSO	1,000	33001
2,6-TDIP	195625-39-9	DMSO	1,000	33000
tebufenozide	112410-23-8	Neat	100mg	PS-2188
tebupirimfos	96182-53-5	Neat	100mg	PS-2192
tebuthiuron	34014-18-1	Neat	1g	PS-407
tefluthrin	79538-32-2	ACN	100	PS-2238
temazepam	896-50-4	PTM	1,000	34056
terbacil	5902-51-2	Neat	1g	PS-375
terbufos	13071-79-9	Neat	500mg	PS-680
terbumeton	33693-04-8	Neat	1g	PS-2120
terbutylazine	5915-41-3	A	1,000	32447
terbutylazine	5915-41-3	Neat	50mg	PS-413
terbutylazine-desethyl	—	Neat	50mg	MET-413A
terbutol	1918-11-2	Neat	100mg	PS-550
terbutryne	886-50-0	Neat	1g	PS-379
<i>p</i> -terphenyl-d14	1718-51-0	D	1,000	31828
<i>o</i> -terphenyl	84-15-1	A	2,000	31066
<i>o</i> -terphenyl	84-15-1	D	10,000	31097
<i>p</i> -terphenyl	92-94-4	D	10,000	31095
α -terpineol	98-55-5	D	2,000	33912
Terrazole®	2593-15-9	Neat	1g	PS-264
2,2',5,5'-tetrachlorobiphenyl (BZ #52)	35693-99-3	I	10	32284
1,1,1,2-tetrachloroethane	630-20-6	PTM	2,000	30411
1,1,2,2-tetrachloroethane	79-34-5	PTM	2,000	30412
tetrachloroethene	127-18-4	PTM	2,000	30413
tetrachlorohydroquinone	87-87-6	Neat	1g	MET-690A
1,2,4,5-tetrachloro- 3-nitrobenzene	117-18-0	Neat	1g	PS-151
2,3,4,6-tetrachlorophenol	58-90-2	M	1,000	31402
2,3,4,6-tetrachlorophenol	58-90-2	Neat	100mg	PS-8
tetrachloroterephthalic acid	2136-79-0	Neat	100mg	MET-33B
2,4,5,6-tetrachloro- <i>m</i> -xylene (5mL)	877-09-8	A	200	32028
2,4,5,6-tetrachloro- <i>m</i> -xylene	877-09-8	A	200	32027
tetrachlorvinphos	961-11-5	Neat	1g	PS-611
<i>n</i> -tetracontane (C40)	4181-95-7	Neat	100mg	31859
tetradifon	116-29-0	Neat	1g	PS-80
tetraethylpyrophosphate	107-49-3	Neat	500mg	PS-601
tetrahydrofuran (THF)	109-99-9	DMSO	3.6mg/mL	36408
tetrahydrofuran (THF)	109-99-9	PTM	2,000	30414
<i>cis</i> -1,2,3,6-tetrahydro- phthalimide	27813-21-4	Neat	5g	MET-260A
tetralin	119-64-2	DMSO	0.5mg/mL	36409
tetramethrin	7696-12-0	Neat	100mg	PS-1042
tetramethylthiuram disulfide	137-26-8	Neat	1g	PS-24
tetrapentyltin	3765-65-9	D	2,000	31475
tetra- <i>n</i> -propyltin	2176-98-9	D	2,000	31474
tefryl	479-45-8	ACN	1,000	31667
Δ^9 -THC	1972-08-3	M	1,000	34067
(\pm)11-nor-9-carboxy- Δ^9 -THC	104874-50-2	M	100	34068
thallium sulfate	7446-18-6	Neat	1g	PS-804
thebaine	115-37-7	PTM	1,000	34009
thiabendazole	148-79-8	Neat	1g	PS-1057
thiamethoxam	153719-23-4	Neat	100mg	PS-2215
thiamylal	337-47-3	PTM	1,000	34040
thiazopyr	117718-60-2	Neat	500mg	PS-2163
thiazuron	51707-55-2	Neat	1g	PS-2044
thiodicarb	59669-26-0	Neat	250mg	PS-2012
thionazin-O-analog	—	Neat	100mg	MET-668A
thiophental	71-73-8	PTM	1,000	34041
thiophanate	23564-06-9	Neat	1g	PS-223
thiophanate-methyl	23564-05-8	Neat	1g	PS-230

Cat.#s starting with PS- & MET- are ChemService pesticides that are subject to size and price changes without notification. They are not available in Europe.

HROMalytic Chromatography Products '08
Australian Distributors ECHnology
 www.chromtech.net.au E-mail : info@chromtech.net.au Tel : +61 3 9762 2034 Fax : +61 3 9761 1169



Individual Compounds: Environmental, Petrochemical, Other

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$ cat.#
Tilt®	60207-90-1	Neat	250mg PS-1075
Tinuvin® P	2440-22-4	Ip	51.8 31629
Tokuthion	34643-46-4	Neat	100mg PS-1094
tolclofos-methyl	57018-04-9	Neat	500mg PS-2156
tolfenpyrad	129588-76-5	Neat	100mg PS-2323
<i>m</i> -tolualdehyde-2,4-DNPH	2880-05-9	ACN	100 33088
<i>o</i> -tolualdehyde-2,4-DNPH	1773-44-0	ACN	100 33087
<i>p</i> -tolualdehyde-2,4-DNPH	2571-00-8	ACN	100 33089
toluene	108-88-3	DMSO	4.45mg/mL 36410
toluene	108-88-3	PTM	2,000 30415
toluene-d8	2037-26-5	PTM	2,000 30224
tolylfluanid	731-27-1	Neat	250mg PS-2060
<i>N-m</i> -tolylphthalamic acid	85-72-3	Neat	1g PS-50
toxaphene	8001-35-2	H	1,000 32005
toxaphene	8001-35-2	I	5,000 32071
toxaphene	8001-35-2	M	2,000 32015
toxaphene	8001-35-2	Neat	1g PS-79
2,4,5-TP (Silvex)	93-72-1	M	1,000 32245
2,4,5-TP (Silvex) methyl ester	4841-20-7	M	1,000 32246
tralkoxydim	87820-88-0	Neat	500mg PS-2135
tralomethrin	66841-25-6	Neat	100mg PS-2115
transformer oil (PCB-free)	64742-53-6	Neat	50mL 32425
transformer oil (PCB-free)	64742-53-6	Neat	5mL 32424
<i>n</i> -triacontane-d62 (C30)	638-68-6	D	500 31816
triadimefon	43121-43-3	Neat	500mg PS-1013
triasulfuron	82097-50-5	Neat	100mg PS-2042
triazolam	28911-01-5	PTM	1,000 34057
1,2,4-triazole	288-88-0	Neat	1g MET-1075A
triazophos	24017-47-8	Neat	100mg PS-2139
triazophos (5mL)	24017-47-8	T	100 PS-1213
tribenuron methyl	101200-48-0	Neat	250mg PS-2099
2,4,6-tribromophenol	118-79-6	M	1,000 31401
tributylphosphate	126-73-8	A	1,000 32280
S,S,S-tributylphosphoro-trithioate	78-48-8	Neat	1g PS-562
tributyltin chloride	1461-22-9	D	2,000 31478
tricaprin (1,2,3-tricaprinoylglycerol)	621-71-6	pyridine	8,000 33025
tricaprin (1,2,3-tricaprinoylglycerol) (5mL)	621-71-6	pyridine	8,000 33033
trichloroacetic acid	76-03-9	Neat	1g PS-34
1,2,3-trichlorobenzene	87-61-6	PTM	2,000 30416
1,2,4-trichlorobenzene	120-82-1	M	1,000 31439
2,3,6-trichlorobenzoic acid	50-31-7	Neat	1g PS-330
2,4,4'-trichlorobiphenyl (BZ #28)	7012-37-5	I	10 32283
1,1,1-trichloroethane	74552-83-3	DMSO	50mg/mL 36411
1,1,1-trichloroethane	71-55-6	PTM	2,000 30418
1,1,2-trichloroethane	79-00-5	PTM	2,000 30419
trichloroethene	79-01-6	DMSO	0.4mg/mL 36412
trichloroethene	79-01-6	PTM	2,000 30420
trichlorofluoromethane (CFC-11)	75-69-4	PTM	2,000 30421
trichloronate	327-98-0	Neat	100mg PS-2007
trichloronitromethane	76-06-2	Neat	1g PS-4
2,4,5-trichlorophenol	95-95-4	A	1,000 32017
2,4,5-trichlorophenol	95-95-4	M	1,000 31299
2,4,5-trichlorophenol	95-95-4	Neat	1g PS-9
2,4,6-trichlorophenol	88-06-2	M	1,000 31400
2,4,6-trichlorophenol	88-06-2	Neat	1g PS-10
(2,4,5-trichlorophenoxy) acetic acid butoxypropylene glycol ester	—	Neat	1g PS-334
(2,4,5-trichlorophenoxy) acetic acid methyl ester	—	Neat	250mg PS-339
1,2,3-trichloropropane	96-18-4	MTBE	1,000 31648
1,2,3-trichloropropane	96-18-4	PTM	2,000 30429
3,5,6-trichloro-2-pyridinol	6515-38-4	Neat	250mg MET-674A
1,1,2-trichlorotrifluoroethane (CFC-113)	76-13-1	PTM	2,000 30462
triclopyr	55335-06-3	Neat	500mg PS-417
triclopyr butoxyethyl ester	64700-56-7	Neat	250mg PS-2308
triclopyr methyl ester	—	Neat	100mg PS-1110

Compound Packaged 1mL/ampul*	CAS#	Solvent Code	Individual $\mu\text{g/mL}^*$ cat.#
tricyclohexyltin hydroxide	13121-70-5	Neat	250mg PS-107
trietazine	1912-26-1	Neat	100mg PS-2107
trifloxystrobin	141517-21-7	Neat	100mg PS-2196
trifloxysulfuron-sodium	199119-58-9	Neat	100mg PS-2317
triflumizole	68694-11-1	Neat	1g PS-2121
triflurumuron	64628-44-0	Neat	100mg PS-2176
α,α,α -trifluoro-3'-hydroxy- O-toluanilide	—	Neat	250mg MET-2057C
α,α,α -trifluorotoluene	98-08-8	PTM	2,000 30048
α,α,α -trifluorotoluene	98-08-8	PTM	2,500 30068
α,α,α -trifluorotoluene	98-08-8	PTM	10,000 30083
α,α,α -trifluoro-O-toluic acid	433-97-6	Neat	500mg MET-2057B
trifluralin	1582-09-8	M	1,000 32238
trifluralin	1582-09-8	Neat	1g PS-393
triflurosulfuron-methyl	126535-15-7	Neat	500mg PS-2147
triforine	26644-46-2	Neat	100mg PS-2014
2,3,5-triiodobenzoic acid	88-82-4	Neat	1g PS-561
1,2,4-trimethylbenzene	95-63-6	PTM	2,000 30422
1,3,5-trimethylbenzene	108-67-8	PTM	2,000 30423
1,2,4-trimethyl-5-nitrobenzene	610-91-3	M	2,000 33903
2,3,5-trimethylphenyl methyl carbamate	2655-15-4	Neat	50mg PS-541
trinexapac-ethyl	95266-40-3	Neat	250mg PS-2158
1,3,5-trinitrobenzene	99-35-4	ACN	1,000 31668
2,4,6-trinitrotoluene	118-96-7	ACN	1,000 31669
tri- <i>n</i> -propyltin chloride	995-25-5	D	2,000 31476
triolein (1,2,3-Tri[<i>cis</i> -octadecenoyl] glycerol)	122-32-7	pyridine	5,000 33023
tripentyltin chloride	3342-67-4	D	2,000 31477
triphenylphosphate	115-86-6	A	1,000 32281
Tycor	64529-56-2	Neat	100mg PS-1092
2-undecanone	112-12-9	Neat	1g PS-951
uniconzole-P	83657-17-4	Neat	10mg PS-2076
unleaded gasoline composite (5mL)	8006-61-9	PTM	50,000 30206
unleaded gasoline composite	8006-61-9	PTM	2,500 30081
unleaded gasoline composite	8006-61-9	PTM	50,000 30205
unleaded gasoline: 25% weathered	8006-61-9	PTM	5,000 30097
unleaded gasoline: 50% weathered	8006-61-9	PTM	5,000 30098
unleaded gasoline: 75% weathered	8006-61-9	PTM	5,000 30099
unleaded gasoline: 99% weathered	8006-61-9	PTM	5,000 30436
unleaded gasoline: unweathered	8006-61-9	PTM	5,000 30096
used motor oil composite	64742-65-0	D	50,000 31465
vacor	53558-25-1	Neat	100mg PS-1059
valeraldehyde-2,4-DNPH	2057-84-3	ACN	100 33090
γ -valerolactone	108-29-2	ACN	1,000 34080
vamidithion	2275-23-2	Neat	100mg PS-2187
Vanicide-20S	—	Neat	1g PS-227
Velpar®	51235-04-2	Neat	1g PS-416
vernolate	1929-77-7	Neat	250mg PS-504
vinclizolin	50471-44-8	Neat	1g PS-1049
vinyl acetate	108-05-4	PTM	2,000 30216
vinyl chloride	75-01-4	PTM	2,000 30089
vinyl chloride	75-01-4	PTM	2,500 30093
warfarin	81-81-2	Neat	1g PS-104
<i>m</i> -xylene	108-38-3	DMSO	6.51mg/mL 36414
<i>m</i> -xylene	108-38-3	PTM	2,000 30424
<i>o</i> -xylene	95-47-6	DMSO	0.97mg/mL 36415
<i>o</i> -xylene	95-47-6	PTM	2,000 30425
<i>p</i> -xylene	106-42-3	DMSO	1.52mg/mL 36416
<i>p</i> -xylene	106-42-3	PTM	2,000 30426
Zectran	315-18-4	Neat	500mg PS-756
zinc phosphate	—	Neat	1g PS-803
Zineb	12122-67-7	Neat	1g PS-23
Zinophos	297-97-2	Neat	100mg PS-668
ziram	137-30-4	Neat	1g PS-21

*Volume is 1mL/ampul unless otherwise noted. Concentration is $\mu\text{g/mL}$ unless otherwise noted.

Solvent code chart available on page 438.

Cat.#s starting with PS- & MET- are ChemService pesticides that are subject to size and price changes without notification. They are not available in Europe.

searching for the perfect **solution?**

Let Restek's Analytical Reference Materials Team
create the right solution for you!

Restek should be your first choice for custom-made reference materials. Maximum convenience, maximum value, minimum time spent blending calibration mixtures in your laboratory.

- Quick quotations.
- Mixtures made to your EXACT specifications.
- Most reference materials shipped within 5-10 working days after receipt of your order*

We have more than 2,000 pure, characterized, neat compounds in our inventory!
If you do not see the **EXACT** mixture you need listed on any of the following pages, contact us.

* Availability of raw materials and final product testing requested may affect delivery of some mixtures.



Analytical Reference Materials Group

Mike Maye, Rebecca Lauver, Randy Fisher, John Lidgett, Joe Tallon,
Joe Moodler, Andrea Gill, Marty Rockey



Jessie Emel, Diane Shaffer, Staci Bodle, Sara Eyster, Mary Ellen Wood,
Brendan Conway, Matt Hepfer (not pictured: Sam Morehead)