

Pesticides Calibration Mixtures

Components of these products are at 16x the Contract Required Quantitation Level (CRQL) and can be used to prepare calibration mixes at 4x CRQL and at 1x CRQL by serial dilution.

Pesticide Standard Mix A w/Surrogates

(11 components)

α -BHC	8 μ g/mL	endosulfan I	8
γ -BHC (lindane)	8	endrin	16
4,4'-DDD	16	heptachlor	8
4,4'-DDT	16	methoxychlor	80
decachlorobiphenyl (SUR)	16	2,4,5,6-tetrachloro-	
dieldrin	16	<i>m</i> -xylene (SUR)	8

In hexane, 1mL/ampul

cat. # 32003 (ea.)

Pesticide Standard Mix B w/Surrogates

(13 components)

aldrin	8 μ g/mL	endosulfan sulfate	16
β -BHC	8	endrin aldehyde	16
δ -BHC	8	endrin ketone	16
α -chlordane	8	heptachlor epoxide	
γ -chlordane	8	(isomer B)	8
4,4'-DDE	16	2,4,5,6-tetrachloro-	
decachlorobiphenyl (SUR)	16	<i>m</i> -xylene (SUR)	8
endosulfan II	16		

In hexane, 1mL/ampul

cat. # 32004 (ea.)

Technical Chlordane, Toxaphene Solutions

Compound	cat.# (ea.)
1,000 μ g/mL in hexane, 1mL/ampul	
chlordane (technical)	32021
toxaphene	32005
2,000 μ g/mL in methanol, 1mL/ampul	
chlordane (technical)	32016
toxaphene	32015
5,000 μ g/mL in isooctane, 1mL/ampul	
chlordane (technical)	32072
toxaphene	32071

Aroclor Solutions

Compound	cat.# (ea.)
1,000 μ g/mL in hexane, 1mL/ampul	
Aroclor 1016	32006
Aroclor 1221	32007
Aroclor 1232	32008
Aroclor 1242	32009
Aroclor 1248	32010
Aroclor 1254	32011
Aroclor 1260	32012
Aroclor 1016/1260	32039
Aroclor 1262	32409
Aroclor 1268	32410
200 μ g/mL in isooctane, 1mL/ampul	
Aroclor 1016	32064
Aroclor 1221	32065
Aroclor 1232	32066
Aroclor 1242	32067
Aroclor 1248	32068
Aroclor 1254	32069
Aroclor 1260	32070
Aroclor 1016/1260	32299
400 μ g/mL in acetone, 1mL/ampul	
Aroclor 1016/1260	32456
50mg/kg in transformer oil (PCB-free)	
Aroclor 1016	32075
Aroclor 1221	32077
Aroclor 1232	32079
Aroclor 1242	32081
Aroclor 1248	32083
Aroclor 1254	32085
Aroclor 1260	32087
500mg/kg in transformer oil (PCB-free)	
Aroclor 1016	32076
Aroclor 1221	32078
Aroclor 1232	32080
Aroclor 1242	32082
Aroclor 1248	32084
Aroclor 1254	32086
Aroclor 1260	32088

please **note**

We test our transformer oil solvent to ensure that it is PCB-free.