

Viva C8 Columns (USP L7)

Chromatographic Properties:

Highly base-deactivated wide pore packing that exhibits excellent peak shape for a wide range of compounds. Less retention in reversed phase assays than Viva C18.

Length	1.0mm ID		2.1mm ID		3.2mm ID		4.6mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
5µm Columns								
30mm	9513531		9513532		9513533		9513535	
50mm	9513551		9513552		9513553		9513555	
100mm	9513511		9513512		9513513		9513515	
150mm	9513561		9513562		9513563		9513565	
200mm	9513521		9513522		9513523		9513525	
250mm	9513571		9513572		9513573		9513575	

Viva C8 Guard Cartridges

Guard Cartridges	3-pk. (10 x 2.1mm)	3-pk. (10 x 4.0mm)	2-pk. (20 x 2.1mm)	2-pk. (20 x 4.0mm)	price
Viva C8 Guard Cartridge	951350212	951350210	951350222	951350220	

Viva C8 HPLC Prep Columns

Length	10mm ID		21.2mm ID		30mm ID		50mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
5µm Columns								
50mm	9513557		9513558		9513559		9513550	
100mm	9513517		9513518		9513519		9513510	
150mm	9513567		9513568		9513569		9513560	
250mm	9513577		9513578		9513579		9513570	

Viva C4 Columns (USP L26)

Chromatographic Properties:

Highly base-deactivated wide pore packing that exhibits excellent peak shape for a wide range of compounds. Less retention in reversed phase assays than Viva C18 or Viva C8.

Length	1.0mm ID		2.1mm ID		3.2mm ID		4.6mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
5µm Columns								
30mm	9512531		9512532		9512533		9512535	
50mm	9512551		9512552		9512553		9512555	
100mm	9512511		9512512		9512513		9512515	
150mm	9512561		9512562		9512563		9512565	
200mm	9512521		9512522		9512523		9512525	
250mm	9512571		9512572		9512573		9512575	

Viva C4 Guard Cartridges

Guard Cartridges	3-pk. (10 x 2.1mm)	3-pk. (10 x 4.0mm)	2-pk. (20 x 2.1mm)	2-pk. (20 x 4.0mm)	price
Viva C4 Guard Cartridge	951250212	951250210	951250222	951250220	

Viva C4 HPLC Prep Columns

Length	10mm ID		21.2mm ID		30mm ID		50mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
5µm Columns								
50mm	9512557		9512558		9512559		9512550	
100mm	9512517		9512518		9512519		9512510	
150mm	9512567		9512568		9512569		9512560	
250mm	9512577		9512578		9512579		9512570	



Physical Characteristics:

particle size: 5µm, spherical
pore size: 300Å
carbon load: 5%
endcap: yes
pH range: 2.5 to 8
temperature limit: 80°C

also available

3µm particles are available for all Viva phases—please inquire.




Physical Characteristics:

particle size: 5µm, spherical
pore size: 300Å
carbon load: 3.5%
endcap: yes
pH range: 2.5 to 8
temperature limit: 80°C

ordering note

Other column dimensions and phases are available. Please call for a quote.



free literature

Viva Wide Pore HPLC Columns

Download your free copy from www.restek.com

lit. cat# 59939