

### Dual Vespel® Ring Cross-Disk Inlet Seals

for Agilent GCs

- Ideal for high-flow split applications.
- Washerless, leak-tight seals.

0.8mm ID Dual Vespel Ring Cross-Disk Inlet Seal	2-pk./price	10-pk./price
Gold-Plated	22083	22084
Siltek Treated	22085	22086



### Replacement Inlet Seals with Washers

- Siltek® treatment provides inertness similar to fused silica.
- All seals include washers.

#### Replacement Inlet Seals

for Agilent GCs

The inlet seal at the base of the Agilent 5890/6890/7890 GC injection port contacts the sample and must be changed frequently to prevent adsorption of active compounds. In addition, septum fragments and sample residue accumulate on the disk surface, requiring disk replacement.

To reduce breakdown and adsorption of active compounds, use Siltek® or gold-plated seals.



Now available in economical 50-pks.

Single-Column Installation 0.8mm ID (Opening)	Similar to Agilent part #	2-pk./price	10-pk./price	50-pk./price
Gold-Plated	18740-20885	21317	21318	23415
Siltek Treated		21319	21320	23416
Stainless Steel	18740-20880	21315	21316	23417

0.25/0.32mm ID Dual-Column Installation 1.2mm ID (Opening)	2-pk./price	10-pk./price
Gold-Plated	21305	21306
Siltek Treated	21307	21308
Stainless Steel	20390	20391

0.53mm ID Dual-Column Installation 1/16-inch ID (Opening)	2-pk./price	10-pk./price
Stainless Steel	20392	20393



Note: The 1.2 mm inlet seal is recommended when installing two columns using a two-hole Vespel®/graphite ferrule.

### Replacement Inlet Seal Washers

Description	Similar to Agilent part #	qty.	cat.#	price
Replacement Inlet Seal Washers	5061-5869	15-pk.	21710	



### Cross-Disk Inlet Seals

for Agilent GCs

- Ideal for high-flow split applications on Agilent 5890 GCs.
- All seals include washers.

0.8mm ID Cross-Disk Inlet Seal	Similar to Agilent part #	2-pk./price	10-pk./price
Gold-Plated	5182-9652	20477	20476
Siltek Treated	—	20475	20474

1.2mm ID Cross-Disk Inlet Seal	Similar to Agilent part #	2-pk./price	10-pk./price
Gold-Plated	—	21009	\$71 21010
Siltek Treated	—	21011	\$71 21012



### Silver PTV Seals

for Agilent 6890 GCs

Description	Similar to Agilent part #	qty.	cat.#	price
Silver PTV Seals for Agilent 6890 GCs	5182-9763	5-pk.	21409	

