



FID Collector Assembly Kit

for Agilent 5890 GCs

- Constructed of high-quality stainless steel.
- Meets or exceeds original manufacturer's performance.
- Individual replacement parts available.

Description	Similar to Agilent part #	qty.	cat.#	price
FID Collector Assembly Kit (includes insulators)	19231-60690	kit	23010	
FID Collector Assembly Kit w/Siltek Ignitor Castle		kit	21131	
FID Collector (includes insulators)	19231-20970, 19231-20960, 19231-20950, 19231-20940,	ea.	21138	
FID Collector Nut and Washer	5181-3311	set	21136	
FID Ignitor*	19231-60680	ea.	21001	
FID Ignitor Castle	19231-20910	ea.	21137	
Siltek FID Ignitor Castle		ea.	21135	

*Also fits OI Analytical 4410 detector (similar to OI part # 191833).

FID Collector Assembly Kit

for Agilent 6890/6850/7890 GCs

- Constructed of high-quality stainless steel.
- Meets or exceeds original manufacturer's performance.
- Individual replacement parts available.

Description	Similar to Agilent part #	qty.	cat.#	price
FID Collector Assembly Kit (includes insulator)	G1531-60690	kit	21699	
FID Collector Assembly Kit w/Siltek Ignitor Castle		kit	21132	
FID Collector (includes insulators)	G1531-20690, G1531-20700, 19231-20940,	ea.	21139	
FID Collector Nut and Washer	5181-3311	set	21136	
FID Ignitor*	19231-60680	ea.	21001	
FID Ignitor Castle	19231-20910	ea.	21137	

*Also fits OI Analytical 4410 detector (similar to OI part # 191833).



FID Flow Measuring Adaptor

for Agilent 5890/6890/6850/7890 GCs

- Makes setting flows easy.
- Meets or exceeds original manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#	price
FID Flow Measuring Adaptor	19301-60660	ea.	21000	

FID Gauge Pack

for Agilent 5890 GCs

- Pressure regulators and gauges for air & hydrogen.
- 1/8-inch bulkhead allows easy hookup to instrument.
- Rated for inlet pressures to 250 psi (1,724 kPa).
- Rated for outlet pressures of 0 to 60 psi (0-414 kPa).



Description	qty.	cat.#	price
FID Gauge Pack for Agilent 5890 GCs	ea.	22071	

FID Maintenance Kits

for Agilent GCs

- Include the most common consumable GC supplies and tools.
- All parts meet or exceed performance of instrument manufacturer's parts.
- Parts list makes reordering easy.

FID kits include:

- 1/4", 0.4, 0.5, and 0.8 mm ID graphite ferrules.
- FID/NPD capillary adaptor.
- Capillary nuts.
- Jet reamers/ferrule removers.
- 1/4" nut.
- Scoring wafer.
- Capillary column caps.
- Ignitor for either Agilent 5890 or 6890/6850/7890 GCs.
- FID flow measuring adaptor.
- 1/4" x 5/16" wrench.
- Installation gauge.
- Wire cleaning brush.
- High-performance Siltek® treated FID jet for either Agilent 5890 (adaptable jet) or 6890/6850/7890 (dedicated jet) GCs.
- Spanner wrench.
- FID jet removal tool.



Description	qty.	cat.#	price
FID Maintenance Kit for Agilent 5890 GCs	kit	22180	
FID Maintenance Kit for Agilent 6850/6890/7890 GCs	kit	22179	