




Liners for APEX ProSep 800 & ProSep 800 Plus GCs

APEX Liners for ProSep 800 & ProSep 800 Plus GCs	Benefits/Uses:	ID*	Similar to APEX part #	cat.#/price ea.
		OD x Length		
 Mega IV (4.0mm ID)	injections ≤125µL	4.0mm 6.0mm x 243mm	L-00410	21075
 Micro I (1.0mm ID)	injections ≤5µL	1.0mm 6.0mm x 243mm	L-00110	21073
 MIDI II (2.0mm ID)	injections ≤25µL	2.0mm 6.0mm x 243mm	L-00210	21074

COLUMN INSTALLS THIS END

*Nominal ID at syringe needle expulsion point.

Accessories for APEX GCs

Viton® O-Rings for Apex Liners



Description	Max. temp.	qty.	cat.#	price
O-Rings for APEX liners	250°C	25-pk.	22067	

Deactivated Wool

This deactivated wool is more inert than our traditional fused silica wool. This wool can be used to vaporize a sample in a liner prior to introduction into a capillary column.



Description	qty.	cat.#	price
Deactivated Wool	10 grams	24324	

restek innovation!

Inlet Liner Removal Tool

- Easily remove liner from injector—no more burned fingers.
- Made from high-temperature silicone.
- Won't chip or crack the liner.



20181

Description	qty.	cat.#	price
Inlet Liner Removal Tool	3-pk.	20181	







20117

Septum Puller

- Use hooked end for removing septa and O-rings; pointed end works well for removing stuck ferrules or debris.
- Keep several on hand in your laboratory—can be used in many different ways.

Description	qty.	cat.#	price
Septum Puller	ea.	20117	

Liners for ATAS Injectors

Liners for ATAS Injectors	Benefits/Uses:	ID*	Similar to ATAS part #	ea.	cat.#/price 5-pk.
		OD x Length			
	universal	3.0mm 5.0mm x 80mm		22415	\$25
 ATAS Open Liner, 3mm	trace, active samples <1µL	1.0mm 5.0mm x 80mm	A100049	22417	\$42
 ATAS Open Liner, 1mm		3.0mm 5.0mm x 80mm	A100126	22419	\$40
 ATAS Fritted Gooseneck	dirty samples	3.0mm 5.0mm x 80mm	A100126	22419	\$40

COLUMN INSTALLS THIS END

Liners for ATAS LINEX DMI System OPTIC 2 and 3 inlets



ATAS LINEX DMI Liner

ID*	cat.#/price 5-pk.
OD x Length	
3.3mm 5.0mm x 81mm	22353



*Nominal ID at syringe needle expulsion point.

Accessories for ATAS LINEX DMI System



22354

22355

Description	qty.	cat.#	price
Needle Guide for ATAS LINEX DMI Liner	5-pk.	22354	
Microvial/Insert for ATAS LINEX DMI Liner	100-pk.	22355	

ordering note

For additional liner ordering information, see pages 206-213.