



### Teflon® Tubing

- Ideal for mobile phase inlet lines.
- Chemically inert.
- Use to 500 psi (3,447 kPa) and 80 °C.

Description	OD	ID	Length	qty.	cat.#	price
Teflon Tubing	1/8"	0.063" (1.6mm)	3m	3m	25306	
Teflon Tubing	1/8"	0.094" (2.4mm)	3m	3m	25307	



20188

### Tubing Dressing Tool

Opens stainless steel tubing bore and removes burrs. For 1/16-inch OD tubing or 1/8-inch OD tubing.

Description	Size	qty.	cat.#	price
Tubing Dressing Tool	1/16"	ea.	20188	
Replacement Insert	for 1/16" Tubing Dressing Tool	ea.	20189	
Tubing Dressing Tool	1/8"	ea.	20190	
Replacement Insert	for 1/8" Tubing Dressing Tool	ea.	20191	



### Tubing Clip

Securely holds 1/16-inch or 1/8-inch tubing in beaker, flask, or bottle up to 4 mm thick.

Description	qty.	cat.#	price
Tubing Clip	5-pk.	25310	



25069

### Clean-Cut Tubing Cutter

- Burr-free, perpendicular cuts that will not distort the tubing OD or close the ID.
- Use on PEEK™, Teflon®, Tefzel®, other polymeric tubing.

Description	qty.	cat.#	price
Clean-Cut Tubing Cutter	ea.	25069	
Replacement Blade for Clean-Cut Cutter	ea.	25070	

For stainless steel tubing cutters, see page 324 or visit [www.restek.com](http://www.restek.com).



25308

25309

### PEEK™ Tubing Elbows

Tubing Elbows (90° and 180°) are ideal for routing 1/16-inch PEEK™ tubing through your system. Simply snap the tubing into the elbow. Prevent pinching of PEEK™ tubing which can cause high pressure.

Description	qty.	cat.#	price
PEEK Tubing Elbow, 90°	5-pk.	25308	
PEEK Tubing Elbow, 180°	5-pk.	25309	