

Syringes for LC Autosamplers

SGE for Agilent 1090 & 1100 LC Autosamplers

Volume	Model	SGE cat.#	qty.	Restek cat.#	price
25µL, 1/4-32 UNEF Thread	25D-HP1090-GT	003670	ea.	22290	



22290

SGE for Hitachi LC Autosamplers

Volume	Model	SGE cat.#	qty.	Restek cat.#	price
500µL, M10 X1 Thread	500C-HITACH1	007660	ea.	22292	



22292

SGE for PerkinElmer LC Autosamplers

Volume	Model	SGE cat.#	qty.	Restek cat.#	price
100µL, 1/4-28 UNF Thread	100D-CX-GT	005990	ea.	22296	
250µL, 1/4-28 UNF Thread	250D-CX-GT	006995	ea.	22297	
500µL, 1/4-28 UNF Thread	500D-CX-GT	007995	ea.	22298	
1mL, 1/4-28 UNF Thread	1MD-C-GT	008185	ea.	22299	



22297

Syringes for Waters WISP LC Autosamplers

Teflon® PTFE-tipped plungers

Hamilton

Volume	Model	Hamilton cat.#	qty.	Restek cat.#	price
25µL, 1/4-28 UNF Thread	1702	80020	ea.	24528	
250µL, 1/4-28 UNF Thread	1725	80024	ea.	24529	



24528

SGE

Volume	Model	SGE cat.#	qty.	Restek cat.#	price
25µL, 1/4-28 UNF Thread	25D-WISP	003990	ea.	22293	
250µL, 1/4-28 UNF Thread	250D-WISP	006690	ea.	22294	



22294

Syringes for CTC LC Autosamplers

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	Hamilton cat.#	qty.	Restek cat.#	price
10µL	N	22s	2"/51mm	3	701N	203073	ea.	22743	
25µL	N	22s	2"/51mm	3	1702N	203075	ea.	22744	
100µL	N	22s	2"/51mm	3	1710N Slim Line*	203077	ea.	22745	
100µL	N	22	2"/51mm	3	1710N	203235	ea.	22746	
250µL	N	22	2"/51mm	3	1725N	203079	ea.	22747	
500µL	N	22	2"/51mm	3	1750N	203349	ea.	22748	



22746

*Barrel OD = 6.7mm; all other 25µL and 100µL syringes have a 7.9mm barrel OD.

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	SGE cat.#	qty.	Restek cat.#	price
10µL	F	22s	2"/51mm	LC	10F-CTC-LC	002710	ea.	22737	
100µL*	R	22s	2"/51mm	LC	100R-C/F-GT-LC	005330	ea.	22741	
500µL*	F	22s	2"/51mm	LC	500F-CTC-GT-LC(0.41)	007720	ea.	22742	



22737

*Gas-tight syringe.