



24385

24383



21724

also available

See page 366 for crimpers and decappers.



24673



21776

24518A

24692



24651

2.0 mL Crimp-Top Vials, 12 x 32 mm, 11 mm

White graduated marking spots are a standard feature unless otherwise noted.

2.0 mL Crimp-Top Vials, 12 x 32 mm, 11 mm

Description	Volume	Color	Material	100-pk.	1000-pk.
Vial w/White Graduated Marking Spot	2.0mL	Clear	Glass	24383	\$18.50 24384
Vial w/White Graduated Marking Spot	2.0mL	Amber	Glass	24385	\$19.50 24386
Vial without Graduated Marking Spot	2.0mL	Clear	Glass	21152	\$16.50 21153

2.0 mL, 11 mm Aluminum Crimp Seals with Septa

Seal Color	Septa Material	100-pk.	500-pk.	1000-pk.
Silver	PTFE/Natural Rubber	21174	\$13.50	— 21175
Blue	PTFE/Red Rubber	24351	\$19.50	— 24352
Green	PTFE/Red Rubber	24353	\$19.50	— 24354
Red	PTFE/Red Rubber	24355	\$19.50	— 24356
Yellow	PTFE/Red Rubber	24357	\$19.50	— 24358
Mixed	PTFE/Rubber	—	— 21724	\$50 —
Silver	PTFE/Silicone	24359	\$25.50	— 24360
Blue	PTFE/Silicone	24361	\$28	— 24362
Green	PTFE/Silicone	24363	\$28	— 24364
Red	PTFE/Silicone	24365	\$28	— 24366
Yellow	PTFE/Silicone	24367	\$28	— 24368
Mixed	PTFE/Silicone	—	— 21725	\$119.50 —
Silver	PTFE/Silicone/PTFE	24369	\$28	— 24370

2.0 mL Crimp Vial Convenience Kits (Vials, Caps, & Septa)

Vials packaged in a clear-lid tray. Caps with septa packaged in a plastic bag.

Description	100-pk.	1000-pk.
Clear 2.0mL Vial, Deactivated, Silver Seal, PTFE/Natural Rubber Septa	24671	\$44 24672
Amber 2.0mL Vial, Deactivated, Silver Seal, PTFE/Natural Rubber Septa	24673	\$44 24674
Clear 2.0mL Vial, Untreated, Silver Seal, PTFE/Natural Rubber Septa	21196	\$34 21197
Amber 2.0mL Vial, Untreated, Silver Seal, PTFE/Natural Rubber Septa	21198	\$34 21199
Clear 2.0mL Vial, Untreated, Silver Seal, PTFE/Silicone Septa	24646	\$40 24647
Amber 2.0mL Vial, Untreated, Silver Seal, PTFE/Silicone Septa	24648	\$40 24649

Inserts

for 2.0 mL, 11 mm Crimp-Top & 2.0 mL, 9 mm Short-Cap, Screw-Thread Vials

Description	Volume	Material	100-pk.	1000-pk.
Big Mouth Insert w/Bottom Spring	50µL	Glass	24513	\$48 21782
Big Mouth Insert w/Bottom Spring	250µL	Glass	21776	\$82 21777
Big Mouth Insert w/Glass Flange (Step Design)*	250µL	Glass	24516	\$71.50 21779
Insert, Flat Bottom	350µL	Glass	21780	\$21.50 24517
Insert, Flat Bottom w/ID Ring	350µL	Glass	24692	\$23.50 24693
Big Mouth Insert w/Bottom Spring	250µL	Polypropylene	24518	\$18.50 —
Big Mouth Insert w/Bottom Spring & Graduated Markings	250µL	Polypropylene	24518A	\$24.50 —
Big Mouth Insert, Top Flange	250µL	Polypropylene	24519	\$11.50 —
Big Mouth Insert, No Spring	250µL	Polypropylene	24520	\$11.50 —

*Not to be used with 9mm screw-thread vials.

2.0 mL Crimp-Top Vial with 100 µL Insert

- Plastic vial.
- Choice of glass or polypropylene insert.

Description	100-pk.	1000-pk.
Clear Plastic Vial w/100µL Polypropylene Insert	24651	\$21.50 24652
Clear Plastic Vial w/100µL Glass Insert	24653	\$149 —