



24654

4.0 mL WISP 48 Screw-Thread Step Vials (15 x 45 mm, 13/425)

4.0 mL WISP 48 Screw-Thread Step Vials (vials only)

Description	Volume	Material	100-pk.	1000-pk.
WISP 48 Step Vial w/White Graduated Marking Spot	4.0mL	Clear	24654	24655
WISP 48 Step Vial w/White Graduated Marking Spot	4.0mL	Amber	24656	24657

4.0 mL WISP 48 Caps and Septa

Description	Septa Material	100-pk.	1000-pk.
WISP 48 Caps and Septa, Top Hat	PTFE/Silicone	21745	21746
WISP 48 Caps and Septa, Preassembled	Red PTFE/Silicone, 0.065"	21743	21744
WISP 48 Caps, Black Polypropylene, Open-Hole		21741	21742
Septa	Red PTFE/Silicone/PTFE, 0.040"	21751	21752
Septa	Red PTFE/Silicone, 0.060"	21749	21750

Inserts

for 4.0 mL WISP 48 Screw-Thread Step Vials

Description	Volume	Material	100-pk.	1000-pk.
Insert for WISP 48 Vials, Flat Bottom**	500µL	Glass	21787	21788

Polypropylene inserts available on request (1,000-packs only).

**Also fit 4.0mL Snap Seal Vials. See page 361.



24662

24664

24678

1.0 mL WISP 96 Style Shell Vials, 8 x 40 mm

Available in 3 styles.

1.0 mL WISP 96 Style, 8 x 40 mm Shell Vials

Description	Volume	Color	Material	200-pk.	1000-pk.
WISP 96 Style Vial	1.0mL	Clear	Glass	24662	24663
WISP 96 Style Vial	1.0mL	Amber	Glass	24664	24665
WISP 96 Style Vial	1.0mL	Clear	Polypropylene	24678	24679
WISP 96 Style Vial	700µL Limited Volume	Clear	Polypropylene	24680	24681

8mm Starburst Plugs sold separately.



21769

21773

8 mm Starburst Plugs

for Shell Vials, 1.0 mL, WISP 96 Style, 8 x 40 mm

Description	Color	Material	200-pk.	1000-pk.
Starburst Snap Plugs for Shell Vials	Clear	Polyethylene	21769	21770
Starburst Snap Plugs for Shell Vials, Conical	Clear	Polyethylene	21771	21772
Starburst Snap Plugs for Shell Vials, Silicone-Lined	Clear	Polyethylene	21773	21774

Limited Volume Inserts

for Shell Vials, 1.0 mL, WISP 96 Style, 8 x 40 mm

Description	Volume	Material	100-pk.	1000-pk.
Big Mouth Insert w/Bottom Spring	250µL	Glass	21776	21777



21776

did you know?

Siltek® coated vials create a highly inert surface for your testing purposes, making them inert to most organics. See page 353 for more details.

instrument reference

For 1.0 mL WISP 96 Style Vials

The following instruments are compatible with 1.0 mL WISP 96 style vials. For further information or questions, please contact your Restek sales representative or technical service.

Manufacturer	Instrument/Model #
Hitachi	L-7200
Waters	96-position M-700, 710, 717+