



Universal plug set for battery charger included in all kits.

Electronic Rechargeable Crimpers and Decappers

- Easy to use; comfortable grip.
- Hundreds of operations from one charge.
- Adjustable crimping force.
- New battery allows for faster charge cycle.
- All kits include a universal plug set for the battery charger.
- For 11 mm and 20 mm caps.
- One-year warranty.

The redesigned electronic crimpers and decappers have a 4.8 volt battery installed. The battery allows for quick, powerful operation and a faster charge cycle, minimizing down-time. Charging the battery takes thirty minutes to two hours. The universal plug set includes U.S., European, UK, and Australian plug clips. Each kit includes a tool, a 4.8 volt battery, which allows for a faster charge cycle, and a universal plug set for the battery charger.

Capping

Easily crimp caps on standard crimp vials.



Decapping has never been easier!

one...

two...

three!



Description	Size	qty.	cat.#	price
Electronic Crimper	11mm	kit	22338	
Electronic Crimper	20mm	kit	22339	
Electronic Decapper	11mm	kit	22340	
Electronic Decapper	20mm	kit	22341	
Replacement Battery for Electronic Crimpers and Decappers		ea.	22342	

Safer than removing caps with pliers, knives, or screwdriver blades!



Crimper



Decapper

Easy Grip Manual Crimpers and Decappers

Our new Easy Grip Crimpers and Decappers feature tough, light glass-filled plastic bodies with curved handles, which improve hand comfort during use. The bottom-pull handle allows for a steadier hold than top-push handles. An easily viewed adjustment knob is located at the head of the product. The clearly shown + and - symbols, along with directional arrows, simplify adjusting the intensity of the crimp/decap needed. The knob is also designed to indicate when the crimp setting has been reached, when it fully contacts the body of the product.

Description	Size	qty.	cat.#	price
Crimper	11mm	ea.	23396	
Decapper	11mm	ea.	23397	
Crimper	20mm	ea.	23398	
Decapper	20mm	ea.	23399	