

No more burned fingers!



The Claw and The Claw Holder Kit

- Easily removes hot liners from injection ports.
- 4 mL vials (not included) can be replaced when dirty.

Never again will you have to burn your fingers removing a hot injection port liner. The Claw safely and cleanly removes liners, O-rings, or other small objects from the injection port. You can then place the hot objects in a clean 4 mL vial situated in the Claw holder until ready for reuse.

Description	qty.	cat. #	price
The Claw	ea.	26261	
The Claw Holder Kit (includes The Claw and holder)	kit	26262	
WISP 48 Snap Seal Vial	100-pk.	24658	

2.0 mL Sample Vial Racks

- Racks feature alphanumeric indexing for easier vial identification.
- Racks can be stacked for efficient storage.
- Fits most 2.0 mL, 12 x 32 mm vials (vials sold separately).

Description	Capacity	qty.	cat.#	price	qty.	cat.#	price
Polypropylene Storage Rack for 12 x 32mm Vials	50 vials, 5 x 10	ea.	22856	\$17.50	5-pk.	22857	



Custom Standards

Restek is your #1 source for custom reference standards!

- **Made to your specifications.**
- **Quick quotations.**
- **Most reference standards shipped within 15 working days.**
- **The more you buy, the more you save.**

www.restek.com/solutions