



Autosampler Syringes for Thermo Scientific GCs

- Hamilton and SGE syringes are designed and tested to meet critical autosampler performance needs.
- Needle point styles are designed to withstand multiple, rapid injections through a septum.

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	Hamilton cat.#	Restek cat.#	price
10µL	N	26s	3.15"	AS	701SN	80318	24962	ea.

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	SGE cat.#	Restek cat.#	price
10µL	F	26	80mm	Cone	10F-C/F-8/0.47C	002992	24924	ea.
10µL	R	26	80mm	Cone	10R-C/F-8/0-0.47C	002993	24934	ea.



Replacement Needles for Thermo Scientific Autosampler Syringes

SGE

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	Restek cat.#	price
0.5µL	26	50mm	Cone	033770	24927	ea.
10µL	26	80mm	Cone	031535	24935	3-pk.



also available

Looking for Restek Genuine Replacement Parts for Thermo Scientific Instruments?
See pages 253-258.

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) Teflon® Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip

instrument reference

Autosampler Syringes for Thermo Scientific GCs

Autosampler syringes for TriPlus, AS3000, AS2000, AS200/800 GCs. The following autosamplers are compatible with the respective SGE and Hamilton syringes. For further information or questions, contact your Restek sales representative or technical service.

SGE Syringes

TriPlus	AS3000	AS2000	AS200/800	Restek cat.#	Page
•	•			24928	377
•	•			24926	377
•	•	•	•	24923	377
•	•	•	•	24932	377
•	•	•	•	24924	376
•	•	•	•	24934	376
•	•	•	•	24922	377
•	•	•	•	24930	377

Replacement Needles for SGE Syringes for Thermo Scientific GCs

TriPlus	AS3000	AS2000	AS200/800	Restek cat.#	Page
•	•			24929	377
•	•			24927	376
•	•	•	•	24933	377
•	•	•	•	24931	377
•	•	•	•	24935	376

Hamilton Syringes

TriPlus	AS3000	AS2000	AS200/800	Restek cat.#	Page
			•	24962	376