

Q-sep™ QuEChERS Extraction Salts

- Salt packets eliminate the need for a second empty tube to transfer salts.
- Go green by using packets with reusable tubes.
- Convenient and easy to use.



Description	Material	Methods	qty.	cat#	price
Q110 Kit	4g MgSO ₄ , 1g NaCl, 1g TSCD, 0.5g DHS with 50mL Centrifuge Tube	European EN 15662	50 packets & 50 tubes	26235	
Q110 Packets	4g MgSO ₄ , 1g NaCl, 1g TSCD, 0.5g DHS	European EN 15662	50 packets	26236	
Q150 Kit	6g MgSO ₄ , 1.5g NaOAc with 50mL Centrifuge Tube	AOAC 2007.01	50 packets & 50 tubes	26237	
Q150 Packets	6g MgSO ₄ , 1.5g NaOAc	AOAC 2007.01	50 packets	26238	
Empty 50mL Centrifuge Tube, Polypropylene			50-pk.	26239	
Empty 50mL Centrifuge Tube, Teflon FEP			2-pk.	23997	

TSCD—trisodium citrate dihydrate
 DHS—disodium hydrogen citrate sesquihydrate
 NaOAc—sodium acetate



Q-sep™ QuEChERS dSPE Tubes for Extract Clean-Up

- Fast, Simple Sample Prep for Multi-Residue Pesticide Analysis
- Packaged in mylar subpacks for enhanced protection and storage stability.
 - Individually labeled tubes for easy sorbent identification.

Description	Material	Methods	qty.	cat#	price
2mL Micro-Centrifuge Tubes for dSPE (clean-up of 1mL extract)					
Q210	150mg MgSO ₄ , 25mg PSA	European EN 15662	100-pk.	26215	
Q211	150mg MgSO ₄ , 25mg PSA, 25mg C18		100-pk.	26216	
Q212	150mg MgSO ₄ , 25mg PSA, 2.5mg GCB	European EN 15662	100-pk.	26217	
Q213	150mg MgSO ₄ , 25mg PSA, 7.5mg GCB	European EN 15662	100-pk.	26218	
Q250	150mg MgSO ₄ , 50mg PSA	AOAC 2007.01	100-pk.	26124	
Q251	150mg MgSO ₄ , 50mg PSA, 50mg C18	AOAC 2007.01	100-pk.	26125	
Q253	150mg MgSO ₄ , 50mg PSA, 50mg GCB		100-pk.	26123	
Q252	150mg MgSO ₄ , 50mg PSA, 50mg C18, 50mg GCB	AOAC 2007.01	100-pk.	26219	
15mL Centrifuge Tubes for dSPE (clean-up of 6mL extract)					
Q350	1200mg MgSO ₄ , 400mg PSA	AOAC 2007.01	50-pk.	26220	
Q351	1200mg MgSO ₄ , 400mg PSA, 400mg C18	AOAC 2007.01	50-pk.	26221	
Q352	1200mg MgSO ₄ , 400mg PSA, 400mg C18, 400mg GCB	AOAC 2007.01	50-pk.	26222	
Q370	900mg MgSO ₄ , 150mg PSA	European EN 15662	50-pk.	26223	
Q371	900mg MgSO ₄ , 150mg PSA, 15mg GCB	European EN 15662	50-pk.	26224	
Q372	900mg MgSO ₄ , 150mg PSA, 45mg GCB	European EN 15662	50-pk.	26225	
Q373	900mg MgSO ₄ , 150mg PSA, 150mg C18		50-pk.	26226	
Q374	900mg MgSO ₄ , 300mg PSA, 150mg GCB		50-pk.	26126	

PSA—primary and secondary amine exchange material
 GCB—graphitized carbon black



Sorbent Guide

Sorbent	Removes
PSA	sugars, fatty acids, organic acids, anthocyanine pigments
C18	lipids, nonpolar interferences
GCB	pigments, sterols, nonpolar interferences
PSA	—primary and secondary amine exchange material
GCB	—graphitized carbon black

Try QuEChERS risk-free today!

Call 800-356-1688 or your Restek representative to request a free sample pack of Q-sep™ QuEChERS tubes.

(Sample pack orders cannot be placed on line. I limit one pack per customer.)

