

**Resprep® 12- or 24-Port SPE Manifolds**

- Use with any standard male luer-end SPE cartridges.
- Inert, Teflon® sample guides reduce cross-contamination and carryover.
- Flexible sample collection rack will accommodate a variety of receiving vessels.
- Quick vacuum-release valve for better system control.
- Individual vacuum control for each cartridge—improves reproducibility.



**Resprep® Manifold Replacement Parts**

Description	qty.	cat.#	price
<b>Complete Resprep SPE Manifold, 12-Port</b> (Includes: glass basin with built-in vacuum regulator, polypropylene lid with 12 individual control valves, 12-position collection rack, 12 Teflon sample guides, and waste container)	kit	26077	
<b>Complete Resprep SPE Manifold, 24-Port</b> (Includes: glass basin with built-in vacuum regulator, polypropylene lid with 24 individual control valves, 24-position collection rack, and 24 Teflon sample guides)	kit	26080	

**Resprep® Manifold Replacement Parts**

Description	qty.	cat.#	price
Replacement Waste Container, 12-Port	ea.	24014	
Replacement Vacuum Valve and Gauge Assembly	ea.	24008	
Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Port	ea.	25991	
<b>Collection Rack</b>			
Collection Rack, 12-Port	ea.	26079	
Collection Rack, 24-Port	ea.	26082	
<b>Manifold Lid</b>			
Replacement Manifold Lid (sample guides not included), 12-Port	ea.	26078	
Replacement Manifold Lid (sample guides not included), 24-Port	ea.	26081	
<b>Manifold Lid Replacement Parts</b>			
Valves, 12 or 24-Port	2-pk.	26083	
Valve Attachment	48-pk.	26130	
Needle Attachment	48-pk.	26131	
Sample Guides, 12 or 24-Port	12-pk.	26084	

\*Not included in 24-port manifold (cat.# 26080 or 26082).

**Oil-Free Vacuum/Pressure Station**

for 12-Port Resprep® SPE Systems

Allows switching from pressure to vacuum in a matter of seconds. Quiet, oil-free unit won't contaminate the extraction system or your samples. Provides 20" Hg (68k Pa) vacuum or 18 psi (124 kPa) pressure.

Description	qty.	cat.#	price
Oil-Free Vacuum/Pressure Station, 115VAC, 60Hz, US	ea.	24002	
Oil-Free Vacuum/Pressure Station, 230VAC, 50Hz, Europe (CE certified)	ea.	24003	
Vacuum Tubing (10ft./3m, 1/4" ID)	ea.	24016	

Not recommended for use with 24-port manifold.

Warranty period is one-year from date of purchase. Evaluation fee is charged for repairs out of warranty.

