

General Purpose Columns



Chemically bonded capillary columns

- Allow for direct solvent injection onto column.
- Columns can be solvent rinsed.

Extensive GC column selection

- Available in many dimensions, including variations in length, internal diameter, and film thickness.
- Internal diameters include 0.10mm and 0.18mm for faster analysis time and greater resolution.

Broad range of stationary phases

- Columns based on polysiloxane backbone; functional groups added to the polymers to vary selectivity:



Rtx[®]-1, Rtx[®]-5, Rtx[®]-5MS, Rtx[®]-20, Rtx[®]-35, Rtx[®]-50, Rtx[®]-65, Rtx[®]-440, Rtx[®]-200, Rtx[®]-200MS, Rtx[®]-1301, Rtx[®]-624, Rtx[®]-1701, Rtx[®]-225, Rtx[®]-2330, Rtx[®]-Wax, Stabilwax[®]

visit www.restek.com for complete product listings