

SGE Syringe Needles

Replacement Needles

Needle Selection

- needle length refers to the exposed length of the needle when fitted to syringe
- choose a different replacement needle by selecting the syringe volume, needle length, OD and tip style
- for example, a 10 μ L MicroVolume syringe (Part No. 002050) is supplied with a 50mm bevel tip needle (Part No. 037110). This can be replaced with any needle designed for 10 μ L syringes

[5 \$\mu\$ L Syringes](#)
 [10 \$\mu\$ L Syringes](#)
 [25-500 \$\mu\$ L Syringes](#)
 [1-2.5mL Syringes](#)
 [5-10mL Syringes](#)
 [Luer Lock](#)
 [Valve](#)
[Plunger-in-Needle Syringes](#)

5 μ L Syringes

Part no. ▲▼	Needle Length (mm) ▲▼	Gauge ▲▼	OD (mm) ▲▼	ID (mm) ▲▼	Tip Style ▲▼	Pack Size ▲▼	
036011	50	23	0.63	0.11	Cone	2	<input type="button" value="Buy"/>
036110	50	26	0.47	0.11	Bevel	5	<input type="button" value="Buy"/>

[Top](#)

10 μ L Syringes

Part no. ▲▼	Needle Length (mm) ▲▼	Gauge ▲▼	OD (mm) ▲▼	ID (mm) ▲▼	Tip Style ▲▼	Pack Size ▲▼	
037010	50	26	0.47	0.11	Cone	2	<input type="button" value="Buy"/>
037011	50	23	0.63	0.11	Cone	2	<input type="button" value="Buy"/>
037030	70	26	0.47	0.11	Cone	2	<input type="button" value="Buy"/>



This website uses cookies

This website uses cookies to enhance the functionality of our website, to personalise content, to provide social media features, and to analyse our traffic. You can click on the "Show details" button to read more about the cookies used by this site. Please indicate whether you consent to the placement of the different types of cookies by ticking the boxes below (other than Necessary cookies which are necessary for the operation of our site).

<input type="checkbox"/> Necessary	<input type="checkbox"/> Preferences	<input type="checkbox"/> Statistics	<input type="checkbox"/> Marketing	<input type="button" value="Show details"/>	<input type="button" value="OK"/>
------------------------------------	--------------------------------------	-------------------------------------	------------------------------------	---	-----------------------------------

037610	100		0.17	0.1	On-column	2	<input type="button" value="Buy"/>
037675	75		0.23	0.1	On-column	2	<input type="button" value="Buy"/>

[Top](#)

25-500 μ L Syringes








Part no. ▲▼	Needle Length (mm) ▲▼	Gauge ▲▼	OD (mm) ▲▼	ID (mm) ▲▼	Tip Style ▲▼	Pack Size ▲▼	
038010	50	25	0.5	0.2	Cone	2	<input type="button" value="Buy"/>
038030	70	25	0.5	0.2	Cone	2	<input type="button" value="Buy"/>
038038	180	23	0.63	0.24	Cone	2	<input type="button" value="Buy"/>
038060	115	25	0.5	0.2	Cone	2	<input type="button" value="Buy"/>

SGE Syringe Needles

038110	50	25	0.5	0.2	Bevel	5	Buy
038111	50	23	0.63	0.24	Bevel	5	Buy
038130	70	25	0.5	0.2	Bevel	5	Buy
038131	70	23	0.63	0.24	Bevel	5	Buy
038138	180	23	0.63	0.24	Bevel	2	Buy
038160	115	25	0.5	0.2	Bevel	5	Buy
038161	115	23	0.63	0.24	Bevel	5	Buy
038250	51(2")	22	0.028"	0.37	LC	5	Buy
038260	50	25	0.5	0.2	LC	5	Buy
038270	70	22	0.028"	0.37	LC	5	Buy
038410	50	25	0.5	0.2	Side Hole/Dome	2	Buy
038430	70	25	0.5	0.2	Side Hole/Dome	2	Buy
038736	50	23	0.63	0.24	Side Hole/Dome	2	Buy








[Top](#)

1-2.5mL Syringes

Part no. 	Needle Length (mm) 	Gauge 	OD (mm) 	ID (mm) 	Tip Style 	Pack Size 	
039038	180	23	0.63	0.32	Cone	2	Buy
039060	115	23	0.63	0.32	Cone	5	Buy
039110	50	23	0.63	0.32	Bevel	5	Buy
039115	50	22	0.028"	0.37	Bevel	2	Buy
039116	50	23	0.63	0.32	Dome	2	Buy
039120	50	23	0.63	0.32	Side Hole/Dome	2	Buy
039130	70	23	0.63	0.32	Bevel	5	Buy
039138	180	23	0.63	0.32	Bevel	2	Buy
039160	115	23	0.63	0.32	Bevel	5	Buy
039161	70	23	0.63	0.32	Cone	2	Buy
039250	51(2")	22	0.028"	0.37	LC	5	Buy
039260	50	25	0.5	0.2	LC	5	Buy

[Top](#)

5-10mL Syringes

Part no. 	Needle Length (mm) 	Gauge 	OD (mm) 	ID (mm) 	Tip Style 	Pack Size 	
031516	50	23	0.63	0.32	Bevel	5	Buy
0315233	51(2")	22	0.028"	0.37	LC	5	Buy
0315234	50	25	0.5	0.2	LC	5	Buy

[Top](#)

Luer Lock Needles

SGE Syringe Needles

Part no.	Needle Length (mm)	Gauge	OD (mm)	ID (mm)	Tip Style	Pack Size	<input type="button" value="Buy"/>
039802	50	23	0.63	0.32	Bevel	5	<input type="button" value="Buy"/>
039803	50	23	0.63	0.32	Side Hole/Dome	2	<input type="button" value="Buy"/>
039807	70	23	0.63	0.32	Bevel	5	<input type="button" value="Buy"/>
039808	70	23	0.63	0.32	Side Hole/Dome	2	<input type="button" value="Buy"/>
039811	115	23	0.63	0.32	Bevel	5	<input type="button" value="Buy"/>
039822	50	19	1.07	0.65	Bevel	5	<input type="button" value="Buy"/>
039823	50	19	1.07	0.65	Side Hole/Dome	2	<input type="button" value="Buy"/>
039827	70	19	1.07	0.65	Bevel	5	<input type="button" value="Buy"/>
039831	115	19	1.07	0.65	Bevel	5	<input type="button" value="Buy"/>
039842	50	18	1.27	0.8	Bevel	5	<input type="button" value="Buy"/>
039847	70	18	1.27	0.8	Bevel	5	<input type="button" value="Buy"/>
039851	115	18	1.27	0.8	Bevel	5	<input type="button" value="Buy"/>
039862	50	16	1.57	1.19	Bevel	5	<input type="button" value="Buy"/>
039867	70	16	1.57	1.19	Bevel	5	<input type="button" value="Buy"/>
039871	115	16	1.57	1.19	Bevel	5	<input type="button" value="Buy"/>
039880	50	14	2.1	1.6	Bevel	5	<input type="button" value="Buy"/>
039891	115	14	2.1	1.6	Bevel	5	<input type="button" value="Buy"/>
039895	51(2")	22	0.028"	0.37	LC	2	<input type="button" value="Buy"/>
039897	50	25	0.5	0.2	LC	2	<input type="button" value="Buy"/>

[Top](#)

Valve Needles

Part no.	Needle Length (mm)	Gauge	OD (mm)	ID (mm)	Tip Style	Pack Size	<input type="button" value="Buy"/>
038810	50	23	0.63	0.2	Cone	1	<input type="button" value="Buy"/>
038815	50	23	0.63	0.2	Side Hole/Dome	1	<input type="button" value="Buy"/>

[Top](#)

Plunger-in-Needle kits

- kits are supplied with both a plunger and a needle
- plunger and needle must be replaced together

MicroVolume Syringes, 0.5uL, 6.5mm OD Barrel

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style	<input type="button" value="Buy"/>
033010	50	23	0.63	Cone	<input type="button" value="Buy"/>
033011	50	23	0.63	Bevel	<input type="button" value="Buy"/>
033012	50	26	0.47	Cone	<input type="button" value="Buy"/>

MicroVolume Syringes, 0.5uL, 8mm OD Barrel

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style
----------	--------------------	-------	---------	-----------

SGE Syringe Needles

033010	50	23	0.63	Cone	Buy
033057	70	23	0.63	Cone	Buy
033060	70	23	0.63	Bevel	Buy
033610	100	0	0.23	Sheathed / LC	Buy
033620	75	0	0.23	LC	Buy
033630	70	26	0.47	Cone	Buy

MicroVolume Syringes, 1uL, 8mm OD Barrel

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style	
034055	50	23	0.63	Cone	Buy
034056	50	23	0.63	Bevel	Buy
034057	70	23	0.63	Cone	Buy
034059	115	23	0.63	Cone	Buy
034060	70	23	0.63	Bevel	Buy
034610	70	26	0.47	Cone	Buy

MicroVolume Syringes, 5uL, 8mm OD Barrel

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style	
035055	50	23	0.63	Cone	Buy
035056	50	23	0.63	Bevel	Buy
035057	70	23	0.63	Cone	Buy
035058	70	23	0.63	Bevel	Buy

Autosampler Plunger-in-Needle Syringes

Agilent Technologies - 7673, 7683 & 6850 ALS

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style	
033708	42	26	0.47	Cone	Buy
033715	42	23	0.63	Cone	Buy
033730	42	23-26	0.63/0.47	Cone	Buy
034715	42	23	0.63	Cone	Buy

CTC

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style	
033770	50	26	0.47	Cone	Buy
033772	50	23	0.63	Cone	Buy
034905	50	23	0.63	Cone	Buy

Perkin Elmer - AutoSystem & Clarus 500

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style	
----------	--------------------	-------	---------	-----------	--

SGE Syringe Needles

033750	70	26	0.47	Cone/Bevel	<input type="button" value="Buy"/>
033765	70	23	0.63	Cone	<input type="button" value="Buy"/>

Shimadzu - AOC 14, 17, 20 & 20i

Part no.	Needle Length (mm)	Gauge	OD (mm)	Tip Style	<input type="button" value="Buy"/>
033738	42	26	0.47	Cone	<input type="button" value="Buy"/>
033745	42	23	0.63	Cone	<input type="button" value="Buy"/>
033747	42	23-26	0.63/0.47	Cone	<input type="button" value="Buy"/>

[Top](#)

Replacement PTFE Seals and Spring

Part no.	Description	<input type="button" value="Buy"/>
032030	Spring for removable needles (pack of 10)	<input type="button" value="Buy"/>

© Trajan Scientific Australia Pty Ltd 2016 [Company Profile](#) [Locations](#) [Careers](#) [Portal](#) [Login](#) ·