

**PEEK one-piece fingertight fittings – color-coded**

10-32

**SPECS**


- For 1/16" OD tubing
- Material: PEEK
- 10-32 threads
- Tolerance: ±0.05 mm (.002")
- Pressure rating: <350 bar (5000 psi) (Varies with tubing material and ID)

- Fingertight up to 350 bar (5000 psi)
- Excellent chemical and mechanical properties
- One-piece – no ferrule required
- Biocompatible

Our PEEK fingertight fittings are available in six colors, for easy identification of tubing lines. They are rated to 350 bar (5000 psi) and can be used in virtually any type of 10-32 thread HPLC fitting detail on the market.

Sold in packages of 5.

| Color   | Product No. |
|---------|-------------|
| Natural | JR-55020-5  |
| Black   | JR-55021-5  |
| Red     | JR-55022-5  |
| Yellow  | JR-55023-5  |
| Blue    | JR-55024-5  |
| Green   | JR-55025-5  |



**PEEK one-piece fingertight fittings**

10-32


**SPECS**

- For 1/16" OD tubing
- Material: PEEK
- 10-32 threads
- Tolerance: ±0.05 mm (.002")
- Pressure rating: <350 bar (5000 psi) (Varies with tubing material and ID)

- Natural PEEK versions fingertight up to 350 bar (5000 psi)
- One-piece – no ferrule required
- Biocompatible
- Four styles, machined for extra precision


Sold in packages of 5.

| Description                 | Product No. |
|-----------------------------|-------------|
| Standard                    | JR-5502-5   |
| Short                       | JR-5504-5   |
| Long                        | JR-5507-5   |
| Narrow hex-head/fingertight | JR-5508-5   |



Standard

Short



Long

Narrow hex-head/fingertight

**TECH TIP**

Coloration pigments are carefully chosen so as to contain no heavy metals or hazardous chemicals.

**MORE INFORMATION**

About PEEK ..... 99  
 About CTFE ..... 99  
 Chemical resistance 96-98