



2009 #60

Low Pressure Unions

Unions *Cheminert to Cheminert* 1/4-28 to 1/4-28

Includes flangeless 1/4-28 fittings for tubing OD indicated. Polypropylene unions are for use with flanged tubing only.


Tubing OD	Bore	PEEK		CTFE		Polypropylene	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CUCPK		CUCKF		-	
1/16"	0.50 mm	CUPK		CUKF		-	
1/16"	0.75 mm	CUMPK		CUMKF		-	
1/8"	1.50 mm	CULPK		CULKF		-	
1/8"	Butt connection	CUTPK		CUTKF		CUTPP *	(* pkg/5)



Unions *Cheminert to 1/16" ZDV* 1/4-28 to 10-32

Includes flangeless 1/4-28 and ZDV 10-32 fittings for 1/16" tubing.

Tubing OD	Bore	PEEK		CTFE		316 Stainless	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CZUCPK		CZUCKF		CZUCS6	
1/16"	0.50 mm	CZUPK		CZUKF		CZUS6	
1/16"	0.75 mm	CZUMPK		CZUMKF		CZUMS6	



Unions *Cheminert to 1/4" tubing* 1/4-28 to 1/2-20

Includes flangeless 1/4-28 and 1/2-20 fittings.

Tubing OD	Bore	PEEK		CTFE	
		Prod No	Price	Prod No	Price
1/8" to 1/4"	1.50 mm	CU4LPK		CU4LKF	

Components	Prod No	Price
1/2-20 nut, CTFE	CFL-4KF	
1/2-20 nut, Delrin	CFL-4D	
CTFE ferrule	CFL-CB4KF-S	




0.25 mm = .010"
0.50 mm = .020"
0.75 mm = .030"
1.0 mm = .040"
1.5 mm = .060"
2.0 mm = .080"
4.6 mm = .180"
6.0 mm = .236"
6.4 mm = .253"
7.0 mm = .275"
10.0 mm = .400"
27.0 mm = 1.08"
1/32" = 0.8 mm
1/16" = 1.6 mm
1/8" = 3.2 mm
1/4" = 6.4 mm
3/8" = 9.5 mm
1/2" = 12.7 mm



Bulkhead unions *Cheminert to Cheminert* 1/4-28 to 1/4-28

Includes flangeless 1/4-28 fittings for tubing OD indicated.

Tubing OD	Bore	PEEK		CTFE		316 Stainless	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CBUCPK		CBUCKF		CBUCS6	
1/16"	0.50 mm	CBUPK		CBUKF		CBUS6	
1/16"	0.75 mm	CBUMPK		CBUMKF		CBUMS6	
1/8"	1.50 mm	CBULPK		CBULKF		CBULS6	



Bulkhead unions *Cheminert to 1/16" ZDV* 1/4-28 to 10-32

Includes flangeless 1/4-28 and ZDV 10-32 fittings for 1/16" OD tubing.

Tubing OD	Bore	PEEK		CTFE		316 Stainless	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CZBUCPK		CZBUCKF		CZBUCS6	
1/16"	0.50 mm	CZBUPK		CZBUKF		CZBUS6	
1/16"	0.75 mm	CZBUMPK		CZBUMKF		CZBUMS6	

