

PEEK Flangeless Nuts

PEEK Flangeless Nuts

- Flange-free
- Compatible with all 1/4"-28 flat bottom ports
- Excellent chemical resistance
- Fits 1/16" or 1/8" OD tubing
- True fingertight – no tools required

PEEK Flangeless Nuts in combination with the inverted ferrule do not require flanging tools. A hand tightening procedure makes leak free seals. The inert

ETFE ferrule is the only part, besides the tubing, which is in contact with the fluid – the nut is never in contact with the fluid stream (see illustration).



PEEK Flangeless Nuts

Part No.	Description	Qty/pkg
JR-20124-10	Nut, PEEK, flangeless, natural, 1/16"	10
JR-20116-10	Nut, PEEK, flangeless, natural, 1/8"	10
JR-041-10	Ferrule, ETFE, 1/16"	10
JR-051-10	Ferrule, ETFE, 1/8"	10

SPECS

Material
Nuts: PEEK
Ferrules: ETFE

Dimensions
For 1/16" and 1/8" OD tubing
Others: see illustration

Threads:
1/4"-28

Pressure rating
1/16" up to 105 bar (1500 psi)
1/8" up to 35 bar (500 psi)
Varies with tubing material and ID