



VICI JOUR - Tubing for Chromatography

VICI JOUR PRODUCTS

- Tubing for chromatography
- High pressure fittings
 - Nanovolume
 - Stainless steel
 - Polymeric
- Low pressure fittings
- Filters & mobile phase filters
- Mobile phase handling systems
- Frits
- PEEK columns & guard columns
- Sample loops
- Accessories/miscellaneous

CUSTOMER SUPPORT

VICI Jour offers an assortment of chromatography grade tubing in 1/32", 1/16", 1/8", and 1/4" ODs, in PEEK, PEEK-clad fused silica, electroformed nickel, stainless, PTFE, and more. You'll also find unique tubing accessories and tools.



Tubing

- PEEK tubing
- PEEKclad-FST™
- Premium grade stainless steel tubing
- Premium grade stainless steel tubing, pre-cut
- Stainless steel tubing
- Electroformed nickel tubing
- No-Ox tubing
- PTFE, PFA, FEP, or ETFE tubing
- Thermo-formed polymeric tubing
- Knitted PTFE reaction coils

Tools and Accessories

- PEEK tubing elbows
- Tubing Clip – the LC tubing organizer
- Clean Cut polymer tubing cutter
- Guillotine polymer tubing cutter
- Stainless steel tubing cutter
- Stainless steel tubing pliers
- Easy-Flange and Easy-Flange Combi Kits

Disclaimer

The maximum holding pressure for any type of connection involving tubing and a ferrule varies considerably with the tubing material, the ferrule material, the clearance between tubing OD and ferrule ID, and the shape of the fitting detail.

TOP ▲