

- 001 (various, see part description)
- 06 Orange Nylon
- 07 Yellow Nylon
- 0 Natural Nylon
- 1 White Nylon
- 2 Black Nylon
- 3 Red Nylon
- 4 Green Nylon
- 5 Blue Nylon
- 7 Glass Filled Natural Nylon
- 8 Gray Nylon
- 8B Gray Nylon
- 9 Clear Polycarbonate
- 10 Acetal
- 40 Polysulfone
- 81 Lace White ABS
- 91 White Polycarbonate
- 1006 Natural Acetal
- 1007 Black Acetal
- 6005 Animal-Free Polypropylene
- 6006 White Polypropylene
- 8003 Natural ABS
- 8012 White ABS
- 9002 Rad. Stable Polycarbonate
- 9010 Rad. Stable Polycarbonate
- B Buna-N O-ring
- CM030 Dark Blue Nylon
- E PTEE/EPDM
- EP EPDM O-ring
- J1A Natural Kynar PVDF
- J1A2 Black Kynar PVDF
- PX Platinum-Cured Silicone
- U71 White Class VI Glass Filled Nylon
- U75 Blue Class VI Glass Filled Nylon
- VP1 Animal-free Polyethylene
- X2 Buna-N Rubber
- XO Stainless Steel

**CF104 (-6005)**  
Male Luer Integral Lock Ring Plug, Ø.033 Orifice at Luer Tip



**CF226 (-6005)**  
High Flow Connector to 500 Series Barb, 1/4" (6.4 mm) ID Tubing (May be used with separate rotating lock ring; CF227)



**CF015 (-1)**  
Female Luer Lug Style Flare to M6x.75 Thread



**CF131 (-1)**  
Male Luer Integral Lock Ring Plug, Closed at Grip (Short Luer)




**CF141 (-9)**  
Male Luer to Female Luer Thread Style Coupler (Ø.045 Thru Bore with Ø.055 Male Side Bond-in Port)



**CF148 (-9002)**  
Male Luer to Female Luer Thread Style Coupler (Ø.045 Thru Bore with Ø.080 Male Side Bond-in Port)



**CF210 (-6005)**  
Male Luer to Female Luer Thread Style Coupler (Ø.125 Center Bond-in Port)



**CF213 (-3)**  
Male Luer Integral Lock Ring Plug, Closed at Grip (Ø.055 to Ø.060 Core)



**CF224 (-9)**  
Male Luer Integral Lock Ring to 500 Series Barb, 5/32" (4 mm) ID Tubing (Ø.076 Bore)



**CF233 (-6005)**  
Female Extended Luer Thread Style to 200 Series Barb, 1/16" (1.6 mm) ID Tubing



**CF235 (-9 -91)**  
Male Luer to Female Luer Thread Style Coupler (Ø.023 Thru Bore with Ø.033 Male Side Bond-in Port)



**CF239 (-6005)**  
Male Luer Integral Lock Ring to Female Luer Thread Style Extended Coupler



**CF322 (-1)**  
Female Luer Lug Style Flare to 1/4-32 UNF Thread



- Pg 82** Chemical Resistance
- Pg 83** Regulatory Status
- Pg 84** Sterilization Stability
- Pg 85** Typical Properties